


PCT

WORLD INTELLECTUAL PROPERTY ORGANIZATION
International Bureau



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁶ : G06F 15/42		A1	(11) International Publication Number: WO 99/35588
			(43) International Publication Date: 15 July 1999 (15.07.99)
(21) International Application Number: PCT/US98/22830			(81) Designated States: AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CZ, DE, DK, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IS, JP, KE, KG, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, US, UZ, VN, YU, ZW, ARIPO patent (GH, GM, KE, LS, MW, SD, SZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG).
(22) International Filing Date: 29 October 1998 (29.10.98)			
(30) Priority Data: 60/071,107 12 January 1998 (12.01.98) US 09/143,483 28 August 1998 (28.08.98) US			
(71) Applicant (for all designated States except US): SANGSTAT MEDICAL CORPORATION [US/US]; 1505 Adams Drive, Menlo Park, CA 94025 (US).			
(72) Inventors; and (75) Inventors/Applicants (for US only): POULETTY, Philippe [FR/US]; 25 Oakhill Drive, Woodside, CA 94062 (US). HAMILTON, Richard, G. [US/US]; 4670 Amiens Avenue, Fremont, CA 94555 (US). ROSSI, Stephen, J. [US/US]; 655 Kansas Street #301, San Francisco, CA 94107 (US). McENROE, Debra, L. [US/US]; 566 Cutwater Lane, Foster City, CA 94404 (US).			
(74) Agents: CARLSON, Stephen, C. et al.; McDermott, Will & Emery, Suite 300, 99 Canal Center Plaza, Alexandria, VA 22314 (US).			Published With international search report. Before the expiration of the time limit for amending the claims and to be republished in the event of the receipt of amendments.
(54) Title: SYSTEM AND METHOD FOR MANAGING ADMINISTRATION OF MEDICINE			
(57) Abstract			
<p>A method for managing doses of medication delivered to a patient is described. A computer system (101) receives dosage data and administration data that represent, respectively, times and quantities for taking a drug that are prescribed for a patient, and the times and quantities the drug is delivered to the patient. Based on the dosage and administration data, compliance information is generated and displayed, representing the degree to which a drug has been delivered in accordance with the dosage data. In one aspect, a calendar (126) in the form of a grid comprised of grid elements is displayed. Each grid element represents a period, such as a day in a month, and contains one or more icons. An icon's appearance indicates whether a particular dose was delivered properly, when a grid element is selected by a user, more detail is displayed about the administration of the drug for the respective day.</p>			
			

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AL	Albania	ES	Spain	LS	Lesotho	SI	Slovenia
AM	Armenia	FI	Finland	LT	Lithuania	SK	Slovakia
AT	Austria	FR	France	LU	Luxembourg	SN	Senegal
AU	Australia	GA	Gabon	LV	Latvia	SZ	Swaziland
AZ	Azerbaijan	GB	United Kingdom	MC	Monaco	TD	Chad
BA	Bosnia and Herzegovina	GE	Georgia	MD	Republic of Moldova	TG	Togo
BB	Barbados	GH	Ghana	MG	Madagascar	TJ	Tajikistan
BE	Belgium	GN	Guinea	MK	The former Yugoslav Republic of Macedonia	TM	Turkmenistan
BF	Burkina Faso	GR	Greece			TR	Turkey
BG	Bulgaria	HU	Hungary	ML	Mali	TT	Trinidad and Tobago
BJ	Benin	IE	Ireland	MN	Mongolia	UA	Ukraine
BR	Brazil	IL	Israel	MR	Mauritania	UG	Uganda
BY	Belarus	IS	Iceland	MW	Malawi	US	United States of America
CA	Canada	IT	Italy	MX	Mexico	UZ	Uzbekistan
CF	Central African Republic	JP	Japan	NE	Niger	VN	Viet Nam
CG	Congo	KE	Kenya	NL	Netherlands	YU	Yugoslavia
CH	Switzerland	KG	Kyrgyzstan	NO	Norway	ZW	Zimbabwe
CI	Côte d'Ivoire	KP	Democratic People's Republic of Korea	NZ	New Zealand		
CM	Cameroon	KR	Republic of Korea	PL	Poland		
CN	China	KZ	Kazakhstan	PT	Portugal		
CU	Cuba	LC	Saint Lucia	RO	Romania		
CZ	Czech Republic	LI	Liechtenstein	RU	Russian Federation		
DE	Germany	LK	Sri Lanka	SD	Sudan		
DK	Denmark	LR	Liberia	SE	Sweden		
EE	Estonia			SG	Singapore		

WO 99/35588

PCT/US98/22830

**SYSTEM AND METHOD FOR MANAGING ADMINISTRATION OF MEDICINE
RELATED APPLICATIONS**

SYSTEM AND METHOD FOR MANAGING ADMINISTRATION OF MEDICINE RELATED APPLICATIONS

This application claims priority from prior U.S. provisional application serial number 60/071,107 filed on January 12, 1998, entitled "Method and System for Monitoring Doses," which is incorporated by reference in its entirety as if fully set forth herein.

COPYRIGHT NOTICE

A portion of the disclosure of this patent document contains material that is subject to copyright protection. The copyright owner has no objection to the facsimile reproduction in paper copies by anyone of the patent document or the patent disclosure, as it appears in the Patent & Trademark Office patent file or records, but otherwise reserves all other copyright rights whatsoever.

FIELD OF THE INVENTION

The present invention relates generally to computer systems. The invention relates more specifically to managing administration of medicine, monitoring dosages of drugs given to patients, and the like.

BACKGROUND OF THE INVENTION

Monitoring dosages of drugs or medicines for patients requires communication among several levels. First, a physician must diagnose and prescribe a dosage for a patient. The medication must then be distributed accurately and, finally, the patient or a care provider must ensure that the dosages are properly administered to or taken by the patient.

For many reasons, ensuring that accurate dosages are delivered to a patient in a consistently timely manner can be difficult despite the importance of accurate administration in many instances.

Therefore, it is desirable to provide a method of automating the delivery of medicine and monitoring the delivery of medicine.

Moreover, special challenges are presented in managing patients who are taking more than one medication. Elderly patients on multiple medications may have difficulty keeping track of whether they have taken all their medications, when, and in what quantity. In the clinical setting, proper administration of multiple medications to acutely ill patients is challenging for care providers.

Thus there is a need to track multiple medications and multiple dispensing mechanisms, and to present data for all such dispensers in a report.

To facilitate the proper administration of medication and the tracking of when it is administered, medication dispensing devices are used. Conventional medication dispensing devices typically include a medicine container and an alarm mechanism which notifies a patient at the time intervals the dose(s) are due. Each time the patient opens the container, the device records the event and the time it occurred. One example of a conventional medication dispensing device is a jar with lid which incorporates an alarm mechanism and a recording mechanism. When the lid is removed, the recording mechanism records this event and the time it occurred.

One drawback to conventional dispensing devices is that they do not control access to medicine or the quantities dispensed. Thus, there is little assurance that when a dispensing device is opened, the proper amount is dispensed. Another drawback is that once opened, the dispensing devices may be re-opened immediately. Thus a confused elderly patient, having
5 forgotten the dose they just took, may take another far too soon.

Thus, there is further need for a system that controls dispensing times and amounts and which tracks those times and amounts.

SUMMARY OF THE INVENTION

The foregoing needs, and other needs and objects that will become apparent from the following discussion, are fulfilled by the present invention, which comprises, in one aspect, a
10 method for managing doses of medication delivered to a patient. Generally, a computer system receives dosage data and administration data. The dosage data represents a drug prescription, and includes, but is not limited to, one or more times for taking the drug, the quantities in which the drug is to be taken by the patient, or a combination thereof. The administration data
15 represents when and in what quantities each dose in a set of doses of the drug is actually delivered to the patient. Based on the dosage and administration data, compliance information is generated and displayed. Compliance information indicates the degree to which a drug has been delivered in accordance with the dosage data. The compliance information can be displayed in variety of forms.

20 According to another aspect, a calendar in the form of a grid comprised of grid elements is displayed. Each grid element represents a period, such as a day in a month, and contains one or more icons. An icon's appearance indicates whether a particular dose was delivered properly. For example, a green square icon indicates that a dose was delivered on time, and a triangular red icon indicates that a dose was not delivered. When a user selects a
25 grid element, more detail is displayed about the delivery of the drug for the respective day. In particular, a graphical object is displayed that contains one or more icons for each dose delivered in the day. An icon's position along an axis of the graphical object reflects when a dose was delivered.

According to another aspect, data is generated that specifies what portion of a set of
30 doses was delivered properly. The data includes values that indicate what portions of the doses were delivered, and what proportion of doses were delivered on time.

According to another aspect, dosage data is transmitted to a dosage-dispensing device. The dosage data includes times and quantities to deliver a drug to a patient. In addition, data representing a lockout period may be transmitted. The dosage-dispensing device dispenses the
35 drug to the patient in accordance with the data transmitted to it.

BRIEF DESCRIPTION OF THE DRAWINGS

The present invention is illustrated by way of example, and not by way of limitation, in the figures of the accompanying drawings and in which like reference numerals refer to similar elements and in which:

40 FIG. 1 is a block diagram illustrating a system for monitoring patient dosages.

FIG. 2 is a flow chart illustrating steps for a computer-implemented method for monitoring patient dosages.

FIG. 3 is a flow chart showing steps for retrieving data that is used in a system for monitoring the administration of doses to a patient.

5 FIG. 4 is flow chart showing steps for transmitting dosage information to a dosage-dispensing device.

FIG. 5A is block diagram depicting a calendar in the form of a grid.

FIG. 5B is a block diagram depicting a grid element and icons used to indicate patient compliance.

10 FIG. 5C is a block diagram depicting a graphical object used to graphically represent when doses were delivered.

Fig. 6 is a block diagram depicting a histogram showing dosage scores over period of time.

DETAILED DESCRIPTION

15 OVERVIEW

One embodiment is a system and method for substantially automating the administration of patient dosages, the monitoring of the delivery of doses, whether or not timely and whether or not accurate in amount, and the accumulation of data for individual patients representing administration data over an extended period of time.

20 Another embodiment encompasses accumulation of data for each patient from a plurality of dosage dispensing devices, and the assimilation of such data into reports which may be either specific for the particular patient, or an accumulation of data for an entire range of patients. In this way, more accurate dispensing of doses is achieved, as well as more accurate monitoring to facilitate detection of whether prescribed doses are being properly
25 administered to the patients.

A preferred embodiment provides a computer-implemented method for monitoring patient dosages by retrieving administration data, including times and amounts of medication prescribed for a patient, retrieving patient data, including times and amounts of medication delivered for the patient, determining evaluation data by analyzing the retrieved dosing and
30 patient data to determine compliance of the delivered medication to the prescribed medication, and displaying the evaluation data.

The method may include one or more of the following features. Patient data, including administration data, may be received from an associated device over a communications line, from local memory, or from user input. The data may be accumulated to provide a basis for
35 patient evaluation. The patient data may be transmitted to a dosage-dispensing device, which dispenses doses to the patient in accordance with the received patient data.

The evaluation data may be displayed in a variety of ways, including display in a patient administration report that may indicate compliance of the delivered doses to the prescribed dosages. In one implementation, the data retrieved may be viewed in a scrollable
40 tabular grid, with displayed values for all medication events, and dates, times and dose sizes

dispensed from the dosage dispensing device. In addition, non-medication events may be displayed, including "bottle replaced" or other ancillary but relevant data.

Additionally, evaluation data may be displayed in the form of a patient summary report which may, for example, include all information for a particular patient including name, ID, monitoring dates, drug, brand, and so on. In addition, a histogram may be prepared summarizing the patient's compliance, including calculation of a "compliance index" or similar quantification of the patient's overall compliance with the prescribed dosing plan. The evaluation data may be displayed for varying periods, such as a week, a month, or a shorter or longer period, and may be displayed in graphical form including options for displaying doses delivered, missed, or delivered but not within compliance parameters. The data may also be displayed in calendar form.

In many instances, patients undergoing treatment may have multiple dosage dispensers. In a manner similar to the single dispenser arrangement discussed above, data for each such dispenser can be tracked and presented in a merged patient summary report. Likewise, a summary of all patients may be provided which may provide, in either graphical or tabular form, any of the selected data including name, ID, compliance index, dosage, time of day, or any other field. Histograms may also be developed across the patient class.

Evaluation data may be provided in any suitable format, such as a data file or hard copy. For example, the data may be printed or transmitted to a remote facsimile machine.

According to one embodiment, the delivery of doses of multiple patients is monitored. In this embodiment, a preferred method comprises retrieving dosage data, including times and amounts of medication prescribed for a plurality of patients, retrieving patient data, including times and amounts of medication delivered for the plurality of patients, determining evaluation data by analyzing the retrieved dosing and patient data for the plurality of patients to determine overall compliance of the delivered medication to the prescribed medications, and displaying the evaluation data.

Another embodiment includes a memory device storing computer readable instructions for aiding a computer to implement a method for monitoring patient dosages such as that described above.

Yet another embodiment provides a system for monitoring patient dosages including a computer implementing a method such as that described above.

MEDICINE ADMINISTRATION MANAGEMENT SYSTEM

Embodiments of the invention may be implemented on special purpose electronic or data processing hardware, software applications running on general purpose hardware, or a combination of both. For example, an embodiment may be implemented in a dose administration system that includes a computer system running one or more application programs that provide functions for manipulating dosing and patient data, having access through appropriate communications links to remote devices.

FIG. 1 shows an illustrative system incorporating the present invention, including personal computer 101 running application software. Computer 101 has access to both dosage

data and patient data. For example, as shown in FIG. 1, the computer 101 includes a communications link 105 that couples computer 101 to dosage dispensing device 110. The dosage dispensing device 110 may be, for example, the portable medication administration device described in U.S. Patent Application Serial No. 08/ 867,010 entitled Liquid Medication Dispenser Apparatus, filed on June 2, 1997 and naming as inventors Debra L. McEnroe, Robert A. Britts, Phillippe Pouletty and Ralph Levy, the entire contents of which are hereby incorporated by reference as if fully set forth herein. Dosage dispensing device 110 may be used to dispense, for example, an analgesic drug, opiate agonist or antagonist drug, or an immunosuppressive drug, such as azathioprine, Tacrolimus, Sirolimus, mycophenolate, mofetil, and their chemical derivatives.

A portable medication administration device is a device which may be transported with the patient outside a medical facility such as a hospital or doctor's office, and which delivers multiple doses to the patient without immediate supervision by a registered medical clinician. Such dispensers are typically used by, for example, physicians and pharmacists, to input dosage data.

Communications link 105 enables the dosage data to be recorded at locations remote from the monitoring system, such as at medical facilities where medications are prescribed.

In the illustrated monitoring system, the computer 101 retrieves information relating to the patient data from data stored on diskette 120 or in a mass storage device, such as the computer's hard disk drive 122. This data typically includes a record of doses delivered to the patient and is typically created by the patient or a caretaker. As with the dosage information, this information may be input at remote locations, such as at a patient's home or a location where the medication is administered.

Of course, dosage and patient data may also be provided by alternative methods. For example, the data may be input directly by a user through the computer keyboard 102. The computer 101 can save input and retrieve information by downloading to the diskette 120 or hard drive 122, or if appropriate, may initiate to medication dispenser and monitor 109 a communications link 107. Communications link 107 may use electrical, electromagnetic, optical signals, or other signals that may carry digital data. These signals are exemplary forms of carrier waves transporting information.

Application software running on the computer 101 processes the dosage and patient data to determine monitoring information for patients. The monitoring information is provided to a user in the format of, for example, patient summary reports and graphs 124, event calendars 126, and summaries of groups of patients 128. The monitoring information can also be provided in hard copy via printer 130 or fax 132 through appropriate communication links.

Computer 110 may transmit data to dosage dispensing device 110 via communications link 105. The data may include times and quantities to administer a drug to a patient, and a value representing a lockout period. Dosage dispensing device 110 delivers a drug in accordance with the received data.

In one embodiment, computer 101 is a personal computer having an Intel or AMD-type processor and running the Microsoft® Windows 95 or Windows NT operating system, and equipped with volatile memory such as RAM and non-volatile memory such as a hard disk. A display device such as a CRT also forms part of computer 101.

5 **MONITORING ADMINISTRATION OF MEDICINE TO A PATIENT**

FIG. 2 is a diagram of a method of monitoring the administration of medicine to a patient. In one embodiment, the method of FIG. 2 is implemented in one or more application programs that are executed by computer 101.

At block 202, a computer such as personal computer 101 of the system of FIG. 1,
10 begins execution of the application software. As shown in block 210, computer 101 retrieves dosage and patient data for a patient from stored data. As indicated by block 212, the steps of block 210 may involve retrieving previously stored data files from a mass storage device such as disk drive 122.

Alternatively, computer 101 may establish an appropriate communications link, such
15 as a modem or ISDN line, to retrieve data from a remote device, such as the portable medication administration device illustrated in FIG. 1 and described in the above-referenced U.S. Patent Application Serial No. 867,010, filed June 2, 1997 and entitled Liquid Medication Dispenser Apparatus, previously incorporated by reference. In this alternative case, as indicated in block 214, the dispensing device 110 is connected to the computer 101 and
20 prepared for communication with the computer.

At block 220, dosage and medicine administration information for a patient is reviewed. Specifically, updated patient data is processed by the application software and displayed as requested by a user. The application software may be adapted to manipulate the dosage and patient information as needed. For example, the software may monitor the
25 dosages delivered to patients by recording times and amounts of doses taken by a specified patient, as indicated by the retrieved patient data. With access also to the dosage information for that patient, the software may determine, for example, compliance of a patient's delivered doses with the prescribed doses, either for specified dose times or over a period of time.

Block 220 may involve generating one or more reports, as shown by block 224. For
30 example, the method may be used to generate calendars showing the dosing events indicating, for example, the times of prescribed doses for specific patients and whether the patient complied with those doses. The method may also generate summary reports and graphs reflecting the progress of treatment for specific patients, incorporating, for example, test results. Additionally, the method may generate summary reports for groups of patients, such
35 as groups of patients taking the same medication or groups of patients of a specific physician.

The analyzed results may be stored and may be provided to a user. For example, the method may display the results on a computer monitor. Alternatively, as indicated in block 222, the computer 101 may provide hard copies of reports by printing to a printer or transmitting the results to a remote facsimile machine.

Optionally, as shown by block 230, the data is saved after it is reviewed. As indicated by block 232, the data is saved to the mass storage device from which it was retrieved.

Alternatively, as indicated in block 234, computer 101 may clear the memory of an external device from which the data was received and save a new copy of the data, or modify

5 appropriate parameters of the external device. A pre-defined format is used. For example, data read from the device 110 may be saved as one or more comma-delimited ASCII files on disk 122. Use of such a format enables the data to be human-readable, and allows the data to be imported into commercial, off-the-shelf application programs such as spreadsheets or word processors.

10 In one embodiment, the data is saved with a validation code that is computed at the time the file is saved. Whenever a saved data file is reopened, the code will be used to test and guarantee the validity of the data against corruption of the data or intentional modification by any means outside of the program. In a preferred embodiment, a relational database system such as the Microsoft Access Jet Engine is used for storing and retrieving all data.

15 At block 240, the operational sequence is complete.

RETRIEVING PATIENT DATA - INCLUDING DOSES AND TIME DELIVERED

FIG. 3 illustrates an embodiment of a method of retrieving data. FIG. 3 illustrates substeps involved in block 210 of FIG. 2 in greater detail.

At block 304, the computer system receives a request to read device data. For example, 20 block 304 may involve receiving a request to read "current patient data" that is stored in the dispensing device 110. The request may be generated in response to, for example, a user selecting a program menu option in a graphical user interface ("GUI").

As shown by block 320, the system determines whether dosage or patient data for the requested patient already exists and has not been saved since a prior retrieval operation. If 25 patient data for the requested patient already exists in memory and has not been saved during a prior retrieval, then in block 324, the system displays a prompt message to the user. The prompt message enables the user to select (1) canceling the request to retrieve patient data from the device, or (2) saving the prior data before continuing with the process of retrieving current patient data from the dosage-dispensing device. If the user chooses to cancel the 30 request to retrieve the current patient data, then execution ends. If the user chooses to save the already existing data, then control flows to block 328, where the data is saved in a user specified file. Block 328 may involve displaying a dialog box or prompt to the user that requests the user to enter a file name or pathname. Control then flows to block 330.

At block 330, the current patient data is retrieved from the dosage-dispensing device 35 and stored in a temporary buffer. The temporary buffer may be, for example, a temporary disk file or a buffer area in memory. At block 334, the data is checked to determine whether any transmission or data errors occurred during transmission from the dosage-dispensing device. For example, an 8-bit checksum algorithm can be applied to data received from a dispensing device 110 to detect errors. Such checksums are conventionally included by the dispensing 40 device 110 in data that it transmits to computer 101. If any errors are detected, then at block

338, a message to the user is displayed, informing the user that errors exist in the data, and execution ends. If no transmission errors are detected, then control flows to block 340.

As indicated by block 340, the disk or other storage device is checked to determine whether any prior patient data for the patient has been retrieved and stored. If previous data
5 has been retrieved from the device, then control flows to block 344. In this case, as shown by block 344, data for the patient is updated by merging the current patient data with the prior data. The merged data is stored in memory. A message is displayed informing the user that the merge has occurred.

As shown by block 348, the current data is stored. Alternatively, the merged data is
10 stored, if merged data was created at block 344. The user interface is updated to reflect the addition of current patient data.

At block 360, a device retrieval dialogue is displayed, which is data about the just retrieved patient data. Such data can include patient name, the drug(s), prescribed doses per day, and the administration times.

15 TRANSMITTING DOSAGE DATA TO DOSAGE DISPENSING DEVICE

In one embodiment, computer 101 transmits dosage data to dosage dispensing device
110. The dosage data is used by dosage dispensing device 110 to control the dispensing of medicine. The dosage data may represent medicine to deliver, administration times, quantities, and a lockout period. A lockout period is a period of time that must elapse after dispensing a
20 dose before another dose may be administered or delivered to the patient. The dosage data may specify medicines that include, for example, an analgesic drug, opiate agonist or antagonist drug, or a immunosuppressive drug. An example of a dosage dispensing device that receives data specifying administration times and quantities and a lock out period, and then which operates in accordance to such data, is the portable medication administration device,
25 described in U.S. Patent Application Serial No. 867,010, filed June 2, 1997 and entitled Liquid Medication Dispenser Apparatus, previously incorporated by reference.

The ability to transmit data to a dosage device that dispenses medicine accordingly provides significant advantages. The amounts of medicine that are actually dispensed to the patient may be controlled, and premature administration of doses may be prevented.

30 FIG. 4 is a diagram of a method of collecting dosage data from a user, such as a physician or other clinician, and transmitting the dosage data to a dosage dispensing device.

As shown by block 410, a request is received from a user to enter dosage data. The request may be generated in response to a user selecting a program menu option in a GUI. As indicated by block 420, current dosage data for the patient is retrieved from stored data. At
35 block 430, a data entry screen or dialog box is displayed, showing the current dosage data as the default data.

As indicated by block 440, dosage data is received from the user. The dosage data includes prescribed administration times and quantities and a lockout period. For example, the user enters the following information:

40 Number of Doses

Quantity and Unit
Times for Each Dose
Lock-out Period

As shown by block 450, the dosage data is transmitted to a dosage dispensing device, such as device 110 shown in FIG. 1. At block 460, the dosage data is stored in a mass storage device of a computer system, for example, hard disk 122 of computer 101.

In an embodiment of the present invention, the application software may be adapted to analyze additional data. This may include device monitoring data, such as the time a drug bottle was changed, temperature monitoring data, battery status, times data was downloaded from a dosage dispensing device, data identifying the bottle of the drug, such as data read from a bar code. Patient data may include test results measured at specified times to measure the effect of the administered dosages, or information on multiple drugs dispensed by a dosage dispensing device. Dosage data may include proper dosages of specified medications, as well as an indication of possible side effects and information regarding whether the dosage should be altered should those side effects be detected. In such a case, the application software may be adapted to provide an analysis of the effectiveness of the administered dosage.

EXEMPLARY GENERATION OF COMPLIANCE INFORMATION

To help determine whether a patient is administering a drug properly, compliance information is generated and displayed to a user. The system may display such compliance information in many forms. For example, the system may display a calendar that indicates whether particular doses were delivered properly. As another example, the system may display one or more compliance indexes, such as the percent of daily doses delivered or the percent of doses delivered on time. The compliance information may be generated by, for example, a computer system executing a computer program according to the source code set forth in the Appendix.

CALENDAR SHOWING PATIENT COMPLIANCE

FIG. 5A is a block diagram depicting a calendar 500. In the preferred embodiment, one or more calendars 500 are displayed to graphically convey user compliance information on a computer display, or other output device such as a printer.

Calendar 500 of FIG. 5A comprises a grid 502, which includes one or more grid elements 520. Each grid element 520 represents a particular day of the month, and may contain one or more icons 521 for each dose due on the particular day. The calendar 500 may also include a legend 523 that identifies each icon 521 with a descriptive label. Thus, each calendar 500 provides a snapshot display to the user of the dosages due for a particular patient throughout a particular month.

FIG. 5B shows grid element 520 in greater detail. Grid element 520 of FIG. 5B pertains to the second day of a particular month, as indicated by the numeric day value 540. Grid element 520 includes one or more icons 521 selected from among a new dosage icon 522, wrong time icon 524, on-time icon 526, and missed dose icon 528. The particular icons

521 that appear in a particular grid element 520 depend upon the content of the data previously entered for the patient by the user.

New dosage icon 522 is displayed so that it reflects the day the dosage was changed, as specified by, for example, dosing data retrieved from a dosage dispensing device 110. The
5 new dose size may be displayed within new dosage icon 522. For example, new dosage icon 522 may include text showing that the dosage is "250 mg".

Preferably, wrong time icon 524, and missed dose icon 528 each are displayed with different patterns that indicate whether a dose was delivered properly. For example, wrong
10 time icon 524 is a square shaped icon that is displayed in a first color, such as brown or tan, and is displayed for a dose that was delivered at the wrong time. A dose is delivered at the wrong time if it was delivered to the patient at a time outside the scheduled administration time.

Similarly, on-time icon 526 may be a green colored icon, and is displayed for a dose that was delivered on time. A dose is delivered on time if it was delivered to the patient within
15 the scheduled administration time.

Missed dose icon 528 is a circular icon displayed, for example, in red, and has a thick border. The missed dose icon 528 indicates that a patient failed to take a scheduled dose.

The colors and shapes of the icons 521 disclosed herein are not required and are not important. What is important is that a wrong time dose, on time dose, and missed dose each
20 are represented by a unique icon or symbol. In addition, another row of icons can be displayed in each grid element to indicate the number of doses due, each icon representing a scheduled dose for a day.

In one embodiment, each of the grid elements in grid 520 are graphical user controls. A user may cause the computer to display more information about a particular day reflected in
25 grid 502 by manipulating the day's respective grid element. For example, a user, using mouse 103 as an input device, moves a mouse cursor of calendar 500 onto the day's respective grid element and then clicks the mouse. In response, computer 101 displays a graphical time line with icons positioned to reflect when the drug was delivered.

FIG. 5C depicts an exemplary graphical time line. Time line 550 is a graphical image
30 having a horizontal length that reflects one 24-hour day. One or more icons 562 each represent a dose delivered for a particular day. Each of the icons 562 are displayed along the horizontal axis 564 of time line 550 so that their respective positions along the horizontal axis of time line 550 reflects when they were delivered. One or more hour labels 566 indicate the time at which a dose was delivered. For example, icon 562 represents a dose that was delivered at
35 approximately 8:00 a.m., as indicated by hour label 568.

In one embodiment, icons 562 may include icons for missed doses. Such icons may be displayed using a different pattern than those used to represent doses delivered on time. In addition, icons representing doses delivered at the wrong time can be displayed using a third pattern.

COMPLIANCE INDEXES

Compliance information can also be provided in the form of compliance indexes. A compliance index is a set of one or more values that reflects the degree to which the actual delivery of a drug complies with the prescribed administration. A variety of compliance indexes may used.

For example, the compliance indexes may include a dosage-on-time index. The dosage-on-time index reflects the percent of doses that were delivered to the patient on time in a given period. For example, assume that a drug is prescribed to be administered three times a day, at 7:00 a.m., 3:00 p.m., and 11:00 p.m., plus or minus an hour. If for a given day the drug is in fact delivered twice at 8:00 a.m. and 6:00 p.m., then the dosage-on-time index for the day is thirty-three percent (33%).

A dose-per-day index reflects the percentage of prescribed doses that were at least delivered in a given period. In the previous example, the dose-per-index would be sixty-six percent (66%) because two out of three doses were delivered in the day.

A unit-per-day-index reflects what portion of the amount of a drug prescribed for a day was delivered to the patient. For example, 2000 mg may be prescribed, but 2200 mg may be delivered to the patient. Thus, the unit-per-day-index would be 110%.

The user may specify the period covered by a compliance index in a variety of ways. For example, a graphical user control list box may provide selectable list box items which each represent a period for which to generate a compliance index. One list box item specifies the last week, another the last two weeks, and another the previous month. In addition, the graphical user control text boxes can be configured to accept the beginning and end dates of a period.

Also, various techniques may be used to display compliance indexes to the user. Each index can be displayed as a numeral, or a graphic, such as a horizontal bar. The length of the bar would represent 100 percent, and a position of an indicator along the length would indicate a percent.

One or more compliance indexes may be presented in the form of a weekly dosing graph, as shown in FIG. 5C, or in other graph forms, such as a line, area, and histogram graph. In addition, a GUI may present a graphical user control through which a user may select the form of the graph for displaying compliance indexes. For example, a GUI may display a graphical user control list box containing list box items for each graph form. By selecting one of the list box items, a user specifies a graph form for displaying a compliance index.

Fig. 6 shows a score histogram graph according to an embodiment of the present invention. Score histogram graph 600 displays patient dosing scores in the form of a graph of "Time Span" versus "Score." The time span is selectable for a time range specified by the user. The score value represents a compliance index over, for example, the last 7, 14, 21, or 28 days, or a time span specified by the user.

Score histogram graph 600 contains one or more graphical bars, such as graphical bar 610. Each graphical bar is used to reflect a dosage score for a time period within the time span,

such as a day. To measure the graphical bars, score histogram graph 600 includes graphical score scale 604. The height of the graphical bars together with graphical score scale 604 indicate a dosage score for a particular time period. Graphical bar 610 reflects a score of 66%.

OTHER REPORTS

5 Other reports can be generated based on the foregoing information.

In particular, a Patient Dosing Report is generated based on data retrieved from the dispenser device 110. The data is displayed in a scrollable tabular grid. Displayed values include all medication events, dates, times, and dose sizes that are retrieved from the dispenser. Other non-medication events that are reported from the dispenser device to the
10 computer 101 can be displayed at the option of the user. For example, when a user replaces a bottle in the dispenser, the dispenser device 110 reports a "bottle replaced" event to the computer 101. Such events can appear in the Patient Dosing Report.

As another example, a Patient Summary Report is generated. The report includes a header containing complete patient information such as Name, ID, Monitoring Dates, Drug,
15 Brand, etc.

A Patient Summary Report, based on the merged data created in block 344 of FIG. 3, can be generated. The report summarizes data downloaded from multiple devices for the same patient.

A Summary of All Patients report presents a summary of all patients in grid form. The
20 grid includes Name, ID, and Score for each patient. The grid may be sorted by any column. The Score value may be selected based on Doses Per Day or Time Of Day.

Preferably, the system provides a Print Preview function whereby the user can view any pages on the screen before they are printed.

PROGRAM STRUCTURE

25 Embodiments of the methods described further below may be implemented, for example, in one or more computer programs developed using Microsoft Visual Basic®. Preferably, the programs provide a multi-document interface whereby a user may view multiple documents simultaneously within the program. For example, the calendar dialog and medication event data dialogs described herein may be viewed at the same time.

30 In one embodiment, the program functions and method steps described above are organized in an application program using one or more pull-down menus, each of which has one or more menu options. Table 1 presents a hierarchy of menu options in one embodiment of such a program.

TABLE 1 -- MENU OPTIONS

35	FILE
	New
	Open
	Save ...
	Save As ...
40	Print Setup ...

Print Preview ...
Print ...
Exit ...
DEVICE
5 Retrieve Dispenser Data
Program Dispenser
VIEW
Dosing Data
Dosing Calendar
10 Reports & Graphs ...
HELP
About

The application program may also provide confirmation dialogs that prompt the user to verify various functions, such as dosing, as they are performed and where appropriate.

15 In the foregoing specification, the invention has been described with reference to specific embodiments thereof. It will, however, be evident that various modifications and changes may be made thereto without departing from the broader spirit and scope of the invention. The specification and drawings are, accordingly, to be regarded in an illustrative rather than a restrictive sense.

APPENDIX

CycloTech Medication Monitoring Program

SangStat Medical Corporation

Produced by Glen Hamilton, Cyber Innovations Corporation

Code Listing From
3/19/98

Table Of Contents

General.bas	1
Comm.bas	26
Printing.bas	46
Fax.bas	57
Calendar.bas	61
frmMain.frm	76
frmSplash.frm	85
frmLogin.frm	86
frmOptions.frm	88
frmAbout.frm	92
frmBrowser.frm	95
frmTip.frm	98
frmAllPatients.frm	100
frmRecentDosingGraph.frm	106
frmDosingCalendar.frm	111
frmPatientDosingRpt.frm	117
frmReadDeviceData.frm	122
frmPrint.frm	124
frmDeviceDiagnostics.frm	128
frmFaxStatus.frm	133
frmFaxSend.frm	134
frmFaxLog.frm	141
frmFaxEditGroups.frm	142
frmFaxEditLocations.frm	144
frmDeviceInitialize.frm	146
frmGetDateTime.frm	151

Attribute VB_Name = "modGeneral"
Option Explicit

'Declare DLL calls

Declare Function OSWinHelp% Lib "user32" Alias "WinHelpA" (ByVal hWnd&, ByVal HelpFile\$, ByVal wCommand%, dwData As Any)
Declare Function GetPrivateProfileInt Lib "kernel32" Alias "GetPrivateProfileIntA" (ByVal lpApplicationName As String, ByVal lpKeyName As String, ByVal nDefault As Long, ByVal lpFileName As String) As Long
Declare Function GetPrivateProfileString Lib "kernel32" Alias "GetPrivateProfileStringA" (ByVal lpApplicationName As String, ByVal lpKeyName As Any, ByVal lpDefault As String, ByVal lpReturnedString As String, ByVal nSize As Long, ByVal lpFileName As String) As Long
Declare Function WritePrivateProfileString Lib "kernel32" Alias "WritePrivateProfileStringA" (ByVal lpApplicationName As String, ByVal lpKeyName As Any, ByVal lpString As Any, ByVal lpFileName As String) As Long

'Set up some temporary buffers for getting strings from DLL calls

Public Const gIBufSize1024 = 1024 'set size of input buffer for strings
Public gsTempBuf As String 'input buffer for strings (defined at program start)

Public giLatestOptionsTabSelected As Integer 'keep track of the tab that was last selected by user (goes back to it next time it is opened)
Public gsLastStartDateChosen As String 'a starting plot date last selected by user
Public gsLastEndDateChosen As String 'an ending plot date last selected by user
Public gsLastDateSet As String 'a temp string used to pass dates back and forth to the calendar
Public gsDateDisplayFormat As String 'holds the user's choice for the displayed date format for dialogs and reports
Public gsTimeDisplayFormat As String 'holds the user's choice for the displayed time format for dialogs and reports
Public gsCustomLblPatientLastName As String 'replacement labels for the dialogs if exist in config file.
Public gsCustomLblPatientFirstName As String 'replacement labels for the dialogs if exist in config file.
Public gsCustomLblPatientID As String
Public gsCustomLblTxCenter As String
Public gsCustomLblDrug As String
Public gsCustomLblOrgan As String
Public gsLabelGridColumnCustom1 As String
Public gsLabelGridColumnCustom2 As String
Public gsLabelGridColumnCustom3 As String

This value is stored in device & indicates the version of data structure within the device.
This does not relate directly to the version of the host software because the host software version can change with meaning that the structure of the data in the device has changed.
This value should be increased when any kind of change occurs to the custom areas of the device such as changing the length of strings to accommodate new features. The purpose of this value is to let us read it back from a device and determine if newer host software is being used on a device programmed with another version.

Public Const gsREV_DATA_STRUCTURE = "01"

There are 4 fields in the device containing 16 characters each. In the original device design, this was intended to contain 4 separate pieces of information.
The length of each data type is as follows:

Public Const giLEN_REV_DATA_STRUCTURE = 2
Public Const giLEN_PATIENT_NAME = 26
Public Const giLEN_ID = 11
Public Const giLEN_DRUG = 2
Public Const giLEN_TX_CENTER = 18
Public Const giLEN_ORGAN = 2

Public Const giMaxDoseTimes = 4 'the max number of prescribed dosing time (entry boxes)
Public Const giDosesPerDayDefault = 2
Public gbPatientDataNotSaved As Boolean 'true once the data in memory has been saved (from device)

Public gdTempDateTime As Double 'get, put this var in the form that uses it, can also be done
Public giTempCya As Integer
Public giTempCreatinine As Single
Public gsTempCustomInfo As String

Public gsActiveFormName As String


```

Public giCurrentTip As Integer      'most recent tip number that was shown
Public gsWebStartingAddress As String 'url address and any associated password for web site

```

Public Function ComputeIniSectionChecksum(ByVal sFileSpec As String, ByVal sSection As String)

'Read each line in the section name of an INI file that was passed here.

'Compute a unique value and pass back to caller

On Error GoTo 0 'ign temp

```

Dim iChecksumTally As Long, r As Integer, l As Long, l As Integer, iKey As Integer
Dim sLine As String

```

'Get the names of all of the keys in this section.

'A null key field in above line loads all keys in that section

```

Dim iStrSize As Integer, sTempBuf As String, iBufSize As Integer

```

```

Dim sKeyList(2000) As String 'make room for this many key names in this section

```

```

sTempBuf = Space$(16384)

```

```

iBufSize = 16384

```

```

l = GetPrivateProfileString(sSection, ByVal 0&, "", sTempBuf, iBufSize, sFileSpec)

```

```

r = ParseDelimString(Left$(sTempBuf, l), Chr$(0), sKeyList()) 'put the key names in a list

```

```

For iKey = 1 To r

```

```

    sLine = GetINISetting(sFileSpec, sSection, sKeyList(iKey), "")

```

```

    For l = 1 To Len(sLine)

```

```

        iChecksumTally = iChecksumTally + (Asc(Mid$(sLine, l, 1)) * iKey)

```

```

    Next l

```

```

Next iKey

```

```

iChecksumTally = iChecksumTally Mod 536870912 'a 29 bit number

```

```

ComputeIniSectionChecksum = iChecksumTally 'pass result back to caller

```

End Function

Public Sub EventDelete(DataStruct As DeviceDataStruct, ByVal iIndex As Integer)

*'Remove an event from the data structure. The index to the position is
'passed here.*

```

Dim l As Integer

```

'It is not a valid index

```

If iIndex < 1 Or iIndex > DataStruct.lEventData(0) Then Exit Sub

```

```

For l = iIndex To DataStruct.lEventData(0) 'move all events up one

```

```

    DataStruct.byteEventType(l) = DataStruct.byteEventType(l + 1)

```

```

    DataStruct.dEventData(l) = DataStruct.dEventData(l + 1)

```

```

    DataStruct.lEventData(l) = DataStruct.lEventData(l + 1)

```

```

Next l

```

```

DataStruct.lEventData(0) = DataStruct.lEventData(0) - 1 'decrement event count

```

```

gbPatientDataNotSaved = True 'set flag to indicate that the file has changed but not yet been saved

```

End Sub

General.bas - EventInsert

3

Public Sub EventInsert(DataStruct As DeviceDataStruct, ByVal iIndex As Integer, ByVal dDate As Double)

'Insert a new event into the data structure at the index location passed here. If the index = 0 then it probably indicates that a previous function could not find where to insert the date in the structure. In this case, the event must be inserted at the beginning or the end of the structure depending on the date.

Dim I As Integer

```
If iIndex = 0 Then 'date was not found
  If dDate <= DataStruct.dEventData(1) Then
    iIndex = 1
  Else
    iIndex = DataStruct.iEventData(0) 'insert at last point
  End If
End If
```

```
If iIndex Then 'there are events in the structure
  For I = DataStruct.iEventData(0) To iIndex Step -1 'move all events down to make room for new one
    DataStruct.byteEventType(I + 1) = DataStruct.byteEventType(I)
    DataStruct.dEventData(I + 1) = DataStruct.dEventData(I)
    DataStruct.iEventData(I + 1) = DataStruct.iEventData(I)
    DataStruct.sUserData1(I + 1) = DataStruct.sUserData1(I)
    DataStruct.sUserData2(I + 1) = DataStruct.sUserData2(I)
    DataStruct.sUserData3(I + 1) = DataStruct.sUserData3(I)
  Next I
Else
  iIndex = iIndex + 1
End If
```

'Now insert the new event

```
DataStruct.iEventData(0) = DataStruct.iEventData(0) + 1 'increment event count
DataStruct.byteEventType(iIndex) = giEVENT_USER_DEFINED
DataStruct.dEventData(iIndex) = dDate
DataStruct.sUserData1(iIndex) = gTempCys 'put change in structure
DataStruct.sUserData2(iIndex) = gTempCreatinine 'put change in structure
DataStruct.sUserData3(iIndex) = gTempCustomInfo 'put change in structure
```

```
If iIndex = 1 Then 'there were no previous events until this one
  DataStruct.iEventData(iIndex) = 0
Else
  DataStruct.iEventData(iIndex) = DataStruct.iEventData(iIndex - 1)
End If
gbPatientDataNotSaved = True 'set flag to indicate that the file has changed but not yet been saved
```

End Sub

Public Function FindPrescribedDoseSizeForSpecificDay(DataStruct As DeviceDataStruct, ByVal lDate As Long)

*'Find the prescribed dose for the day that is passed here.
'This is accomplished by looking for the most recent dose change
'event that occurred on or prior to this date.*

```
Dim I As Integer, iIndex As Integer
iIndex = FindClosestDateInArray(DataStruct, lDate)
If iIndex = 0 Then 'all events are occurring after the date requested
  For I = 1 To DataStruct.iEventData(0) 'look through whole array if necessary
    If DataStruct.byteEventType(I) = giEVENT_DOSE_CHANGED Then
      FindPrescribedDoseSizeForSpecificDay = I 'DataStruct.iEventData(I)
    End If
  Next I
Else
```

```
'an event date was found
For I = iIndex To 1 Step -1
  If DataStruct.byteEventType(I) = giEVENT_DOSE_CHANGED Then
```

General.bas - FindPrescribedDoseSizeForSpec Day

4

```

    FindPrescribedDoseSizeForSpecificDay = i 'DataStruct.iEventData(i)
    Exit For
End If

Next i
End If

End Function

```

Public Function CalcDayDoseScore_OnTime(DataStruct As DeviceDataStruct, ByVal IStartingDate As Long) As

*'Compute the dosing score for the day passed here.
 'This score tests to see if the doses taken was within the prescribed time range.
 'Pass the score back to the caller as nearest whole percent.
 'Index is the index in the array where computation is to start.
 'It should already be set to the first event that occurred on that day.*

**Dim i As Long, I As Integer, iTotalDoses As Integer
 Dim iIndex As Integer, r As Integer**

```

iIndex = FindClosestDateInArray(DataStruct, IStartingDate) 'returns 0 if date is not found
If iIndex Then 'an event was found
    Do 'look at all past events for the past iScoreDays
        If Int(DataStruct.dEventDate(iIndex)) = IStartingDate Then
            'date still in range, ok to continue
            If DataStruct.byteEventType(iIndex) = giEVENT_DOSE_TAKEN Then 'this is a medication
                'Now test to see if time is within the daily prescribed range
                r = IsDoseWithinPrescribedTimeRange(DataStruct, iIndex) 'pass index to event time
                If r Then iTotalDoses = iTotalDoses + 1
            End If
            iIndex = iIndex + 1
        Else
            Exit Do
        End If
    Loop

```

CalcDayDoseScore_OnTime = 100 * iTotalDoses / DataStruct.iDosesPerDay

End If

End Function

Public Function CalcDayDoseScore_AllDoses(DataStruct As DeviceDataStruct, ByVal IStartingDate As Long) As

*'Compute the dosing score for the day passed here.
 'Calculate for all doses taken on that day regardless of if they were taken on time or not.
 'Pass the score back to the caller as nearest whole percent.
 'Index is the index in the array where computation is to start.
 'It should already be set to the first event that occurred on that day.*

Dim iTotalDoses As Integer, iIndex As Integer

```

iIndex = FindClosestDateInArray(DataStruct, IStartingDate)
If iIndex Then 'an event was found on this date
    Do 'look at all dosing events for this day
        If Int(DataStruct.dEventDate(iIndex)) = IStartingDate Then
            'date still in range, ok to continue
            'this is a medication
            If DataStruct.byteEventType(iIndex) = giEVENT_DOSE_TAKEN Then iTotalDoses = iTotalDoses + 1
            iIndex = iIndex + 1
        Else
            Exit Do
        End If
    Loop

```

CalcDayDoseScore_AllDoses = 100 * iTotalDoses / (DataStruct.iDosesPerDay)

End If

General.bas - CalcDayDoseScore_AltDose

5

End Function

Public Function CalcDosesSumTakenOnSpecificDay(DataStruct As DeviceDataStruct, ByVal IStartingDate As I

*'Compute the dosing total number of doses taken on a specific date
'Note, this calculation does not take into consideration whether or not the dose
'was taken within the prescribed time. This is all doses for a particular day
'Pass the count back to the caller.*

Dim ITodayDoseCount As Integer, IIndex As Integer

IIndex = FindFirstMatchingDateInArray(DataStruct, IStartingDate)

If IIndex Then *'an event was found on this date*

Do *'Look at all dosing events for this day*

If Int(DataStruct.dEventDate(IIndex)) = IStartingDate Then

'date still in range, ok to continue

'this is a medication

If DataStruct.byteEventType(IIndex) = giEVENT_DOSE_TAKEN Then ITodayDoseCount = ITodayDoseCount + 1

IIndex = IIndex + 1 *'goto next higher event in array*

'exit if at end of array (prevents error)

If UBound(DataStruct.dEventDate()) = IIndex Then Exit Do

'exit if no data in array

If IIndex > Int(DataStruct.iEventData(0)) Then Exit Do

Else

Exit Do

End If

Loop

CalcDosesSumTakenOnSpecificDay = ITodayDoseCount

End If

End Function

Public Sub EraseDataInMemory(DataStruct As DeviceDataStruct)

Dim I As Integer

'clear out any data that may be in memory and initialize the arrays

DataStruct.sPatientLastName = ""

DataStruct.sPatientFirstName = ""

DataStruct.sPatientID = ""

DataStruct.sDrug = ""

DataStruct.sOrgan = ""

DataStruct.sTxCenter = ""

DataStruct.sSerialNumber = ""

DataStruct.sFirmwareVer = ""

DataStruct.sDoseSize = ""

DataStruct.sPatientDataFileName = ""

For i = 0 To giMaxDoseTimes

DataStruct.dPrescribedDoseTime(i) = -1

Next i

DataStruct.iDosesPerDay = 0

DataStruct.sDoseResolution = ""

DataStruct.sMedRemaining = ""

Erase DataStruct.sScoreData

Erase DataStruct.iEventData

Erase DataStruct.dEventData

Erase DataStruct.byteEventType *'erases all elements of a fixed array*

Erase DataStruct.sUserData1

Erase DataStruct.sUserData2

General.bas - EraseDataInMemory

6

```

Erase DataStruct.sUserData3

DataStruct.JDeviceInitDate = 0
DataStruct.sBatteryChangeTimer = ""
DataStruct.sDoseLockoutHours = ""

DataStruct.bErrorFatal = False
DataStruct.bErrorNonFatal = False
DataStruct.bErrorDoseSize = False
DataStruct.bErrorMedRemaining = False
DataStruct.bErrorMemoryFull = False
DataStruct.bErrorsExist = False
DataStruct.bErrorBrownOut = False
DataStruct.dLastDownloadDate = 0

gbPatientDataNotSaved = False

End Sub

Public Sub CreateTxtSummaryFile()
    'This routine creates a temp text file in the "fax" subdirectory
    'This will allow the information to be faxed as a text document.

    Dim l As Integer, r As Integer, sFileSpec As String, lErrorCode As Long
    Dim sLblName As String, sLblID As String, sLbITxCenter As String, sLblDrug As String, sLblOrgan As String

    'Get rid of the previous temporary file.
    sFileSpec = App.Path + "Vaxes\temp.txt"
    sFileSpec = App.Path + "Vaxes\" + PAT_DATA.sPatientLastName + ", " + PAT_DATA.sPatientFirstName + " " + PAT_DATA.sPatientID + ".txt"
    r = FileExists(sFileSpec, lErrorCode)
    If r Then Kill sFileSpec

    sLblName = gsCustomLblPatientLastName
    sLblID = gsCustomLblPatientID
    sLbITxCenter = gsCustomLblTxCenter
    sLblDrug = gsCustomLblDrug
    sLblOrgan = gsCustomLblOrgan

    Open sFileSpec For Output Shared As #1
    Print #1, sLblName + " " + PAT_DATA.sPatientLastName + " " + PAT_DATA.sPatientFirstName
    Print #1, sLblID + " " + PAT_DATA.sPatientID
    Print #1, sLbITxCenter + " " + PAT_DATA.sTxCenter
    Print #1, sLblDrug + " " + PAT_DATA.sDrug
    Print #1, sLblOrgan + " " + PAT_DATA.sOrgan
    Print #1,
    Print #1, "Device Serial Number: " + PAT_DATA.sSerialNumber
    Print #1, "FirmWare Version: " + PAT_DATA.sFirmwareVer
    Print #1, "Last Download Date: " + Format$(PAT_DATA.dLastDownloadDate, gsDateDisplayFormat)
    Close #1

    'gh: ensure that the complete file is printed
End Sub

```

General.bas - FileExists

7

Function FileExists(ByVal sPath As String, IErrorCode As Long) As Integer

** Check for existence of a file by attempting an OPEN.
 * Return true (-1) if exists else return False (0) or error condition
 * Note that since this function tries to open a file, an error could
 * return to caller if file is there but in use by another application.*

```
Dim X As Integer
```

```
X = FreeFile
```

```
On Error Resume Next
```

```
Open sPath For Input As X
```

```
Close X
```

```
If Err = 0 Then
```

```
FileExists = True
```

```
IErrorCode = 0 'clear error code
```

```
Else
```

```
FileExists = False 'set flag for error
```

```
IErrorCode = Err 'pass error back to caller
```

```
End If
```

```
End Function
```

Public Function GetINISetting(sFileSpec As String, sSection As String, sKeyField As String, sDefault As String)

```
Dim IStrSize As Integer, sTempBuf As String, IBufSize As Integer
```

```
sTempBuf = Space$(1024)
```

```
IBufSize = 1024
```

```
IStrSize = GetPrivateProfileString(sSection, sKeyField, sDefault, sTempBuf, IBufSize, sFileSpec)
```

```
If IStrSize Then
```

```
GetINISetting = Trim$(Left$(sTempBuf, IStrSize))
```

```
Else
```

```
GetINISetting = sDefault
```

```
End If
```

```
End Function
```

Public Function GetPatientDataFromDisk(ByVal sFileSpec As String, DataStruct As DeviceDataStruct, IErrorRe

**Get all of the patient data from the file on disk and place into memory.*

**The filename that is passed here must be a valid patient file and verified*

by the calling procedure.

```
Dim sSection As String, I As Integer, sTemp As String, r As Integer
```

```
Dim IFileChecksum As Long, IChecksumTally As Long
```

```
On Error GoTo GetPatientDataFromDisk_Error
```

**Read the file and calculate the checksum.*

```
IFileChecksum = ComputeIniSectionChecksum(sFileSpec, "Device Data")
```

```
ICheckSumTally = GetINISetting(sFileSpec, "General", "Device Data Validation", 0)
```

```
If IFileChecksum <> IChecksumTally Then
```

```
IErrorReturn = ERR_DATA_CHECKSUM
```

```
Exit Function
```

```
End If
```

```
IFileChecksum = ComputeIniSectionChecksum(sFileSpec, "Event Data")
```

```
ICheckSumTally = GetINISetting(sFileSpec, "General", "Event Data Validation", 0)
```

```
If IFileChecksum <> IChecksumTally Then
```

```
IErrorReturn = ERR_DATA_CHECKSUM
```

```
Exit Function
```

```
End If
```

General.bas - GetPatientDataFromDisk

8

```

If FileChecksum = ComputeIniSectionChecksum(sFileSpec, "Device Error Flags")
IChecksumTally = GetINISetting(sFileSpec, "General", "Device Error Flags Validation", 0)
If IFileChecksum <> IChecksumTally Then
    IErrorReturn = ERR_DATA_CHECKSUM
    Exit Function
End If

sSection = "Device Error Flags"
DataStruct.bErrorFatal = CBool(GetINISetting(sFileSpec, sSection, "Fatal", False))
DataStruct.bErrorNonFatal = CBool(GetINISetting(sFileSpec, sSection, "Non Fatal", False))
DataStruct.bErrorDoseSize = CBool(GetINISetting(sFileSpec, sSection, "Dose Size", False))
DataStruct.bErrorMedRemaining = CBool(GetINISetting(sFileSpec, sSection, "Med Remaining", False))
DataStruct.bErrorMemoryFull = CBool(GetINISetting(sFileSpec, sSection, "Memory Full", False))
DataStruct.bErrorBrownOut = CBool(GetINISetting(sFileSpec, sSection, "Brownout", False))

sSection = "Device Data"
EraseDataInMemory DataStruct

sTemp = GetINISetting(sFileSpec, sSection, "Device Init Date", 0)
If IsDate(sTemp) Then
    DataStruct.iDeviceInitDate = DateValue(sTemp)
End If

sTemp = GetINISetting(sFileSpec, sSection, "Events Ref Date Time", 0)
If IsDate(sTemp) Then
    DataStruct.dDeviceRefDate = DateValue(sTemp)
End If

sTemp = GetINISetting(sFileSpec, sSection, "Last Download Date", 0)
If IsDate(sTemp) Then
    DataStruct.dLastDownloadDate = DateValue(sTemp)
End If

DataStruct.sPatientLastName = GetINISetting(sFileSpec, sSection, "Last Name", "")
DataStruct.sPatientFirstName = GetINISetting(sFileSpec, sSection, "First Name", "")
DataStruct.sPatientID = GetINISetting(sFileSpec, sSection, "Patient ID", "")
DataStruct.sTxCenter = GetINISetting(sFileSpec, sSection, "Tx Center", "")
i = CInt(GetINISetting(sFileSpec, sSection, "Organ Reference Number", 0))
If i And i <= UBound(gsOrganNames) Then DataStruct.sOrgan = gsOrganNames(i)
i = CInt(GetINISetting(sFileSpec, sSection, "Drug Reference Number", 0))
If i And i <= UBound(gsDrugNames) Then DataStruct.sDrug = gsDrugNames(i)
DataStruct.sSerialNumber = GetINISetting(sFileSpec, sSection, "Serial Number", "")
DataStruct.sFirmwareVer = GetINISetting(sFileSpec, sSection, "Firmware Version", "")
DataStruct.sDoseSize = GetINISetting(sFileSpec, sSection, "Dose Size", "")
DataStruct.iDosesPerDay = CInt(GetINISetting(sFileSpec, sSection, "Doses Per Day", 0))
DataStruct.sDoseResolution = GetINISetting(sFileSpec, sSection, "Dose Resolution", "")
DataStruct.sMedRemaining = GetINISetting(sFileSpec, sSection, "Medication Remaining", "")
DataStruct.sBatteryChangeTimer = GetINISetting(sFileSpec, sSection, "Battery Change Timer", "")
DataStruct.sDoseLockoutHours = GetINISetting(sFileSpec, sSection, "Lockout Hours Between Doses", "")

For i = 1 To 14
    DataStruct.sScoreData(i) = GetINISetting(sFileSpec, sSection, "Patient Score Data " & CStr(i), "")
Next i

For i = 1 To giMaxDoseTimes
    sTemp = GetINISetting(sFileSpec, sSection, "Prescribed Dose Time " & CStr(i), "-1")
    DataStruct.dPrescribedDoseTime(i) = -1 ' default value
    If IsDate(sTemp) Then DataStruct.dPrescribedDoseTime(i) = CDate(sTemp)
Next i

DataStruct.iEventData(0) = CInt(GetINISetting(sFileSpec, "Event Data", "Event Count", 0))

Dim sTempList(10) As String
For i = 1 To DataStruct.iEventData(0)

```

General.bas - GetPatientDataFromDisk

9

```

sTemp = GetINISetting(sFileSpec, "Event Data", CStr(i), "")
r = ParseDelimString(sTemp, "-", sTempList())
DataStruct.dEventDate(i) = CDate(sTempList(1))
Select Case Trim$(LCase$(sTempList(2)))
    Case "dose taken"
        DataStruct.byteEventType(i) = giEVENT_DOSE_TAKEN
        DataStruct.iEventData(i) = sTempList(3)
    Case "dose change"
        DataStruct.byteEventType(i) = giEVENT_DOSE_CHANGED
        DataStruct.iEventData(i) = sTempList(3)
    Case "custom event"
        DataStruct.byteEventType(i) = giEVENT_USER_DEFINED
        DataStruct.iEventData(i) = sTempList(3)
End Select
DataStruct.sUserData1(i) = sTempList(4)
DataStruct.sUserData2(i) = sTempList(5)
DataStruct.sUserData3(i) = sTempList(6)
Next i

' bErrorFatal As Boolean      True if this flag was set in the returned flags string
' bErrorNonFatal As Boolean   True if this flag was set in the returned flags string
' bErrorDoseSize As Boolean    True if this flag was set in the returned flags string
' bErrorMedRemaining As Boolean True if this flag was set in the returned flags string
' bErrorMemoryFull As Boolean  True if this flag was set in the returned flags string
' bErrorBrownOut As Boolean    True if this flag was set in the returned flags string
' bErrorsExist As Boolean      (1 byte) Bits are set if various errors have occurred and have not
GetPatientDataFromDisk = True 'return success flag to caller

```

GetPatientDataFromDisk_Exit:
Exit Function

GetPatientDataFromDisk_Error:
iErrorReturn = Err
Resume GetPatientDataFromDisk_Exit
End Function

Public Sub GetProgramPreferences()

Load the program and user preferences into the global variables
Dim iStrSize As Integer, i As Integer, sFileSpec As String, r As Integer
Dim sSection As String
sSection = "Preferences"
gsDateDisplayFormat = GetINISetting(gsAppIniFileSpec, sSection, "Date Display Format", "Short Date")
gsTimeDisplayFormat = GetINISetting(gsAppIniFileSpec, sSection, "Time Display Format", "Short Time")
gsngComplianceTimeRange = CSng(GetINISetting(gsAppIniFileSpec, sSection, "Compliance Time Range", "2"))
sSection = "Custom Settings"
Get any custom field labels that may be in the INI file. If none exist the set some defaults here.
gsCustomLblPatientLastName = GetINISetting(gsAppIniFileSpec, sSection, "Last Name Label", "")
If gsCustomLblPatientLastName = "" Then gsCustomLblPatientLastName = "Last Name"
gsCustomLblPatientFirstName = GetINISetting(gsAppIniFileSpec, sSection, "First Name Label", "")
If gsCustomLblPatientFirstName = "" Then gsCustomLblPatientFirstName = "First Name"
gsCustomLblPatientID = GetINISetting(gsAppIniFileSpec, sSection, "Patient ID Label", "")
If gsCustomLblPatientID = "" Then gsCustomLblPatientID = "Patient ID"
gsCustomLblTxCenter = GetINISetting(gsAppIniFileSpec, sSection, "TX Center Label", "")
If gsCustomLblTxCenter = "" Then gsCustomLblTxCenter = "TX Center"
gsCustomLblDrug = GetINISetting(gsAppIniFileSpec, sSection, "Drug Label", "")
If gsCustomLblDrug = "" Then gsCustomLblDrug = "Drug"

General.bas - GetProgramPreferences

10

```

gsCustomLblOrgan = GetINISetting(gsAppIniFileSpec, sSection, "Organ Label", ")
If gsCustomLblOrgan = "" Then gsCustomLblOrgan = "Organ"

gsLabelGridColumnCustom1 = GetINISetting(gsAppIniFileSpec, sSection, "Grid Column 1", ")
If gsLabelGridColumnCustom1 = "" Then gsLabelGridColumnCustom1 = "CYA Level (ng/ml)"

gsLabelGridColumnCustom2 = GetINISetting(gsAppIniFileSpec, sSection, "Grid Column 2", ")
If gsLabelGridColumnCustom2 = "" Then gsLabelGridColumnCustom2 = "Creatinine (mg/dl)"

gsLabelGridColumnCustom3 = GetINISetting(gsAppIniFileSpec, sSection, "Grid Column 3", ")
If gsLabelGridColumnCustom3 = "" Then gsLabelGridColumnCustom3 = "Custom"

'Get the list of most recently used files to the menu
For i = 1 To frmMain.mnuFileMRU.UBound
    frmMain.mnuFileMRU(i).Tag = GetINISetting(gsAppIniFileSpec, "Recent Files", CStr(i), ")
    If frmMain.mnuFileMRU(i).Tag <> "" Then
        frmMain.mnuFileMRU(i).Visible = True
        'strip the filespec away from the tag and put into the caption for display purposes
        r = GetFileNameFromSpec(frmMain.mnuFileMRU(i).Tag, sFileSpec) 'hold the name of the file
        frmMain.mnuFileMRU(i).Caption = sFileSpec
        frmMain.mnuFileBar6.Visible = True
    End If
Next i

'Get last values for the Fax control that was last set by user
sSection = "User Selections"
With FAX_DATA
    .sSenderName = GetINISetting(gsFaxFileSpec, sSection, "Sender Name", ")
    .sSenderCompany = GetINISetting(gsFaxFileSpec, sSection, "Sender Company", ")
    .sSenderVoiceNumber = GetINISetting(gsFaxFileSpec, sSection, "Sender Voice Number", ")
    .sSenderFaxNumber = GetINISetting(gsFaxFileSpec, sSection, "Sender Fax Number", ")
    .sFaxID = GetINISetting(gsFaxFileSpec, sSection, "Fax ID", ")
    .sDialPrefix = GetINISetting(gsFaxFileSpec, sSection, "Dial Prefix", ")
    .iRetries = CInt(GetINISetting(gsFaxFileSpec, sSection, "Retries", "0"))
    .iRetryInterval = CInt(GetINISetting(gsFaxFileSpec, sSection, "Retry Interval", "1"))
    .bFaxResolution = GetINISetting(gsFaxFileSpec, sSection, "Resolution", "0")
End With

'Get the Drug types from file and place in global list
sSection = "Transplant Centers"
TxCenters(0) = GetINISetting(gsAppIniFileSpec, sSection, "Count", "0")
For i = 1 To TxCenters(0)
    TxCenters(i) = GetINISetting(gsAppIniFileSpec, sSection, CStr(i), "0")
Next i

'Get the Drug types from file and place in global list
sSection = "Drugs"
gsDrugNames(0) = GetINISetting(gsAppIniFileSpec, sSection, "Count", "0")
For i = 1 To gsDrugNames(0)
    gsDrugNames(i) = GetINISetting(gsAppIniFileSpec, sSection, CStr(i), "0")
Next i

'Get the Drug types from file and place in global list
sSection = "Organs"
gsOrganNames(0) = GetINISetting(gsAppIniFileSpec, sSection, "Count", "0")
For i = 1 To gsOrganNames(0)
    gsOrganNames(i) = GetINISetting(gsAppIniFileSpec, sSection, CStr(i), "0")
Next i

giCurrentTip = CInt(GetINISetting(gsAppIniFileSpec, "Options", "Current Tip", 1))

'Get settings of calendar form
CAL_DEFAULTS.chkDosesMissed = CByte(GetINISetting(gsAppIniFileSpec, "Calendar Settings", "chkDosesMissed", 1))
CAL_DEFAULTS.chkDosesNotComplied = CByte(GetINISetting(gsAppIniFileSpec, "Calendar Settings", "chkDosesNotComplied", 1))

```

General.bas - GetProgramPreferences

11

```

CAL_DEFAULTS.chkDosesTaken = CByte(GetINISetting(gsAppIniFileSpec, "Calendar Settings", "chkDosesTaken", 0))
CAL_DEFAULTS.chkDoseChanged = CByte(GetINISetting(gsAppIniFileSpec, "Calendar Settings", "chkDoseChanged", 1))

'Get Settings of Patient summary form
PAT_SUM_DEFAULTS.cmboDataToView = CByte(GetINISetting(gsAppIniFileSpec, "Patient Summary Settings", "cmboDataToView", 1))
PAT_SUM_DEFAULTS.cmboChartType = CByte(GetINISetting(gsAppIniFileSpec, "Patient Summary Settings", "cmboChartType", 1))

'Get the web address for the browser
gsWebStartingAddress = "http://"
gsWebStartingAddress = gsWebStartingAddress + GetINISetting(gsAppIniFileSpec, "Web Data", "User Name", "") + ":"
gsWebStartingAddress = gsWebStartingAddress + GetINISetting(gsAppIniFileSpec, "Web Data", "Access", "") + "@"
gsWebStartingAddress = gsWebStartingAddress + GetINISetting(gsAppIniFileSpec, "Web Data", "URL", "")

End Sub

```

Sub Main()

```

Dim l As Long, r As Integer, i As Integer, sMSG As String, sTemp As String, dTime As Double
Dim bBrowserFound As Boolean, iLastAccessDate As Long, iNextReminderDate As Long

'Initialize some application settings
gsAppIniFileSpec = App.Path + "\CycloTech.ini"
gsFaxFileSpec = App.Path + "\Fax.ini"
gsTempBuf = Space$(1024)
gbCommOK = 99 'set to some value other than true or false to properly initialize the dialog

frmSplash.Show
frmSplash.Refresh
dTime = Now 'get the time value of now

Wait 0.75
frmLogin.Show vbModal
If Not frmLogin.OK Then End 'Login Failed so exit app
Unload frmLogin
DoEvents

'rgn note that the debug flag is turned on and the fax icon is hidden
i = Shell(App.Path + "\Faxman32.exe /D /H", 1) 'start the fax server
i = Shell("\Faxman32.exe /H", 1) 'start the fax server
Load frmMain
Set gcFax = frmMain.FaxMan1

'If browser feature is turned on in the ini file, then activate item on the menu.
'See if we should allow access to visit Sangstat on the Internet
r = CBool(GetINISetting(gsAppIniFileSpec, "Web Data", "Active", "False"))
If r = False Then frmMain.mnuAccessWebSite.Visible = False 'no key found in ini file

Comm_InitializeCommPort 'initialize the comm port from INI file settings
GetProgramPreferences
EraseDataInMemory PAT_DATA
EraseDataInMemory TEMP_DATA

'Set some menu items
frmMain.mnuFileSave.Enabled = False
frmSplash.ZOrder

Do
    If CDBl(Now) > dTime + 0.00005 Then Exit Do 'wait for a minimum amount of time before unloading splash screen
    DoEvents
Loop

```

General.bas - Main

12

Unload frmSplash

frmMain.Show

SetPrinterIcon False, "

*'See if we should shown tips at startup*r = CBool(GetINISetting(gsAppIniFileSpec, "Options", "Show Tips at Startup", True))
If r Then frmTip.Show*'See if we should remind user to visit Sangsat on the Internet*If frmMain.mnuAccessWebSite.Visible = True Then *'user must have menu selection on to access web**'Try to find a browser by looking in different locations in the registry*

frmMain.MhIni1.Key = ClassesRoot

frmMain.MhIni1.EntrySection = "HTTP\shell\open\command"

frmMain.MhIni1.EntryItem = "" *'gets the default value*frmMain.MhIni1.Action = 13 *'get registry key*

If Len(Trim(frmMain.MhIni1.EntryValue)) > 0 Then

bBrowserFound = True *'Looks like a value is there*

End If

frmMain.MhIni1.Key = LocalMachine

frmMain.MhIni1.EntrySection = "SOFTWARE\Classes\HTTP\shell\open\command"

frmMain.MhIni1.EntryItem = "" *'gets the default value*frmMain.MhIni1.Action = 13 *'get registry key*

If Len(Trim(frmMain.MhIni1.EntryValue)) > 0 Then

bBrowserFound = True *'Looks like a value is there*

End If

If bBrowserFound Then

iLastAccessDate = 0

sTemp = GetINISetting(gsAppIniFileSpec, "Web Data", "Last Web Visit Date", "")

If IsDate(sTemp) Then iLastAccessDate = DateValue(sTemp)

iNextReminderDate = 0

sTemp = GetINISetting(gsAppIniFileSpec, "Web Data", "Next Web Visit Reminder Date", "")

If IsDate(sTemp) Then iNextReminderDate = DateValue(sTemp)

*'On this line, if iLastAccessDate is negative then it indicates that the user chose to not connect**'the last time he/she was reminded. In this case we wait a much**'shorter period of time before reminding them again.*If DateValue(Now) >= iNextReminderDate Then *'it's been too long since the user was last on the web.*

If iLastAccessDate > 0 Then

sMSG = "You last connected to our internet web site on " & sTemp & ". "

Else

sMSG = "You have not yet connected to our internet web site. "

End If

sMSG = sMSG & "There may be a program update or other valuable information there."

sMSG = sMSG & vbCrLf & vbCrLf & "Would you like to connect to the web site now? (you must already have web access available

)"

Beep

r = MsgBox(sMSG, vbQuestion + vbYesNo + vbDefaultButton2, "Internet Connection Reminder")

*'regardless of the answer to the next question, set a minimum time to ask user again.**'If user actually connects to internet, then this time is overwritten with a longer one by the browser.*

SaveINISetting gsAppIniFileSpec, "Web Data", "Next Web Visit Reminder Date", Format\$(Now + 15, "Medium Date")

If r = vbYes Then

Call LogonToWebSite

End If

End If

End If

End If

End Sub

General.bas - LogonToWebSite

13

Public Sub LogonToWebSite()

```

'Visit Sangstat on the Internet
Dim frmB As New frmBrowser
Load frmBrowser
frmBrowser.StartingAddress = gsWebStartingAddress
DoEvents 'allow time to paint
DoEvents 'allow time to paint
frmBrowser.Refresh
frmBrowser.Show
End Sub

```

Public Function OpenPatientData(ByVal sFileSpec As String) As Integer

```

'Open patient data file and load to memory
'Return a true if load was successful, false if not, and vbCancel if user cancelled

Dim r As Integer, sTemp As String, IErrorCode As Long
On Error GoTo OpenPatientData_Error

r = ValidatePatientDataSaved 'make sure any device data has first been saved
If r = vbCancel Then Exit Function

'Get a filename from the common dialog
'Setup the common dialog control prior to showing it
With frmMain.dlgCommonDialog
.Flags = cdlOFNOverwritePrompt Or cdlOFNPathMustExist Or cdlOFNExplorer Or cdlOFNExtensionDifferent Or
.cdlOFNReadOnlyReturn Or cdlOFNHideReadOnly Or cdlOFNFileMustExist
.CancelError = True 'generate error if CANCEL button is pressed
.InitialDir = App.Path & "\Patient Data"
.Filter = "CycloTech Data File *.*cpd"
.DialogTitle = "Open Patient Data File"
.DefaultExt = "CPD" 'append "Structure" extension when saving.
If sFileSpec <> "" Then .filename = sFileSpec
.ShowOpen 'Open dialog
End With

'Now get the data from file
frmMain.MousePointer = vbHourglass
DoEvents

r = GetPatientDataFromDisk(frmMain.dlgCommonDialog.filename, PAT_DATA, IErrorCode)
If r <> True Then
    If IErrorCode = ERR_DATA_CHECKSUM Then
        sTemp = "The contents of the data file have changed since it was last saved. "
        sTemp = sTemp & "This could be due to a corrupt file, but is more likely that the file was manually changed."
        sTemp = sTemp & vbCrLf & vbCrLf & "The file will not be loaded."
        Beep
        MsgBox sTemp, vbCritical, "File Contents Changed"
    Else
        MsgBox "An error occurred while retrieving data from the file. It was not read.", vbExclamation, "Error In File - " & Error(r)
    End If
End If

PAT_DATA.sPatientDataFileName = frmMain.dlgCommonDialog.filename
OpenPatientData = True
r = GetFileNameFromSpec(PAT_DATA.sPatientDataFileName, sTemp) 'hold the name of the file
UpdateRecentFileMenu sFileSpec
UpdateRecentFileMenu sTemp
frmMain.mnuFileSave.Enabled = True

OpenPatientData_Exit:
On Error GoTo 0
RefreshAllOpenForms
frmMain.MousePointer = vbDefault
Exit Function

```

General.bas - OpenPatientData

14

```

OpenPatientData_Error:
  If Err = cdiCancel Then      'cancel button was pressed in dialog
    OpenPatientData = vbCancel
    Resume OpenPatientData_Exit
  Else
    Beep
    MsgBox "The CycloTech Data file contains invalid data and can not be read.", vbExclamation, "Invalid Data File - " & Error
    PAT_DATA.sPatientDataFileName = ""
    frmMain.mnuFileSave.Enabled = False
  End If
  OpenPatientData = False
  Resume OpenPatientData_Exit  'exit anyway for now

```

End Function

Public Sub PopulateDeviceDiagDialog(DataStruct As DeviceDataStruct, SourceForm As Form)

'There are two possible dialogs that have the same controls on them.

'This common procedure will populate both

Dim I As Integer

With SourceForm

'Show custom labels from config file if there were any

.Label1(3) = gsCustomLblPatientLastName

.Label1(1) = gsCustomLblPatientFirstName

.Label1(5) = gsCustomLblPatientID

.Label1(6) = gsCustomLblTxCenter

.Label1(7) = gsCustomLblDrug

.Label1(0) = gsCustomLblOrgan

If TypeOf .btDrug Is SSPanel Then

.btDrug = DataStruct.sDrug

ElseIf TypeOf .btDrug Is ComboBox Then 'It is a list box

.btDrug.Clear

For i = 1 To gsDrugNames(0) 'fill the drugs list box with available choices

.btDrug.AddItem gsDrugNames(i)

Next i

For i = 0 To .btDrug.ListCount - 1

If .btDrug.List(i) = DataStruct.sDrug Then

.btDrug.ListIndex = i

Exit For

End If

Next i

End If

If TypeOf .btOrgan Is SSPanel Then

.btOrgan = DataStruct.sOrgan

ElseIf TypeOf .btOrgan Is ComboBox Then 'It is a list box

.btOrgan.Clear

For i = 1 To gsOrganNames(0) 'fill the drugs list box with available choices

.btOrgan.AddItem gsOrganNames(i)

Next i

For i = 0 To .btOrgan.ListCount - 1

If .btOrgan.List(i) = DataStruct.sOrgan Then

.btOrgan.ListIndex = i

Exit For

End If

Next i

End If

.btPatientLastName = DataStruct.sPatientLastName

General.bas - PopulateDeviceDiagDialog

15

```

        .btPatientFirstName = DataStruct.sPatientFirstName
        .btPatientID = DataStruct.sPatientID
        .btTxCenter = DataStruct.sTxCenter
        .btSerialNumber = DataStruct.sSerialNumber
        .btDoseSize = DataStruct.sDoseSize

        .btPatientLastName.SetFocus

        If DataStruct.iEventData(0) Then
            .btEventCount = "" + CStr(DataStruct.iEventData(0)) + ""
        Else
            .btEventCount = ""
        End If

        For i = 1 To giMaxDoseTimes
            If DataStruct.dPrescribedDoseTime(i) >= 0 Then .txtDoseTime(i) = Format$(DataStruct.dPrescribedDoseTime(i),
                gsTimeDisplayFormat)
        Next i

        .bdDosesPerDay = CStr(DataStruct.iDosesPerDay)
        .btDoseResolution = DataStruct.sDoseResolution
        .btDoseLockoutHours = DataStruct.sDoseLockoutHours
        If DataStruct.iDeviceInitDate Then
            .btDeviceStarted = "" + Format$(CDate(DataStruct.iDeviceInitDate), "Medium Date")
        End If

        .btMedicationRemaining = "" + DataStruct.sMedRemaining
        .btBatteryChangeTimer = "" + DataStruct.sBatteryChangeTimer
        .btFirmwareVer = "" + DataStruct.sFirmwareVer

        'set indicators for error flags
        If DataStruct.bErrorFatal Then
            .imgFatalPicture = .imgError.Picture
        Else
            .imgFatalPicture = .imgNoError.Picture
        End If

        If DataStruct.bErrorNonFatal Then
            .imgNonFatalPicture = .imgError.Picture
        Else
            .imgNonFatalPicture = .imgNoError.Picture
        End If

        If DataStruct.bErrorDoseSize Then
            .imgDoseSizePicture = .imgError.Picture
        Else
            .imgDoseSizePicture = .imgNoError.Picture
        End If

        If DataStruct.bErrorMedRemaining Then
            .imgMedRemainingPicture = .imgError.Picture
        Else
            .imgMedRemainingPicture = .imgNoError.Picture
        End If

        If DataStruct.bErrorMemoryFull Then
            .imgMemoryFullPicture = .imgError.Picture
        Else
            .imgMemoryFullPicture = .imgNoError.Picture
        End If

        If DataStruct.bErrorBrownOut Then
            .imgBrownOutPicture = .imgError.Picture
        Else
            .imgBrownOutPicture = .imgNoError.Picture
        End If

        End With

```

General.bas - PopulateDeviceDiagDialog

16

End Sub

Public Sub RefreshAllOpenForms()

Dim r As Integer

*If any of these forms are open at the time a new file is loaded,
then refresh them.*

For r = 0 To Forms.Count - 1

Select Case Forms(r).Name

Case "frmPatientDosingReport"

frmPatientDosingReport.UpdatefrmPatientDosingReportHeader

frmPatientDosingReport.UpdatePatientGridDisplay

Case "frmDosingCalendar"

*If PAT_DATA.dEventData(PAT_DATA.iEventData(0)) > 0 Then frmDosingCalendar.Calendar.Date = CDate(PAT_DATA.
dEventData(PAT_DATA.iEventData(0)))*

UpdateCalendar

Case "frmPrint"

RefreshPreview

Case "frmPatientSummary"

frmPatientSummary.UpdatefrmPatientSummaryHeader

frmPatientSummary.cmboDateSelection_Click

frmPatientSummary.UpdatePatientDosingGraph

Case "frmDeviceInitialize"

PopulateDeviceCommDialog PAT_DATA, frmDeviceInitialize

Case "frmReadDeviceData"

PopulateDeviceCommDialog PAT_DATA, frmReadDeviceData

End Select

Next r

End Sub

Public Sub SetPrinterIcon(bEnable As Boolean, sTip As String)

On Error Resume Next

frmMain.mnuFilePrint.Enabled = bEnable

If sTip = "" Then

frmMain.mnuFilePrint.Caption = "Print..."

Else

frmMain.mnuFilePrint.Caption = sTip

End If

frmMain.tbToolBar.Buttons.Item(5).Enabled = bEnable

frmMain.tbToolBar.Buttons.Item(5).ToolTipText = sTip

*If the active form is not the print form, then keep the name of the
form in the key property of the icon. This is so that the print
form will know what kind of information to display and print.*

If frmMain.ActiveForm.Name <> "frmPrint" Then gsActiveFormName = frmMain.ActiveForm.Name

On Error GoTo 0

End Sub

General.bas - UpDateRecentFileMenu

17

Public Sub UpDateRecentFileMenu(ByVal sFileSpec As String)

'Add the newest FileName to the menuist and move the other ones down.
On Error GoTo 0

Dim bDuplicateFound As Boolean, i As Integer, r As Integer, sFileName As String
r = GetFileNameFromSpec(sFileSpec, sFileName) 'hold the name of the file

With frmMain

For i = 1 To .mnuFileMRU.UBound - 1

If LCase\$(.mnuFileMRU(i).Caption) = LCase\$(sFileName) Then 'remove any duplicates that might appear
.mnuFileMRU(i).Caption = "
.mnuFileMRU(i).Tag = "
bDuplicateFound = True
End If
Next i

For i = .mnuFileMRU.UBound - 1 To 1 Step -1

If .mnuFileMRU(i).Caption <> "" Then 'contains a filename ok to shift down
.mnuFileMRU(i + 1).Caption = .mnuFileMRU(i).Caption 'holds filename only for display purposes
.mnuFileMRU(i + 1).Tag = .mnuFileMRU(i).Tag 'holds the filespec
.mnuFileMRU(i + 1).Visible = True
Else
.mnuFileMRU(i + 1).Visible = False
End If
Next i

.mnuFileMRU(1).Tag = sFileSpec
.mnuFileMRU(1).Caption = sFileName
.mnuFileMRU(1).Visible = True
.mnuFileBar6.Visible = True

End With

End Sub

Public Function GetFileNameFromSpec(ByVal sFileSpec As String, sFileName As String) As Integer

'Strip the filename and extension from the filespec (drive\path\filename)

Dim r As Integer

ReDim sList(50) As String

On Error GoTo GetFileNameFromSpec_Error

If Len(sFileSpec) > 0 Then
r = ParseDelimString(sFileSpec, "\", sList()) 'delimit all subpaths
sFileName = LCase\$(sList(r)) 'the name is last item in list

'something was returned

If Len(sFileName) > 0 Then GetFileNameFromSpec = True 'return success to caller

End If

GetFileNameFromSpec_Exit:
Exit Function

GetFileNameFromSpec_Error:
'Resume 0
Resume GetFileNameFromSpec_Exit

End Function

General.bas - ParseDelimString

18

Public Function ParseDelimString(ByVal sParse As String, ByVal sDelim As String, sFieldStrings() As String) As Integer

'Parse "sParse" passed here. Put resulting parsed names in a list called sFieldStrings.

'Use sDelim as the delimiter to parse string.

'Trim any leading and trailing spaces from each field.

'sFieldString list must pre-exist before calling here, and should be big enough to hold all delimited strings.

'The list contains fields in the order they appeared from left to right.

'Function returns number of fields found.

Dim i As Integer 'loop counter
Dim iDelim1 As Integer, iDelim2 As Integer 'marks beginning and end of a field

If Len(sParse) = 0 Then Exit Function 'exit if no chars in string

iDelim1 = 0 'set first delim marker to beginning of line

If Right\$(sParse, Len(sDelim)) <> sDelim Then 'see if a delim is already at the end of the string

sParse = sParse + sDelim 'put a delim at end of line

End If

'Note: an Erase method can not be used as it reduces the array to only a few elements

For i = 0 To UBound(sFieldStrings) 'clear out old data from the array

sFieldStrings(i) = ""

Next i

i = 0

Do While iDelim1 < Len(sParse) 'keep looking til all delims are found

iDelim2 = InStr(iDelim1 + 1, sParse, sDelim) 'look for delim in string

'get field from string, trim off spaces and put field into list

sFieldStrings(i + 1) = Trim\$(Mid\$(sParse, iDelim1 + 1, iDelim2 - iDelim1 - 1))

iDelim1 = iDelim2 'reset first delim marker to latest one found

i = i + 1 'increment field counter

Loop 'repeat search

ParseDelimString = i 'put parsed items count in element 0 of list

End Function

Public Function SaveDataToNewFile() As Integer

'Get a filename from the common dialog

'Setup the common dialog control prior to showing it

On Error GoTo SaveDataToNewFile_Error

With frmMain.dlgCommonDialog

.Flags = cdiOFNOverwritePrompt Or cdiOFNCreatePrompt Or cdiOFNPathMustExist Or cdiOFNExplorer Or cdiOFNExtensionDifferent

Or cdiOFNNoReadOnlyReturn Or cdiOFNHideReadOnly

.CancelError = True 'generate error if CANCEL button is pressed

.InitDir = App.Path + "\Patient Data"

*.Filter = "CydoTech Data File *.cpd|.cpd"*

.DialogTitle = "Save Patient Data As..."

.DefaultExt = ".CPD" 'append "Structure" extension when saving.

If PAT_DATA.sPatientDataFileName = "" Then

.filename = PAT_DATA.sPatientLastName + " " + PAT_DATA.sPatientFirstName + " " + PAT_DATA.sPatientID + ".cpd" 'set a

default file name

Else

.filename = PAT_DATA.sPatientDataFileName

End If

.ShowSave 'save as dialog

End With

PAT_DATA.sPatientDataFileName = frmMain.dlgCommonDialog.filename

SavePatientData PAT_DATA.sPatientDataFileName

SaveDataToNewFile = True

SaveDataToNewFile_Exit:

Exit Function

SaveDataToNewFile_Error:

General.bas - SaveDataToNewFile

19

```

If Err = cdCancel Then      'cancel button was pressed in dialog
    SaveDataToNewFile = vbCancel
    Resume SaveDataToNewFile_Exit
End If
SaveDataToNewFile = False

```

End Function

Public Function SavePatientData(ByVal sFileSpec As String) As Integer

'Save all of the patient data currently in memory to a disk file
'Return a true if save was successful, false if not, and vbCancel if user cancelled

```

Dim sTemp As String, r As Integer, i As Integer, sSection As String
Dim iCheckSumTally As Long

```

On Error GoTo SavePatientData_Error

' r = GetFileNameFromSpec(frmMain.dlgCommonDialog.filename, sTemp) 'save the dir was selected by user

'We need to confirm with user that it is desired to save these file modifications under the same name as the one that was just loaded.

```

If Len(frmMain.dlgCommonDialog.filename) And frmMain.dlgCommonDialog.filename = UCase$(Proj.sStructFileName) Then
    MSGS = "You are about to save changes to the same file they were loaded from!"
    MSGS = MSGS + " Are you sure you want to do this?"
    Beep
    r = MsgBox(MSGS, MB_YES_NO, "Confirm Over Write")
    If r = ID_NO Then Exit Function      'oops, user almost made mistake, exit sub
End If

```

'Now save the data to the file

```

r = GetFileNameFromSpec(sTemp, sFileSpec)      'hold the name of the file
sFileSpec = App.Path & "\Patient Data"

```

```

frmMain.MousePointer = vbHourglass
DoEvents

```

sSection = "Device Data"

SaveINISetting sFileSpec, sSection, "Date Saved To File", Now

SaveINISetting sFileSpec, sSection, "Host Software Version", CStr(App.Major & "." & App.Minor & "." & App.Revision)

SaveINISetting sFileSpec, sSection, "Firmware Version", PAT_DATA.sFirmwareVer

SaveINISetting sFileSpec, sSection, "Last Download Date", CDate(PAT_DATA.dLastDownloadDate) 'short date must be used to prevent error when loading back

SaveINISetting sFileSpec, sSection, "Device Init Date", CDate(PAT_DATA.iDeviceInitDate)

SaveINISetting sFileSpec, sSection, "Events Ref Date Time", CDate(PAT_DATA.dDeviceRefDateTime)

SaveINISetting sFileSpec, sSection, "Last Name", PAT_DATA.sPatientLastName

SaveINISetting sFileSpec, sSection, "First Name", PAT_DATA.sPatientFirstName

SaveINISetting sFileSpec, sSection, "Serial Number", PAT_DATA.sSerialNumber

SaveINISetting sFileSpec, sSection, "Patient ID", PAT_DATA.sPatientID

SaveINISetting sFileSpec, sSection, "Organ", PAT_DATA.sOrgan

SaveINISetting sFileSpec, sSection, "Organ Reference Number", CStr(GetOrganRefNumber())

SaveINISetting sFileSpec, sSection, "Tx Center", PAT_DATA.sTxCenter

SaveINISetting sFileSpec, sSection, "Drug", PAT_DATA.sDrug

SaveINISetting sFileSpec, sSection, "Drug Reference Number", CStr(GetDrugRefNumber())

SaveINISetting sFileSpec, sSection, "Dose Size", PAT_DATA.sDoseSize

SaveINISetting sFileSpec, sSection, "Doses Per Day", CStr(PAT_DATA.iDosesPerDay)

SaveINISetting sFileSpec, sSection, "Dose Resolution", PAT_DATA.sDoseResolution

SaveINISetting sFileSpec, sSection, "Medication Remaining", PAT_DATA.sMedRemaining

SaveINISetting sFileSpec, sSection, "Battery Change Timer", PAT_DATA.sBatteryChangeTimer

SaveINISetting sFileSpec, sSection, "Lockout Hours Between Doses", PAT_DATA.sDoseLockoutHours

For i = 1 To 14

SaveINISetting sFileSpec, sSection, "Patient Score Data " & CStr(i), PAT_DATA.sScoreData(i)

Next i

For i = 1 To giMaxDoseTimes

If PAT_DATA.dPrescribedDoseTime(i) >= 0 Then

General.bas - SavePatientData

20

```

SaveINISetting sFileSpec, sSection, "Prescribed Dose Time " + CStr(i), Format$(PAT_DATA.dPrescribedDoseTime(i),
    gsTimeDisplayFormat)
Else
    SaveINISetting sFileSpec, sSection, "Prescribed Dose Time " + CStr(i), "None"
End If
Next i

```

This section is finished. Go compute the checksum and save it.
 ICheckSumTally = ComputeIniSectionChecksum(sFileSpec, sSection)
 SaveINISetting sFileSpec, "General", "Device Data Validation", CStr(ICheckSumTally)

Before saving new event data, clear out the old strings
 sSection = "Event Data"
 r = WritePrivateProfileString(sSection, ByVal 0&, ByVal 0&, sFileSpec)
 SaveINISetting sFileSpec, sSection, "Event Count", CStr(PAT_DATA.iEventData(0))
 For j = 1 To PAT_DATA.iEventData(0) *'total number of events*
 sTemp = Format\$(PAT_DATA.dEventData(i), "General Date") + ", "
 Select Case PAT_DATA.byteEventType(i)
 Case giEVENT_DOSE_TAKEN
 sTemp = sTemp + "Dose Taken, "
 Case giEVENT_DOSE_CHANGED
 sTemp = sTemp + "Dose Change, "
 Case giEVENT_USER_DEFINED
 sTemp = sTemp + "Custom Event, "
 End Select
 sTemp = sTemp + CStr(PAT_DATA.iEventData(i))
 sTemp = sTemp + " " + PAT_DATA.sUserData1(i)
 sTemp = sTemp + " " + PAT_DATA.sUserData2(i)
 sTemp = sTemp + " " + PAT_DATA.sUserData3(i)
 SaveINISetting sFileSpec, sSection, CStr(i), sTemp
 Next j

This section is finished. Go compute the checksum and save it.
 ICheckSumTally = ComputeIniSectionChecksum(sFileSpec, sSection)
 SaveINISetting sFileSpec, "General", "Event Data Validation", CStr(ICheckSumTally)

sSection = "Device Error Flags"
 SaveINISetting sFileSpec, sSection, "Fatal", CStr(PAT_DATA.bErrorFatal)
 SaveINISetting sFileSpec, sSection, "Non Fatal", CStr(PAT_DATA.bErrorNonFatal)
 SaveINISetting sFileSpec, sSection, "Dose Size", CStr(PAT_DATA.bErrorDoseSize)
 SaveINISetting sFileSpec, sSection, "Med Remaining", CStr(PAT_DATA.bErrorMedRemaining)
 SaveINISetting sFileSpec, sSection, "Memory Full", CStr(PAT_DATA.bErrorMemoryFull)
 SaveINISetting sFileSpec, sSection, "Brownout", CStr(PAT_DATA.bErrorBrownOut)

This section is finished. Go compute the checksum and save it.
 ICheckSumTally = ComputeIniSectionChecksum(sFileSpec, sSection)
 SaveINISetting sFileSpec, "General", "Device Error Flags Validation", CStr(ICheckSumTally)

gbPatientDataNotSaved = False
 r = GetFileNameFromSpec(sFileSpec, PAT_DATA.sPatientDataFileName) *'hold the name of the file*
 UpdateRecentFileMenu sFileSpec
 frmMain.mnuFileSave.Enabled = True
 SavePatientData = True

SavePatientData_Exit:
 On Error GoTo 0
 frmMain.MousePointer = vbDefault
 Exit Function

SavePatientData_Error:
 SavePatientData = False

General.bas - SavePatientData

21

Resume SavePatientData_Exit 'exit anyway for now

End Function

Public Sub PopulateDeviceCommDialog(DataStruct As DeviceDataStruct, SourceForm As Form)

'There are two possible dialogs that have the same controls on them.

'This common procedure will populate both

On Error Resume Next 'not all text boxes will appear on every form

Dim i As Integer

With SourceForm

'Show custom labels from config file if there were any

.Label1(3) = gsCustomLblPatientLastName

.Label1(1) = gsCustomLblPatientFirstName

.Label1(5) = gsCustomLblPatientID

.Label1(6) = gsCustomLblTxCenter

.Label1(7) = gsCustomLblDrug

.Label1(0) = gsCustomLblOrgan

If TypeOf .txtDrug Is SSPanel Then

.txtDrug = DataStruct.sDrug

ElseIf TypeOf .txtDrug Is ComboBox Then 'it is a list box

.txtDrug.Clear

For i = 1 To gsDrugNames(0) 'fill the drugs list box with available choices

.txtDrug.AddItem gsDrugNames(i)

Next i

For i = 0 To .txtDrug.ListCount - 1

If .txtDrug.List(i) = DataStruct.sDrug Then

.txtDrug.ListIndex = i

Exit For

End If

Next i

End If

If TypeOf .txtOrgan Is SSPanel Then

.txtOrgan = DataStruct.sOrgan

ElseIf TypeOf .txtOrgan Is ComboBox Then 'it is a list box

.txtOrgan.Clear

For i = 1 To gsOrganNames(0) 'fill the drugs list box with available choices

.txtOrgan.AddItem gsOrganNames(i)

Next i

For i = 0 To .txtOrgan.ListCount - 1

If .txtOrgan.List(i) = DataStruct.sOrgan Then

.txtOrgan.ListIndex = i

Exit For

End If

Next i

End If

.txtPatientLastName = DataStruct.sPatientLastName

.txtPatientFirstName = DataStruct.sPatientFirstName

.txtPatientID = DataStruct.sPatientID

.txtTxCenter = DataStruct.sTxCenter

.txtSerialNumber = DataStruct.sSerialNumber

.txtDoseSize = "" & DataStruct.sDoseSize

General.bas - PopulateDeviceCommDia

22

```

If DataStruct.iEventData(0) Then
    .bEventCount = "" + CStr(DataStruct.iEventData(0)) + ""
Else
    .bEventCount = ""
End If

If PAT_DATA.dLastDownloadDate Then
    .bLastRetrievalDate = "" + Format$(CDate(DataStruct.dLastDownloadDate), "Short Date") + "" + Format$(CDate(DataStruct.
dLastDownloadDate), "Medium Time") + ""
Else
    .bLastRetrievalDate = ""
End If

.bPatientLastName Take focus away from list box after it was set

For i = 1 To giMaxDoseTimes
    If DataStruct.dPrescribedDoseTime(i) >= 0 Then
        .bDoseTime(i) = "" + Format$(DataStruct.dPrescribedDoseTime(i), gsTimeDisplayFormat) + ""
    End If
Next i

.bDosesPerDay = "" + CStr(DataStruct.iDosesPerDay) + ""
.bDoseResolution = "" + DataStruct.sDoseResolution + ""
.bDoseLockoutHours = "" + DataStruct.sDoseLockoutHours + ""
.bMedicationRemaining = "" + DataStruct.sMedRemaining + ""

If DataStruct.iDeviceInitDate Then
    .bDeviceStarted = "" + Format$(CDate(DataStruct.iDeviceInitDate), "Medium Date")
End If
.bBatteryChangeTimer = "" + DataStruct.sBatteryChangeTimer + ""

'set indicators for error flags
With DataStruct
    If .bErrorFatal Or .bErrorNonFatal Or .bErrorDoseSize Or .bErrorMedRemaining Or .bErrorMemoryFull Or .bErrorBrownOut Then
        SourceForm.imgErrorReceived.Visible = True 'errors were found
        SourceForm.lblErrorsReceived.Visible = True
    Else
        SourceForm.imgErrorsReceived.Visible = False 'no errors exist
        SourceForm.lblErrorsReceived.Visible = False
    End If
End With

End With
On Error GoTo 0

End Sub

```

Public Sub SaveProgramPreferences()

```

Dim i As Integer, sSection As String, sFileSpec As String
sSection = "Preferences"
SaveINISetting gsAppIniFileSpec, sSection, "Date Display Format", gsDateDisplayFormat
SaveINISetting gsAppIniFileSpec, sSection, "Time Display Format", gsTimeDisplayFormat
SaveINISetting gsAppIniFileSpec, sSection, "Compliance Time Range", CStr(gsngComplianceTimeRange)

```

'Save the names of the most recently used files from the menu

```

For i = 1 To frmMain.mnuFileMRU.UBound
    SaveINISetting gsAppIniFileSpec, "Recent Files", CStr(i), frmMain.mnuFileMRU(i).Caption
Next i

```

```

SaveINISetting gsAppIniFileSpec, "Options", "Current Tip", CStr(giCurrentTip)

```

'Save Settings of Calendar Form

```

SaveINISetting gsAppIniFileSpec, "Calendar Settings", "chkDosesMissed", CStr(CAL_DEFAULTS.chkDosesMissed)
SaveINISetting gsAppIniFileSpec, "Calendar Settings", "chkDosesNotComplied", CStr(CAL_DEFAULTS.chkDosesNotComplied)
SaveINISetting gsAppIniFileSpec, "Calendar Settings", "chkDosesTaken", CStr(CAL_DEFAULTS.chkDosesTaken)
SaveINISetting gsAppIniFileSpec, "Calendar Settings", "chkDoseChanged", CStr(CAL_DEFAULTS.chkDoseChanged)

```

General.bas - SaveProgramPreferences

```

'Save Settings of Patient Summary Form
SaveINISetting gsAppIniFileSpec, "Patient Summary Settings", "cmboDataToView", CStr(PAT_SUM_DEFAULTS.cmboDataToView)
SaveINISetting gsAppIniFileSpec, "Patient Summary Settings", "cmboChartType", CStr(PAT_SUM_DEFAULTS.cmboChartType)

```

```
End Sub
```

```

Public Function FindFirstMatchingDateInArray(DataStruct As DeviceDataStruct, ByVal IBeginDate As Long)
'Find the earliest event date in the global structure that starts on the same day
'as the date passed here. Return 0 if not found or return the index to the date
'if one is found.

```

```

'Note that this date is not necessarily a dosing event date. It could be any kind of event

```

```

'Conduct a successive approximation lookup of the data in the array
Dim I As Integer, iLowIndex As Integer, iHighIndex As Integer, iTestIndex As Integer

```

```

iLowIndex = 1 'start at begin of array
iHighIndex = DataStruct.iEventData(0) 'stop at end of array
iTestIndex = (iHighIndex + iLowIndex) / 2
For i = 1 To 7 'this number of tries is all that is necessary to find the date
    If IBeginDate <= Int(DataStruct.dEventDate(iTestIndex)) Then
        iHighIndex = iTestIndex
    ElseIf IBeginDate > Int(DataStruct.dEventDate(iTestIndex)) Then
        iLowIndex = iTestIndex
    ElseIf IBeginDate = Int(DataStruct.dEventDate(iTestIndex)) Then
        iHighIndex = iTestIndex
        FindFirstMatchingDateInArray = iTestIndex
    End If
    iTestIndex = (iHighIndex + iLowIndex + 0.5) / 2
Next i

```

```
Next I
```

```
End Function
```

```

Public Function FindClosestDateInArray(DataStruct As DeviceDataStruct, ByVal IFromDate As Long) As Long
'Find the latest event date in the global structure that starts on the same day
'as the date passed here. Return 0 if not found or return the index to the date
'if one is found.

```

```

'If a 0 value is passed here then find the most recent date in the array.
'Note that this date is not necessarily a dosing event date. There is a separate
'procedure to find that date.

```

```

'If necessary for faster speed, this procedure can be recoded to do a successive approximation

```

```

Dim I As Integer
If IFromDate = 0 Then IFromDate = 99999
For i = 1 To DataStruct.iEventData(0) 'find the latest date in the array
    If IFromDate <= Int(DataStruct.dEventDate(i)) Then
        FindClosestDateInArray = i
    End If
Next i

```

```
End Function
```

General.bas - SaveINISetting

24

```
Public Sub SaveINISetting(ByVal sFileName As String, ByVal sSection As String, ByVal sKeyField As String, s
```

```
Dim i As Long
```

```
    i = WritePrivateProfileString(sSection, sKeyField, sValue, sFileName)
```

```
End Sub
```

```
Public Function ValidateDoseNumbers(frmTarget As Form)
```

```
    'Ensure that there are at least as many dose times as there are
```

```
    'for the number of doses per day.
```

```
    Dim i As Integer, iDailyDoseCounts As Integer, iDosesPerDay As Integer
```

```
    If Len(PAT_DATA.sPatientDataFileName) = 0 And Len(PAT_DATA.sSerialNumber) = 0 Then
```

```
        ValidateDoseNumbers = True
```

```
        Exit Function
```

```
    End If
```

```
    With frmTarget
```

```
        iDosesPerDay = Val(.txtDosesPerDay)
```

```
        For i = 1 To 4
```

```
            If IsDate(.txtDoseTime(i)) Then
```

```
                iDailyDoseCounts = iDailyDoseCounts + 1
```

```
            End If
```

```
        Next i
```

```
        If iDailyDoseCounts = iDosesPerDay Then
```

```
            ValidateDoseNumbers = True
```

```
        Else
```

```
            Beep
```

```
            MsgBox "You have indicated " + CStr(iDosesPerDay) + " Doses Per Day, yet " + CStr(iDailyDoseCounts) + " Dose Times were
```

```
            entered. They must match.", vbExclamation, "Mis-matched Dosing Values"
```

```
            .txtDosesPerDay.SetFocus
```

```
            .UpDownDoseTime(4).SetFocus
```

```
        End If
```

```
    End With
```

```
End Function
```

```
Public Function ValidatePatientDataSaved()
```

```
    'Ensure that the patient data in memory is saved before proceeding to load new data from device
```

```
    'Return true if successful, else return vbCancel if user cancelled
```

```
    Dim r As Integer
```

```
    ValidatePatientDataSaved = True    'this is the default condition unless set otherwise below
```

```
    If gbPatientDataNotSaved Then
```

```
        Beep
```

```
        r = MsgBox("The patient data currently in memory has not been saved. Do you want to save it?", vbYesNoCancel + vbQuestion, "
```

```
        Patient Data Not Saved")
```

```
        If r = vbYes Then
```

```
            r = SaveDataToNewFile()
```

```
            If r = vbCancel Then ValidatePatientDataSaved = vbCancel    'cancelled from the save as dialog
```

```
        ElseIf r = vbCancel Then
```

```
            ValidatePatientDataSaved = vbCancel    'cancelled from message box
```

```
        End If
```

```
    End If
```

Comm.bas - File Declarations

26

```

Attribute VB_Name = "modComm"
Option Explicit
'Global definitions for device communication.
'by Glen Hamilton 10/5/97

'for RS232 communication
Public gbCommTimerExpired As Integer 'this flag is set when the comm timer expires
Public giCommPort As Integer 'Communication Port #
Public gsCommDeviceSettings As String 'speed settings (ie 2400,N,8,2)
Public gbCommReplyPending As Boolean 'a command was just sent, and reply is pending
Public gbCommBusy As Boolean 'a command is in progress. Get's cleared when reply is received or times out
Public gbCommOK As Integer 'needs to be an integer (no boolean) keeps current status of communications. false = no
comm, true = comm ok, any other value for simulation
Public giDeviceResponseWait As Integer 'milliseconds to wait for next char before assuming end of received string
Private Const ERR_COMM_BADRESPONSE = 31001
Private Const ERR_COMM_TIMEOUT = 30998
Private Const ERR_COMM_STRINGLENGTH = 30997
Private Const ERR_COMM_BUSY = 30996
Private Const ERR_COMM_CHECKSUM = 30995

Public Const ERR_DATA_CHECKSUM = 99997
Public Const ERR_NEWER_HOST_SOFTWARE = 99998 'set when device returns custom data that was saved with a newer revision
level

Device communications
Public gbKeepPollingDevice As Boolean 'when true, continuous polling of device is done

Define some application specific variables & constants
Public gsAppIniFileSpec As String

Type DeviceDataStruct
sPatientLastName As String '(16 bytes) uses 1st 16 byte block of the patient/pharmacy ID & Names
sPatientFirstName As String '(16 bytes) uses 1st 16 byte block of the patient/pharmacy ID & Names
sPatientID As String
sDrug As String '(16 bytes) uses 2nd 16 byte block of the patient/pharmacy ID & Names
sOrgan As String '(16 bytes) uses 3rd 16 byte block of the patient/pharmacy ID & Names
sTxCenter As String '(16 bytes) uses 3rd 16 byte block of the patient/pharmacy ID & Names
sSerialNumber As String '(10 bytes) device serial number
sFirmwareVer As String 'Rev version and date of firmware
sDoseSize As String '(1 byte) stored here in mg. The device uses "ml" (100mg = 1 ml)
'Device Dose size is in optical ticks (0 to 200) max dose = 5 ml.

sPatientDataFileName As String 'file path and filename of the data in memory

Note: the daily prescribed dosing times below are stored in fractional days. This is done
to speed display operations and reduce the amount of memory needed. The device actually stores
these values as intervals relative to 1:00 in the morning. Thus, the times are converted to
intervals when communicating with the device.
dPrescribedDoseTime(4) As Double 'doses due during the day (prescribed) usually a max of four
iDosesPerDay As Integer '(1 byte) # of doses per day (1 to 4)
sDoseResolution As String '(1 byte) Called "Dose Conversion" in firmware.
'Optical ticks to mg multiplier. (IE 2 ticks = 10 mg.)
'Optical ticks are fixed at 0.05 ml per tick.
sMedRemaining As String '(2 bytes) Medication "Supply volume" remaining (in optical ticks)
sScoreData(14) As String 'Today's score(14 bytes for all scores) of last 14 days doses taken. Circular buffer.
'valid data is value from 0-4 representing number of doses taken each day.
'Note: The "Score pointer" points to the current day.

Note that the following arrays can not be larger than 1500 events or else the space limit
of 64K will be exceeded. If necessary in the future to have more events than this for
a single file then make a separate array for the diagnostic data or temp data.
iEventData(1400) As Integer 'the data occurring for the event data. Might be a dose size, error flags, etc.
dEventData(1400) As Double 'list of dose days in order of first taken to most recent
byteEventType(1400) As Byte 'value =0 if it is a dose taken
'value = 1 if data event is a dose command change
'value = 2 if user entered entered
sUserData1(1400) As String 'user entered data in the first column of the gnd
sUserData2(1400) As String 'user entered data in the first column of the gnd

```


Comm.bas - File Declarations

27

```

sUserData3(1400) As String      'user entered data in the first column of the gnd

sClock As String                '(2 bytes) 10 minute resolution "0000" = 1 am on first dose day.
!DeviceInitDate As Long         'date the device started
!DeviceRefDateTime As Double    'date and time that all events are referenced to.
sBatteryChangeTimer As String   '(2 bytes) Battery change timer, in 10 minute increments
sDoseLockoutHours As String     '(1 byte) Hours to lockout dosing after a dose is taken
bErrorFatal As Boolean          'true if this flag was set in the returned flags string
bErrorNonFatal As Boolean       'true if this flag was set in the returned flags string
bErrorDoseSize As Boolean       'true if this flag was set in the returned flags string
bErrorMedRemaining As Boolean   'true if this flag was set in the returned flags string
bErrorMemoryFull As Boolean     'true if this flag was set in the returned flags string
bErrorBrownOut As Boolean       'true if this flag was set in the returned flags string
bErrorsExist As Boolean         '(1 byte) Bits are set if various errors have occurred and have not
                                'been corrected. A value of "0" is normal (no errors). Errors
                                'are corrected by either correcting the specific situation or
                                'resetting & reloading the dosing parameters.
                                'B0=1 if fatal system failure
                                'B1=1 if non-fatal system failure has occurred
                                'B2=1 if error has occurred in Dose Size Volume
                                'B3=1 if error has occurred in Supply Volume value
                                'B4=1 if compliance memory is near full
                                'B5=1 if brownout (low voltage) occurred

!LastDownloadDate As Double     'date of last data retrieval from device

End Type
Public PAT_DATA As DeviceDataStruct
Public TEMP_DATA As DeviceDataStruct

Public gsDrugNames(25) As String 'names of drugs used to populate the list boxes on dialogs
Public gsOrganNames(25) As String 'names of gsOrganNames used to populate the list boxes on dialogs

Public Const g!EVENT_DOSE_TAKEN = 0
Public Const g!EVENT_DOSE_CHANGED = 1
Public Const g!EVENT_USER_DEFINED = 2

'These values indicate the string position (returned from the device) where each element
'begins. This is the string that is returned when a request for "all memory" is sent.
'See above structures for more detail information about format.
Public Const DATA_BEGIN_DOSE_SIZE = 1 '1 byte
Public Const DATA_BEGIN_DOSE_INTERVAL1 = 1 * 2 + 1 '1 byte
Public Const DATA_BEGIN_DOSE_INTERVAL2 = 2 * 2 + 1 '1 byte
Public Const DATA_BEGIN_DOSE_INTERVAL3 = 3 * 2 + 1 '1 byte
Public Const DATA_BEGIN_DOSE_INTERVAL4 = 4 * 2 + 1 '1 byte
Public Const DATA_BEGIN_DOSES_PER_DAY = 5 * 2 + 1 '1 byte
Public Const DATA_BEGIN_DOSE_CONVERSION = 6 * 2 + 1 '1 byte
Public Const DATA_BEGIN_DOSE_LOCKOUT_HOURS = 7 * 2 + 1 '1 byte
Public Const DATA_BEGIN_DOSE_SCORE_DAY_POINTER = 8 * 2 + 1 '1 byte
Public Const DATA_BEGIN_MED_REMAINING = 9 * 2 + 1 '2 bytes
Public Const DATA_BEGIN_CLOCK = 11 * 2 + 1 '2 bytes clock starts at 1am on first dosing day (10 min increments)
Public Const DATA_BEGIN_BATTERY_CHANGE_TIMER = 13 * 2 + 1 '2 bytes
Public Const DATA_BEGIN_ERROR_FLAGS = 15 * 2 + 1 '1 byte
Public Const DATA_BEGIN_PREV_DOSE_PARAMS = 16 * 2 + 1 '16 bytes of copy of prev dosing params
                                'Used for error checking internal to dispenser!
Public Const DATA_BEGIN_KEY_BITS = 32 * 2 + 1 'activation of keys on device
Public Const DATA_BEGIN_LIFE_COUNT = 33 * 2 + 1 'LSB in 33, MSB in 34
Public Const DATA_BEGIN_LIFE_COMPLETION = 35 * 2 + 1 '=0 when life cycle is programmed, =1 when life test completes
                                successfully
Public Const DATA_BEGIN_COMPENSATION_FACTOR = 36 * 2 + 1 'values from 64-192 (128= 1.0 factor)
Public Const DATA_BEGIN_SERIAL_NUMBER = 38 * 2 + 1 '10 bytes
Public Const DATA_BEGIN_CUSTOM1 = 48 * 2 + 1 '16 bytes of patient/pharmacy ID & names
Public Const DATA_BEGIN_CUSTOM2 = 64 * 2 + 1 '16 bytes of patient/pharmacy ID & names
Public Const DATA_BEGIN_CUSTOM3 = 80 * 2 + 1 '16 bytes of patient/pharmacy ID & names
Public Const DATA_BEGIN_CUSTOM4 = 96 * 2 + 1 '16 bytes of patient/pharmacy ID & names
Public Const DATA_BEGIN_SCORE = 112 * 2 + 1 '14 bytes
Public Const DATA_BEGIN_COMPLIANCE_CHECKSUM = 128 * 2 + 1 '2 bytes Includes compliance pointer and data
                                'Up to data word 1 before data pointed to by comp pointer

```

Comm.bas - File Declarations

28

```
Public Const DATA_BEGIN_COMPLIANCE_POINTER = 130 * 2 + 1    '2 bytes points to next location after end of current compliance
    Date
    'base value = 132 (0x0084)
Public Const DATA_BEGIN_COMPLIANCE_DATA = 132 * 2 + 1    '~1900 bytes max. Array of 2 byte values for dose compliance
    history.
    'Clock time values (in 10 minutes resolution from start) when each dose
    'was taken. Represented by values 0-65279 (0-0x1fff). Each dose time
    'is changed via the "Set Mode", a value between 0x1f00 and 0x1fff is written
    'with the LSD byte representing the dose size. When compliance memory is
    'cleared, the current dose size is always written as the first location in
    'the compliance memory.
```

Public Sub ChangeBatteriesRequest()

```
Dim r As Integer, IErrorCode As Long, sMSG As String
```

```
sMSG = "You should continue only if you are replacing the batteries in the device."
sMSG = sMSG + " This ensures that the battery time counter will be accurate." + vbCrLf + vbCrLf
sMSG = sMSG + " Did you just replace the batteries or are you about to change them now?"
r = MsgBox(sMSG, vbQuestion + vbYesNo + vbDefaultButton2, "Change Batteries")
```

```
If r = vbNo Then Exit Sub
```

```
gbKeepPollingDevice = False    'stop polling for now
Wait 0.25
```

```
On Error GoTo ChangeBatteriesRequest_Error
```

```
r = Comm_SendResetClockAndBattery(IErrorCode)
```

```
If IErrorCode Then
```

```
    Error IErrorCode    'error number
```

```
Else
```

```
    sMSG = "Replace the device batteries now and retrieve data from the device again when complete."
    r = MsgBox(sMSG, vbExclamation, "Change Batteries")
End If
```

```
ChangeBatteriesRequest_Exit:
gbKeepPollingDevice = True    'continue polling device
Exit Sub
```

```
ChangeBatteriesRequest_Error:
DisplayErrorMessage IErrorCode
Resume ChangeBatteriesRequest_Exit
```

```
End Sub
```

Public Function Comm_CheckComm(IErrorCode As Long) As Integer

*'Check the device communication by sending a command and waiting for a reply.
'If no reply is received, then return a "false" flag to caller.*

*'Important Note: Due to the way the firmware was designed for the device, it seems not
'to return anything if the command is in error. This is not good because
'we would not know whether or not a failed reply is due to a back cable, incorrect comm port
'or settings, etc. Hopefully in a future version, the comm check can return some sort of
'character to indicate that a common byte was received, but could not be interpreted correctly.*

```
Dim sOut As String, sChecksum As String, sIn As String
```

```
sOut = "Pp"    'this is the code for checking communication with device
CreateChecksum sOut, sChecksum    'calculate a checksum
sOut = sOut + sChecksum + "T"    'append checksum and ending string identifier
```

```
If Not frmMain.CommDevice.PortOpen Then frmMain.CommDevice.PortOpen = True
```

```
gbCommBusy = True    'prevent other procedures from communicating with device
frmMain.CommDevice.InputLen = 0    'clear input buffer
```

Comm.bas - Comm_CheckComm

29

```

frmMain.CommDevice.Output = sOut      'send string to device
gbCommReplyPending = True             'prevent other procedures from communicating with device
SetCommTimer giDeviceResponseWait     'set timer to wait for response

IErrorcode = 0                        'reset error code
Do
    If gbCommTimerExpired Then         'timer event sets this to true
        IErrorcode = ERR_COMM_TIMEOUT 'no response, get the error code
        GoTo Comm_CheckComm_Exit      'return to calling procedure
    End If
    DoEvents
Loop Until frmMain.CommDevice.InBufferCount > 0 'loop till a reply is received or timeout occurs

sIn = frmMain.CommDevice.Input        'Recd response from serial port
'comm is ok
If sIn = "$" Then Comm_CheckComm = True 'return success to caller

Comm_CheckComm_Exit:
If frmMain.CommDevice.PortOpen Then frmMain.CommDevice.PortOpen = False 'Close the serial port
gbCommReplyPending = False          'reset flag
gbCommBusy = False                  'reset flag

End Function

```

Public Function Comm_GetDeviceReply(sReply As String, IErrorcode As Long) As Integer

'A command should have been just sent to the device from another procedure and a reply is pending.
 'Get the reply into 'sReply' and return to caller.
 'Return false if no reply and set IErrorcode to reason.
 'Return ERR_COMM_TIMEOUT if no reply from the device.
 'error code = 0 if comm is already busy.
 'If reply, then return number of characters received.
 'Close comm port once a reply is received or if an error occurs.

Dim iLastBufferCount As Integer, r As Integer

On Error GoTo Comm_GetDeviceReply_Error
 gbCommReplyPending = True 'set busy flag
 frmMain.MousePointer = vbHourglass

'Open comm port in case it is closed

'prevent device unavailable error
 If frmMain.CommDevice.PortOpen = False Then frmMain.CommDevice.PortOpen = True 'open port
 sReply = "" 'init reply

'Wait for first char to arrive

SetCommTimer giDeviceResponseWait '20 milliseconds is normally sufficient

Do Until frmMain.CommDevice.InBufferCount > 0

DoEvents

'timer event sets this to true

If gbCommTimerExpired Then Error ERR_COMM_TIMEOUT 'return message to caller. No response

Loop

'First char has been received

'Wait for all data to arrive

iLastBufferCount = -1

'init buffer count

Do While frmMain.CommDevice.InBufferCount > iLastBufferCount 'characters are still arriving

iLastBufferCount = frmMain.CommDevice.InBufferCount 'remember intermediate count

SetCommTimer giDeviceResponseWait

'this value works as low as 25 milliseconds

Do Until gbCommTimerExpired = True

'wait for timer to expire

DoEvents

DoEvents

Loop

Loop

'loop if characters are still coming in

'All data has arrived or time has been too long

Comm.bas - Comm_GetDeviceReply

30

```

r = frmMain.CommDevice.InBufferCount      'get character length
sReply = frmMain.CommDevice.Input         'read string from buffer
Comm_GetDeviceReply = Len(sReply)         'return buffer length to caller (- checksum)

```

```

Comm_GetDeviceReply_Exit:
'prevent device unavailable error
If frmMain.CommDevice.PortOpen = True Then frmMain.CommDevice.PortOpen = False 'close port if open
frmMain.MousePointer = vbDefault
On Error GoTo 0 'clear error status
gbCommReplyPending = False 'reset pending flag
gbCommBusy = False 'reset busy flag
Exit Function

```

```

Comm_GetDeviceReply_Error:
IErrorCode = Err 'return error code to caller
'Resume 0 'for testing only
Resume Comm_GetDeviceReply_Exit
End Function

```

Public Function Comm_ReadFirmwareVersion(DataStruct As DeviceDataStruct, IReturnError As Long) As Integer
 Dim sOut As String, sChecksum As String, sIn As String, IErrorCode As Long, r As Integer

```

sOut = "" 'this is the code for version number
CreateChecksum sOut, sChecksum 'calculate a checksum
sOut = sOut + sChecksum + "I" 'append checksum and ending string identifier

```

```

If Not frmMain.CommDevice.PortOpen Then frmMain.CommDevice.PortOpen = True

```

```

frmMain.CommDevice.InputLen = 0 'clear input buffer
frmMain.CommDevice.Output = sOut 'send string to device
r = Comm_GetDeviceReply(sIn, IErrorCode)

```

```

If IErrorCode = 0 Then 'comm was received
  r = ValidateChecksum(sIn)
  If r Then
    DataStruct.sFirmwareVer = Left$(sIn, Len(sIn) - 5) 'put string in global array
    Comm_ReadFirmwareVersion = True 'return success to caller
  End If

```

```

Else
  IReturnError = IErrorCode
  DisplayErrorMessage IErrorCode
End If
End Function

```

Public Sub DisplayCommError(SourceForm As Form)

```

gbCommOK = False
SourceForm.imgCommStatus.Picture = SourceForm.imgRedLight
SourceForm.lblCommStatus = "No Device Found"

```

```

'Play disconnect sound and show status visually
'Set properties needed by MCI to open
With SourceForm.MMControl1
  .Notify = False
  .Wait = False
  .Shareable = False
  .filename = App.Path + "\ProbDetectVoice.wav"

```

```

'Open the MCI WaveAudio device
.Command = "Open"
.Command = "sound"
.Command = "close"
End With

```

Comm.bas - DisplayCommError

31

End Sub

Public Sub DisplayCommOk(SourceForm As Form)

*Play connect sound and show status visually**Set properties needed by MCI to open*

gbCommOK = True

With SourceForm.MMControl1

.Notify = False

.Wait = False

.Shareable = False

.filename = App.Path & "\morsecode.wav"

Open the MCI WaveAudio device

.Command = "Open"

.Command = "sound"

.Command = "sound"

.Command = "close"

End With

SourceForm.imgCommStatus.Picture = SourceForm.imgGreenLight

SourceForm.lblCommStatus = "Device Ready"

End Sub

Public Sub DisplayErrorMessage(ErrorCode As Long)

Dim sMSG As String

Select Case ErrorCode

Case ERR_COMM_TIMEOUT

sMSG = "No response was received from the device to the command just issued. Remove the device from the communicator and re-insert it to ensure that it is seated properly."

Case ERR_COMM_CHECKSUM

sMSG = "Data retrieved from the device is corrupted. This probably occurred during transmission. Please read the device again."

Case ERR_COMM_BADRESPONSE

sMSG = "The device did not interpret the command properly."

Case ERR_NEWER_HOST_SOFTWARE

sMSG = "The device was previously programmed with a newer version of this software. The data can not be retrieved." + vbCrLf + "Please obtain an updated version this software."

Case Else

sMSG = "An error was detected while communicating with the device. Please try again." + vbCrLf + vbCrLf + Error\$(ErrorCode)

End Select

MsgBox sMSG, , App.ProductName & " Comm Error - " & CStr(ErrorCode)

End Sub

Comm.bas - GetDrugRefNumber

32

Public Function GetDrugRefNumber() As Integer*'Find the index to the organ name being using in the global structure*

Dim i As Integer

For i = 1 To UBound(gsDrugNames)

If LCase(PAT_DATA.sDrug) = LCase(gsDrugNames(i)) Then Exit For

Next i

GetDrugRefNumber = i *'return ref number to caller*

End Function

Public Function GetOrganRefNumber() As Integer*'Find the index to the organ name being using in the global structure*

Dim i As Integer

For i = 1 To UBound(gsOrganNames)

If LCase(PAT_DATA.sOrgan) = LCase(gsOrganNames(i)) Then Exit For

Next i

GetOrganRefNumber = i *'return ref number to caller*

End Function

Private Function InterpretDosingData(DataStruct As DeviceDataStruct, ByVal sData As String, ByVal dChecksumTally As Double, ByVal iChecksum As Long, iErrorCode As Long) As Integer*'Parse apart the dosing data that is passed here and put into global structure.**'Each dosing event is 2 bytes in length.**'The checksum is passed here for comparison to the string.**'The checksum includes the pointer bytes which is why it is passed in.*

Dim sTemp As String, iLowByte As Integer, iHiByte As Integer, iCurrentDoseAmount As Long

Dim r As Integer, iPosition As Integer, iCount As Integer, iTemp As Long

On Error GoTo InterpretDosingData_Error

For iPosition = 1 To Len(sData) Step 4 *'there are 2 hex bytes to every dose event*

iCount = iCount + 1

iLowByte = CInt("&H" + Mid\$(sData, iPosition, 2)) *'get first of the 2 byte hex value*dChecksumTally = dChecksumTally + iLowByte *'add to checksum*iHiByte = CInt("&H" + Mid\$(sData, iPosition + 2, 2)) *'get second of the 2 byte hex value*dChecksumTally = dChecksumTally + iHiByte *'add to checksum*If iHiByte = 255 Then *'indicates a dose change*

DataStruct.byteEventType(iCount) = giEVENT_DOSE_CHANGED

iCurrentDoseAmount = CLng(iLowByte) / 40 * 100 *'convert from ml to mg; this is a new dose size change*

DataStruct.iEventData(iCount) = iCurrentDoseAmount

DataStruct.dEventData(iCount) = CLng(DataStruct.dEventData(iCount - 1)) *'get date from last dose**'This is the very first dose change*

If DataStruct.dEventData(iCount) = 0 Then

DataStruct.dEventData(iCount) = 1

End If

Else *'this is a dose*

DataStruct.byteEventType(iCount) = giEVENT_DOSE_TAKEN

DataStruct.iEventData(iCount) = iCurrentDoseAmount *'convert from ml to mg; this is a new dose size change*iTemp = CLng(iHiByte) * 256 + iLowByte *'date and time is 10 minute intervals since first dose day*

DataStruct.dEventData(iCount) = DateAdd("n", iTemp * 10, DataStruct.dDeviceRefDateTime)

End If

Next iPosition

DataStruct.iEventData(0) = iCount *'put the total number of events in the 0 element of list*

dChecksumTally = dChecksumTally Mod 65536

If dChecksumTally = iChecksum Then

Comm.bas - InterpretDosingData

33

```

    InterpretDosingData = True 'return success to caller
Else
    IErrorCode = ERR_COMM_CHECKSUM
End If

InterpretDosingData_Exit:
    On Error GoTo 0
    Exit Function

InterpretDosingData_Error:
    IErrorCode = Err
    Resume InterpretDosingData_Exit:

End Function

Public Function ValidateChecksum(ByVal sData As String) As Integer
    'Look at the data string passed here and get the checksum from the
    'end of the string.
    'sData should be a string that was returned from the device.
    'The last char in the string is a termination char preceded by 4 bytes of checksum.

    Dim sTemp As String, iByte As Integer, r As Integer, iPosition As Integer
    Dim iChecksum As Long, iChecksumTally As Long
    On Error GoTo ValidateChecksum_Error

    r = Len(sData)
    For iPosition = 1 To r - 5
        iByte = Asc(Mid$(sData, iPosition, 1))
        iByte = CInt("&H" + Mid$(sData, iPosition, 2))
        iChecksumTally = iChecksumTally + iByte 'add to checksum
    Next iPosition

    iChecksumTally = iChecksumTally Mod 65536

    sTemp = "&H" + "0" + Mid$(sData, r - 2, 2) + Mid$(sData, r - 4, 2)
    iChecksum = CLng(sTemp)
    If iChecksumTally = iChecksum Then ValidateChecksum = True 'pass success flag back to caller

ValidateChecksum_Exit:
    On Error GoTo 0
    Exit Function

ValidateChecksum_Error:
    Resume ValidateChecksum_Exit

End Function

Private Sub InterpretErrorFlags(DataStruct As DeviceDataStruct, ByVal iFlagsByte As Integer)
    'Break out the bits of the flags bytes passed here.
    'Put the results into the global arrays

    'If any flags exist, then set this to true
    If iFlagsByte Then DataStruct.bErrorsExist = True

    'Parse out flags separately
    DataStruct.bErrorFatal = (iFlagsByte And 2)
    DataStruct.bErrorNonFatal = (iFlagsByte And 4)
    DataStruct.bErrorDoseSize = (iFlagsByte And 8)
    DataStruct.bErrorMedRemaining = (iFlagsByte And 16)
    DataStruct.bErrorMemoryFull = (iFlagsByte And 32)
    DataStruct.bErrorBrownOut = (iFlagsByte And 64)
    'remaining upper 3 bits not used at present

End Sub

```

Comm.bas - Comm_ReadEntireMemoryContents

34

Public Function Comm_ReadEntireMemoryContents(DataStruct As DeviceDataStruct, IReturnError As Long) As String
 Dim sOut As String, sChecksum As String, sin As String, IErrorCode As Long, r As String

On Error GoTo Comm_ReadEntireMemoryContents_Error

EraseDataInMemory DataStruct

sOut = "R" *'this is the code for reading entire memory*

CreateChecksum sOut, sChecksum *'calculate a checksum*

sOut = sOut + sChecksum + "I" *'append checksum and ending string identifier*

If Not frmMain.CommDevice.PortOpen Then frmMain.CommDevice.PortOpen = True

frmMain.CommDevice.InputLen = 0 *'clear input buffer*

frmMain.CommDevice.Output = sOut *'send string to device*

r = Comm_GetDeviceReply(sin, IErrorCode)

If r Then *'comm was received. Should be at least this many bytes*

r = ValidateChecksum(sin)

If r = False Then

IReturnError = ERR_COMM_CHECKSUM

Exit Function

End If

r = ParseMemoryContents(DataStruct, sin, IErrorCode) *'parse out the string*

If IErrorCode Then

IReturnError = IErrorCode

Exit Function

End If

Comm_ReadEntireMemoryContents = True *'return success to caller*

DataStruct.LastDownloadDate = Now

gbPatientDataNotSaved = True *'set this flag to true*

Else

IReturnError = IErrorCode

Exit Function

End If

r = Comm_ReadFirmwareVersion(DataStruct, IErrorCode)

If IErrorCode Then

IReturnError = IErrorCode

End If

Comm_ReadEntireMemoryContents_Exit:

On Error GoTo 0

Exit Function

Comm_ReadEntireMemoryContents_Error:

IReturnError = Err

Resume Comm_ReadEntireMemoryContents_Exit

End Function

Comm.bas - Comm_SendResetClockAndBattery

35

Public Function Comm_SendResetClockAndBattery(IReturnError As Long)

*'Resets the device clock to an offset that represents 1:00am
'and resets the battery timer to zero.*

Dim sData As String, sOut As String, sReply As String, sChecksum As String
Dim r As Integer, IErrorCode As Long, ITemp As Integer

ITemp = CInt(Format(Now, "hh"))
If ITemp = 0 Then ITemp = 24 *'midnight*
ITemp = (ITemp - 1) * 6 *'calc number of 10-minute period ~ hours*
ITemp = ITemp + CInt((Format\$(Now, "nn") - 5) / 10) *'calc number of 10-min periods in this hour*
sData = CStr(Hex(ITemp))
If Len(sData) < 2 Then sData = "0" + sData *'ensure string is always 2 bytes long*

sOut = "0x" + sData *'add data string to command*
CreateChecksum sOut, sChecksum *'calculate a checksum*
sOut = sOut + sChecksum + "T" *'append checksum and ending string identifier*
Comm_SendDataToDevice(sOut) *'send string to comm port*
r = Comm_GetDeviceReply(sReply, IErrorCode)
If sReply = "\$" Then *'string was successfully interpreted by device*
Comm_SendResetClockAndBattery = True *'return success to caller*
ElseIf sReply = "T" Then *'string was not interpreted properly*
IReturnError = ERR_COMM_BADRESPONSE
Else
IReturnError = IErrorCode
End If

End Function

Public Function Comm_SendCustomData(DataStruct As DeviceDataStruct, ByVal sLocation As String, IReturn

*'There are 4 locations in the device, each containing a 16 byte string.
'There first location is usually reserved for Patient ID.
'Any string can be contained in any location.
'Data is taken from the global structure*

Dim sData As String, sOut As String, sReply As String, sChecksum As String
Dim r As Integer, IErrorCode As Long, sCustomData As String
Dim l As Integer, sTemp As String

*'Determine the appropriate command for the location of the data to be stored.
'There are 4 fields in the device containing 16 characters each. In the original
'device design, this was intended to contain 4 separate pieces of information.
'The client has now decided that some fields are too short and others are too long.
'Thus, the fields are combined to be one string of 64 characters.
'Data Structure Rev level = gLEN_REV_DATA_STRUCTURE
'Patient name = gLEN_PATIENT_NAME
'ID = gLEN_ID
'Drug = gLEN_DRUG
'TX Center = gLEN_TX_CENTER
'Organ = gLEN_ORGAN*

*'Create a 64 byte string from the various data elements to be saved
'This string identifies the format of custom information. If the format
'changes in a future version, then this ID can be used to determine which
'version of the software saved the info to the device.
sCustomData = gsREV_DATA_STRUCTURE*

*'Save the 2 digit number that represents this drug.
'To save spec, a numerical index of the Organ name is stored in the device
l = GetDrugRef(Number)
sTemp = CStr(l)
If Len(sTemp) < 2 Then sTemp = "0" + sTemp *'force code to be 2 digits*
sCustomData = sCustomData + sTemp *'concatenate result to outbound string**

*'Save the 2 digit number that represents this organ
'To save spec, a numerical index of the Organ name is stored in the device*

Comm.bas - Comm_SendCustomData

36

```

I = GetOrganRefNumber()
sTemp = CStr(I)
If Len(sTemp) < 2 Then sTemp = "0" + sTemp 'force code to be 2 digits
sCustomData = sCustomData + sTemp

sTemp = DataStruct.sPatientID
If Len(sTemp) > giLEN_ID Then
    sTemp = Left$(sTemp, giLEN_ID) 'name is too long, trim letters of first name
ElseIf Len(sTemp) < giLEN_ID Then
    sTemp = sTemp + Space$(giLEN_ID - Len(sTemp)) 'name is too short to fill the designated length
End If
sCustomData = sCustomData + sTemp 'concatenate result to outbound string

sTemp = DataStruct.sTxCenter
If Len(sTemp) > giLEN_TX_CENTER Then
    sTemp = Left$(sTemp, giLEN_TX_CENTER) 'name is too long, trim letters of first name
ElseIf Len(sTemp) < giLEN_TX_CENTER Then
    sTemp = sTemp + Space$(giLEN_TX_CENTER - Len(sTemp)) 'name is too short to fill the designated length
End If
sCustomData = sCustomData + sTemp 'concatenate result to outbound string

sTemp = DataStruct.sPatientLastName + " " + DataStruct.sPatientFirstName
If Len(sTemp) > giLEN_PATIENT_NAME Then
    sTemp = Left$(sTemp, giLEN_PATIENT_NAME) 'name is too long, trim letters of first name
ElseIf Len(sTemp) < giLEN_PATIENT_NAME Then
    sTemp = sTemp + Space$(giLEN_PATIENT_NAME - Len(sTemp)) 'name is too short to fill the designated length
End If
sCustomData = sCustomData + sTemp 'concatenate result to outbound string

'Assemble the string to be sent
Select Case sLocation
    Case DATA_BEGIN_CUSTOM1
        sOut = "Vw"
        sData = Mid$(sCustomData, 1, 16)
    Case DATA_BEGIN_CUSTOM2
        sOut = "Xx"
        sData = Mid$(sCustomData, 17, 16)
    Case DATA_BEGIN_CUSTOM3
        sOut = "Yy"
        sData = Mid$(sCustomData, 33, 16)
    Case DATA_BEGIN_CUSTOM4
        sOut = "Zz"
        sData = Mid$(sCustomData, 49, 16)
End Select

'Since the device string is limited to 16 chars, ensure that only the 1st 16 chars
'of the string are sent.
sOut = sOut + sData 'add data string to command
CreateChecksum sOut, sChecksum 'calculate a checksum
sOut = sOut + sChecksum + "T" 'append checksum and ending string identifier
Comm_SendDataToDevice (sOut) 'send string to comm port
r = Comm_GetDeviceReply(sReply, IErrorCode)
If sReply = "S" Then 'string was successfully interpreted by device
    Comm_SendCustomData = True 'return success to caller
ElseIf sReply = "T" Then 'string was not interpreted properly
    IReturnError = ERR_COMM_BADRESPONSE
Else
    IReturnError = IErrorCode
End If

End Function

```

Comm.bas - Comm_SendCustomData

37

Public Function Comm_SendDosingParams(DataStruct As DeviceDataStruct, IReturnError As Long)
Sends the dosing parameters from the global structure to the device

'Structure: "Dduuvvwwwxyzmmss"
'Response: "3"
'"d" is Dose Size in pump ticks. Note that pump ticks per milliliter is fixed at 40. Hex value from 0 to 0xFF. (Maximum dose size is currently 5 ml or 200 decimal)
'"uu", "vv", "ww", "xx" are four Dose Interval values in hours between doses. Hex values between 1 and 0xFF.
'"yy" is Number of Doses per day. Hex value from 1 to 4.
'"zz" is Pump Ticks per 10 mg Conversion value. Hex value. Typical value with present medication is 4 ticks per 10 milligrams.
'"mm" is Lockout Hours value. Number of hours to prevent dosing after a dose is taken.
'"ss" are the two checksum digits (hex) equal the one's complement value of the two's complement sum of the command characters and the data. The ASCII values are simply added together in an 8 bit sum, then one is subtracted (modulo 255).

Dim IData As Integer, sData As String, sOut As String, sReply As String
 Dim r As Integer, IErrorCode As Long, sChecksum As String, i As Integer
 Dim iLastIntervalSet As Integer

On Error GoTo 0

sOut = "Dd" *'put command in string*

'Get Dose Size in pump ticks
 IData = Val(DataStruct.sDoseSize) * 40 / 100 *'get dose size from global struct & convert to pump ticks (convert from mg to ml)*
 sData = CStr(Hex(IData)) *'convert value to a hex string*
 If Len(sData) < 2 Then sData = "0" + sData *'ensure string is always 2 bytes long*
 sOut = sOut + sData

iLastIntervalSet = 1 *'1:00 am is the ref time for the first dose*

'Get Dose intervals in hours
 For i = 1 To gMaxDoseTimes *'max doses per day*
 IData = 0 *'in case of conversion error, reset temp value*
 If DataStruct.dPrescribedDoseTime(i) > 0 Then *'a negative number indicates no time was set*
 IData = Format\$(DataStruct.dPrescribedDoseTime(i), "hh") *'convert fractional day to hours*
 IData = IData - iLastIntervalSet *'this time is relative to the last interval that was set*
 'change time to midnight
 If IData < 0 Then IData = 24 - Abs(IData)
 iLastIntervalSet = Format\$(DataStruct.dPrescribedDoseTime(i), "hh") *'the next interval is relative to the last one that is set*
 End If
 sData = CStr(Hex(IData)) *'convert value to a hex string*
 If Len(sData) < 2 Then sData = "0" + sData *'ensure string is always 2 bytes long*
 sOut = sOut + sData
 Next i

'Get number of doses per day
 sData = CStr(Hex(DataStruct.iDosesPerDay)) *'convert value to a hex string*
 If Len(sData) < 2 Then sData = "0" + sData *'ensure string is always 2 bytes long*
 sOut = sOut + sData

'Get conversion value
 IData = 0 *'in case of error, reset temp value*
 IData = CInt(DataStruct.sDoseResolution) *'get dose resolution from global struct*
 sData = CStr(Hex(IData)) *'convert value to a hex string*
 If Len(sData) < 2 Then sData = "0" + sData *'ensure string is always 2 bytes long*
 sOut = sOut + sData

'Get Dose lockout hours
 IData = 0 *'in case of error, reset temp value*
 IData = CInt(DataStruct.sDoseLockoutHours) *'get lockout hour from global struct*
 sData = CStr(Hex(IData)) *'convert value to a hex string*
 If Len(sData) < 2 Then sData = "0" + sData *'ensure string is always 2 bytes long*
 sOut = sOut + sData

CreateChecksum sOut, sChecksum *'calculate a checksum*
 sOut = sOut + sChecksum + "T" *'append checksum and ending string identifier*

Comm.bas - Comm_SendDosingParams

38

```

Comm_SendDataToDevice (sOut)      'send string to comm port
r = Comm_GetDeviceReply(sReply, IErrorCode)
If sReply = "S" Then              'string was successfully interpreted by device
    Comm_SendDosingParams = True  'return success to caller
ElseIf sReply = "T" Then          'string was not interpreted properly
    IReturnError = ERR_COMM_BADRESPONSE
Else
    IReturnError = IErrorCode
End If

```

End Function

Public Function Comm_SendSerialNumber(DataStruct As DeviceDataStruct, IReturnError As Long) As Integer
'Sends the serial number from the global structure to the device

```

Dim sData As String, sOut As String, sReply As String, sChecksum As String
Dim r As Integer, IErrorCode As Long

sData = Left$(Trim(DataStruct.sSerialNumber), 10) 'trim leading spaces and use first 16 chars
sData = sData + Space$(10 - Len(sData)) 'pad the string with spaces
sOut = "Nn" + sData 'add data string to command
CreateChecksum sOut, sChecksum 'calculate a checksum
sOut = sOut + sChecksum + "T" 'append checksum and ending string identifier
Comm_SendDataToDevice (sOut) 'send string to comm port
r = Comm_GetDeviceReply(sReply, IErrorCode)
If sReply = "S" Then          'string was successfully interpreted by device
    Comm_SendSerialNumber = True 'return success to caller
ElseIf sReply = "T" Then      'string was not interpreted properly
    IReturnError = ERR_COMM_BADRESPONSE
Else
    IReturnError = IErrorCode
End If

```

End Function

Private Sub ConvertHexStringToAscii(ByVal sData, sConverted As String)

*'String Data from the device is usually returned as Hex characters.
'Convert the string (passed in here) to an ASCII string and return to caller.
'Such strings are items like patient name, serial number, etc.*

```

Dim sTemp As String, i As Integer, iTemp As Integer
On Error Resume Next
sConverted = "" 'clear out any old string
For i = 1 To Len(sData) Step 2
    sTemp = "&H" + Mid$(sData, i, 2) 'get a 2 char hex byte from string
    sTemp = Chr$(sTemp) 'convert value to ASCII
    sConverted = sConverted + sTemp 'concatenate to existing string being built
Next i

```

End Sub

Comm.bas - CreateChecksum

39

Private Sub CreateChecksum(sOut As String, sChecksum As String)

*'The string "sOut" will be sent to the device by another procedure. Before it is sent, this procedure calculates a checksum and returns it to the caller
'Return the ASCII representation of the checksum value.*

Dim i As Integer, iChecksumTally As Long, iChecksumByteLow As Integer

*For i = 1 To Len(sOut) 'calculate checksum
iChecksumTally = iChecksumTally + Asc(Mid\$(sOut, i, 1))
Next i*

iChecksumByteLow = iChecksumTally Mod 256

iChecksumByteHigh = iChecksumTally \ 256 'not being using by the device

*sChecksum = Hex(iChecksumByteLow - 1) 'checksum is the "one's complement value of a two's complement checksum"
'value must always be 2 chars*

If Len(sChecksum) < 2 Then sChecksum = "0" + sChecksum 'place a leading "0" in front of checksum

End Sub

Public Sub EstablishDeviceComm()

*'This procedure continues to try and establish communication with the Device until it succeeds. When successful, control is returned to the calling procedure.
'The purpose of this procedure is to allow the user to try cable changes, device movement, etc. without having to continue pressing keys on the keyboard.*

Dim r As Integer, iErrorCode As Long

QueryDevice:

r = Comm_CheckComm(iErrorCode)

If r <> True Then

DoEvents 'allow other Windows events to be processed, so we don't lock up the computer

Wait 1 'wait an additional amount of time before trying

GoTo QueryDevice 'try comm again

End If

End Sub

Function Comm_InitializeCommPort() As Integer

*'Get the initial values from INI file and
'Initialize device comm port settings*

Dim iReply As Long

Const sSection = "Communications"

'Get the amount of time needed for a reply to be received from the device

giDeviceResponseWait = CInt(GetINISetting(gsAppIniFileSpec, sSection, "Device Response Wait", "50"))

If giDeviceResponseWait < 25 Then giDeviceResponseWait = 25 'set a minimum delay time

If giDeviceResponseWait > 500 Then giDeviceResponseWait = 500 'clamp upper limit

frmMain.CommTimer.Interval = giDeviceResponseWait 'set up the timer

'Get the comm port speed settings

giCommPort = CInt(GetINISetting(gsAppIniFileSpec, sSection, "Port", "1"))

If giCommPort = 0 Then giCommPort = 2 'set a default of comm 2 if nothing is in the file

gsCommDeviceSettings = GetINISetting(gsAppIniFileSpec, sSection, "Settings", "2400,n,8,2")

'prevent device unavailable error

If frmMain.CommDevice.PortOpen = True Then frmMain.CommDevice.PortOpen = False

'close port if open

frmMain.CommDevice.Settings = gsCommDeviceSettings

frmMain.CommDevice.CommPort = CStr(giCommPort)

Comm.bas - Comm_InitializeCommPort

40

```

frmMain.CommDevice.InBufferSize = 1024
frmMain.CommDevice.InputLen = 0
Comm_InitializeCommPort = True           'return success to caller

```

```

Comm_InitializeCommPort_Exit:
Exit Function

```

```

Comm_InitializeCommPort_Error:
Comm_InitializeCommPort = Err           'return error to caller
On Error GoTo 0
Resume Comm_InitializeCommPort_Exit

```

```

End Function

```

Sub Device_OnComm()

This procedure is called by the OnComm event of the comm control located on frmMain. This is so that the code can be shared between applications.

```

Dim r As Integer, sTemp As String

```

```

r = frmMain.CommDevice.CommEvent

```

```

If r = MSComm_ER_RXOVER Then
    'An over run error occurred. Usually happens when getting events.
    Exit Sub
End If

If r = MSComm_EV_EOF Then
    'end of file flag received.
    Exit Sub
End If

If r = MSComm_ER_BREAK Then
    'break signal received
    Exit Sub
End If

If r = 3 Or r = 4 Or r = 5 Then
    'rts, xon xoff, CD error
    Exit Sub
End If

MsgBox "Unexpected error occurred with the device. Please try again.", "Comm Event - " & Str$(r)

```

Comm.bas - ParseMemoryContents

42

```

    If (Temp Then InterpretScoreData DataStruct, sTemp, CInt(Temp) 'parse out the scores and place in global structure
    End If

    'Get Error Flags
    sTemp = Mid$(sAllData, DATA_BEGIN_ERROR_FLAGS, 2)
    If Len(sTemp) > 0 Then InterpretErrorFlags DataStruct, Val(sTemp) 'parse out the flags and save in global structure

    'Get Medication remaining in device
    sTemp = "&H0" + Mid$(sAllData, DATA_BEGIN_MED_REMAINING + 2, 2) + Mid$(sAllData, DATA_BEGIN_MED_REMAINING, 2)
    DataStruct.sMedRemaining = CStr(CSng(sTemp / 40) * 100) + " mg" 'pump pcks are fixed at 40 per milliliter (100 mg per ml)

    'Get Dose Lockout Hours
    sTemp = "&H0" + Mid$(sAllData, DATA_BEGIN_DOSE_LOCKOUT_HOURS, 2)
    DataStruct.sDoseLockoutHours = CStr(CSng(sTemp))

    'Get Doses per day
    sTemp = "&H0" + Mid$(sAllData, DATA_BEGIN_DOSES_PER_DAY, 2)
    DataStruct.iDosesPerDay = CInt(sTemp)

    'Get Dose Resolution
    sTemp = "&H0" + Mid$(sAllData, DATA_BEGIN_DOSE_CONVERSION, 2)
    DataStruct.sDoseResolution = CStr(CSng(sTemp))

    'Get Dose Intervals
    sLastIntervalTime = "1:00"

    'Get Dose Interval 1 (alarm time)
    sTemp = "&H0" + Mid$(sAllData, DATA_BEGIN_DOSE_INTERVAL1, 2)
    If Val(sTemp) Then
        sTemp = DateAdd("H", CDBl(sTemp), sLastIntervalTime) 'the first dose is relative to 1:00 am
        sLastIntervalTime = sTemp
        DataStruct.dPrescribedDoseTime(1) = TimeValue(sTemp)
    Else
        DataStruct.dPrescribedDoseTime(1) = -1 'this value indicates that no time was received
    End If

    'Get Dose Interval 2 (alarm time)
    sTemp = "&H0" + Mid$(sAllData, DATA_BEGIN_DOSE_INTERVAL2, 2)
    If Val(sTemp) Then
        sTemp = DateAdd("H", CDBl(sTemp), sLastIntervalTime) 'the first dose is relative to 1:00 am
        sLastIntervalTime = sTemp
        DataStruct.dPrescribedDoseTime(2) = TimeValue(sTemp)
    Else
        DataStruct.dPrescribedDoseTime(2) = -1 'this value indicates that no time was received
    End If

    'Get Dose Interval 3 (alarm time)
    sTemp = "&H0" + Mid$(sAllData, DATA_BEGIN_DOSE_INTERVAL3, 2)
    If Val(sTemp) Then
        sTemp = DateAdd("H", CDBl(sTemp), sLastIntervalTime) 'the first dose is relative to 1:00 am
        sLastIntervalTime = sTemp
        DataStruct.dPrescribedDoseTime(3) = TimeValue(sTemp)
    Else
        DataStruct.dPrescribedDoseTime(3) = -1 'this value indicates that no time was received
    End If

    'Get Dose Interval 4 (alarm time)
    sTemp = "&H0" + Mid$(sAllData, DATA_BEGIN_DOSE_INTERVAL4, 2)
    If Val(sTemp) Then
        sTemp = DateAdd("H", CDBl(sTemp), sLastIntervalTime) 'the first dose is relative to 1:00 am
        sLastIntervalTime = sTemp
        DataStruct.dPrescribedDoseTime(4) = TimeValue(sTemp)
    Else
        DataStruct.dPrescribedDoseTime(4) = -1 'this value indicates that no time was received
    End If

    'Get Dose Size

```

Comm.bas - ParseMemoryContents

43

```

sTemp = "&H0" + Mid$(sAllData, DATA_BEGIN_DOSE_SIZE, 2)
If IsNumeric(sTemp) Then
    DataStruct.sDoseSize = CStr(CSng(sTemp) / 40 * 100) 'convert from mg to ml
Else
    DataStruct.sDoseSize = ""
End If

```

'There are 4 fields in the device containing 16 characters each. In the original device design, this was intended to contain 4 separate pieces of information. The client has now decided that some fields are too short and others are too long. Thus, the fields are combined to be one string of 64 characters.'

```

'Patient name = giLEN_PATIENT_NAME
ID = giLEN_ID
Drug = giLEN_DRUG
TX Center = giLEN_TX_CENTER
Organ = giLEN_ORGAN

```

'Send a message to the user if the revision level is higher than the one this software is using to send custom data to the device. The user must upgrade to the current version in order to get accurate custom data. This code should also handle any previous versions that saved data to the device.'

```

'Get Custom string 1
sTemp = Mid$(sAllData, DATA_BEGIN_CUSTOM1, 32)
ConvertHexStringToAscii sTemp, sConverted
sCustomData = sConverted

```

```

'Get Custom string 2
sTemp = Mid$(sAllData, DATA_BEGIN_CUSTOM2, 32)
ConvertHexStringToAscii sTemp, sConverted
sCustomData = sCustomData + sConverted

```

```

'Get Custom string 3
sTemp = Mid$(sAllData, DATA_BEGIN_CUSTOM3, 32)
ConvertHexStringToAscii sTemp, sConverted
sCustomData = sCustomData + sConverted

```

```

'Get Custom string 4
sTemp = Mid$(sAllData, DATA_BEGIN_CUSTOM4, 32)
ConvertHexStringToAscii sTemp, sConverted
sCustomData = sCustomData + sConverted

```

*'Pull apart the 64 char string into its sub-components
'Get the custom data structure revision level that was previously saved to the device.
'Note: this is not the same as the major and minor versions of the host software.'*

```

iStartingLocation = 1
sTemp = Mid$(sAllData, iStartingLocation, giLEN_REV_DATA_STRUCTURE)
ConvertHexStringToAscii sTemp, sConverted
'The device custom data was apparently saved with a newer version of software than this one.
If Val(sConverted) > gsREV_DATA_STRUCTURE Then
    iErrorCode = ERR_NEWER_HOST_SOFTWARE
    GoTo ParseMemoryContents_Exit
End If

```

'Determine the real name of the Drug by the reference number received from the device'

```

iStartingLocation = iStartingLocation + giLEN_REV_DATA_STRUCTURE
sTemp = Trim(Mid$(sCustomData, iStartingLocation, giLEN_DRUG))
r = Val(sTemp)
If r > 0 And r < UBound(gsDrugNames) Then DataStruct.sDrug = gsDrugNames(r)

```

'Determine the real name of the Organ by the reference number received from the device'

```

iStartingLocation = iStartingLocation + giLEN_DRUG
sTemp = Trim(Mid$(sCustomData, iStartingLocation, giLEN_ORGAN))
r = Val(sTemp)
If r > 0 And r < UBound(gsOrganNames) Then DataStruct.sOrgan = gsOrganNames(r)

```

```

iStartingLocation = iStartingLocation + giLEN_ORGAN

```


Comm.bas - ParseMemoryContents

44

```

DataStruct.sPatientID = Trim(Mid$(sCustomData, IStartingLocation, giLEN_ID))

IStartingLocation = IStartingLocation + giLEN_ID
DataStruct.sTxCenter = Trim(Mid$(sCustomData, IStartingLocation, giLEN_TX_CENTER))

IStartingLocation = IStartingLocation + giLEN_TX_CENTER
r = ParseDelimString(Trim(Mid$(sCustomData, IStartingLocation, giLEN_PATIENT_NAME)), ".", sTempList())
DataStruct.sPatientLastName = Trim$(sTempList(1))
DataStruct.sPatientFirstName = Trim$(sTempList(2))

'Get Serial Number
sTemp = Mid$(sAllData, DATA_BEGIN_SERIAL_NUMBER, 20)
ConvertHexStringToAscii sTemp, sConverted
DataStruct.sSerialNumber = Trim(sConverted)
ParseMemoryContents = True 'send success to caller

ParseMemoryContents_Exit:
On Error GoTo 0
Exit Function

ParseMemoryContents_Error:
'force a checksum error here because any type of error is likely due to a checksum problem
IErrorCode = ERR_COMM_CHECKSUM
Resume ParseMemoryContents_Exit

End Function

```

Public Sub PollDeviceContinually(SourceForm As Form)

This procedure continues to try and establish communication with the Device until it succeeds. When successful, control is returned to the calling procedure. The purpose of this procedure is to allow the user to try cable changes, device movement, etc. without having to continue pressing keys on the keyboard.

Dim r As Integer, bProcedureInProgress As Boolean, IErrorCode As Long

```

If bProcedureInProgress Then Exit Sub
If gbCommBusy Or gbCommReplyPending Then Exit Sub
bProcedureInProgress = True 'prevent recursive calls to this procedure

```

QueryDevice:

```

DoEvents 'allow other Windows events to be processed, so we don't lock up the computer
If gbCommOK = True Then 'no need to poll as often if device was working the last time we checked
    Wait 5 'wait an additional amount of time before trying
Else
    Wait 0.05 'poll faster until a good comm is made
End If

```

```

If Not gbCommBusy And Not gbCommReplyPending Then 'poll only if port not busy

```

```

    If gbKeepPollingDevice = False Then Exit Sub
    SourceForm.imgPolling.Visible = True
    SourceForm.imgPolling.Refresh
    r = Comm_CheckComm(IErrorCode)
    If gbKeepPollingDevice = False Then Exit Sub
    If r = True Then 'comm is working
        'display status has not yet been updated
        If Not gbCommOK Then DisplayCommOk SourceForm

```

```

    Else 'comm is NOT working
        'display status has not yet been updated
        If gbCommOK Then DisplayCommError SourceForm

```

```

    End If
End If

```

```

Wait 0.05 'allow polling icon to be viewed
If gbKeepPollingDevice = False Then Exit Sub

```

Comm.bas - PollDeviceContinually

45

SourceFormingPolling.Visible = False

If r = ERR_COMM_TIMEOUT Then 'send an additional error message that no reply was received
MsgBox "No response was received to a wake up command that was sent to the DosPro device. ", "Communication Error"
End If

'flag has not been reset yet

If gbKeepPollingDevice Then GoTo QueryDevice 'try comm again

bProcedureInProgress = False 'allow future calls to this procedure now that we are finished

End Sub

Public Sub Comm_SendDataToDevice(ByVal sOut As String)

Dim dGoAheadTime As Double

'A data string should have been assembled by another procedure and is now
'ready to sent to the device.

'If another command is in progress then wait till it is done

'If the flags have not been reset after this delay, then exit loop anyway

'This prevents lockups inside this loop in case there is a problem elsewhere

dGoAheadTime = DateAdd("s", 5, CDBl(Now))

Do While gbCommBusy Or gbCommReplyPending

DoEvents

DoEvents

DoEvents

If CDBl(Now) > dGoAheadTime Then Exit Do

Loop

gbCommBusy = True

'set busy flag (gets reset if timeout or reply not received)

'If comm port is not open then open it

If Not frmMain.CommDevice.PortOpen Then frmMain.CommDevice.PortOpen = True

frmMain.CommDevice.InputLen = 0 'clear input buffer

frmMain.CommDevice.Output = sOut 'send string to device

End Sub

Private Sub SetCommTimer(ITime As Integer)

'The comm timer determines whether or not a reply has come back from the device.

'The timer fires an event if the ITime has passed without the timer being reset.

'Reset the timer to the interval passed in, then start it.

'Set the comm busy flag, then return to caller

frmMain.CommTimer.Enabled = False

frmMain.CommTimer.Interval = ITime

gbCommTimerExpired = False

frmMain.CommTimer.Enabled = True

'disable timer while resetting it

'set interval

'reset timer expiration flag

'start timer

End Sub

Comm.bas - SetCommTimef

46

Attribute VB_Name = "modPrinting"
Option Explicit

Public gbPrintFormLoading As Integer
Public gbPrinterErrorDetected As Integer
Public giTotalPrintPages As Integer
Public gbPreventPreviewUpdates As Integer *'track the number of pages being previewed*
Public giPrintedPageNumber As Integer
Public gbPageNumberSuspend As Integer *'if true, don't print page number for the active page (probably cover picture)*
Public giFontOptSel As Integer
Public gsSelectPrintPages As String
Public gbPrinterErrorReceived As Integer *'tells other procs that error occurred. Proc must reset flag*
Public gbPrintSpoolingInProgress As Integer *'prevent crashes during spooling*

Private Sub btnPrinter_Preview_Click_Proc()

if giTotalPrintPages > 1 Then *'there is more than one page to print*
frmSelectPages.Show MODAL
Select Case gsSelectPrintPages
Case "All"
frmPrint.MousePointer = vbHourglass
gbPrintSpoolingInProgress = True
frmPrint.vsPrinter1.Action = paPrintAll *'print all pages*
gbPrintSpoolingInProgress = False

Case "Page"
frmPrint.MousePointer = vbHourglass
gbPrintSpoolingInProgress = True
frmPrint.vsPrinter1.Action = paPrintPage *'print current page only*
gbPrintSpoolingInProgress = False

Case ""
'nothing to do
gbPrintSpoolingInProgress = False
End Select

Else *'printing a single page that is not a picture*
frmPrint.MousePointer = vbHourglass
gbPrintSpoolingInProgress = True
frmPrint.vsPrinter1.Action = paPrintAll *'print all pages*
gbPrintSpoolingInProgress = False

End If

Printing.bas - binPrinter_Preview_Clit .oc

47

```

frmPrint.btnClose.Enabled = True
frmPrint.btnPrintNow.Enabled = True
frmPrint.MousePointer = vbDefault
DoEvents

```

**allow button to show*

End Sub

Private Sub DrawHorizontalLine(cPrinter As Control, lPenColor As Long)

```

    frmPrint.vsPrinter1.FontSize = ogBodyTextSize
    frmPrint.vsPrinter1 = ""
    Exit Sub

```

**skip a line from above*

*Draw a horizontal divider on the page
Usually divides the header or topic from the rest of the page*

```

cPrinter.FontSize = 12
cPrinter = ""
cPrinter.PenStyle = 0
cPrinter.PenWidth = 10
cPrinter.PenColor = lPenColor
cPrinter.BrushColor = lPenColor

```

**skip a line from above*

**0=solid 2=dot*

**set pen width*

**set pen color*

```

Print line only across a portion of page
cPrinter.X1 = (cPrinter.PageWidth / 2) - (cPrinter.PageWidth * 0.25)
cPrinter.X2 = (cPrinter.PageWidth / 2) + (cPrinter.PageWidth * 0.25)

```

```

cPrinter.Y1 = cPrinter.CurrentY
cPrinter.Y2 = cPrinter.CurrentY + 50

```

```

cPrinter.Draw = 2

```

**1=line, 2=rectangle*

End Sub

Private Sub PrintAllPatientsSummary()

```

Dim l As Integer, lErrorCode As Long, sTableFormat As String, sTable As String, sList As String
Dim fFontSize As Single, lCount As Integer

```

On Error Resume Next

Prepare progress gauge

With frmPrint

.pnlProgress.FloodPercent = 0

.pnlProgressContainer.Visible = True

.pnlProgressContainer.Refresh

End With

InitPageProperties

fFontSize = 10

Print the logo on the first page

With frmPrint.vsPrinter1

.X1 = 700

.X2 = frmPrint.vsPrinter1.X1 + frmPrint.picLogo.Width

.Y1 = 500

.Y2 = frmPrint.vsPrinter1.Y1 + frmPrint.picLogo.Height

.Picture = frmPrint.picLogo.Picture

End With

Print information header

gbPageNumberSuspend = False

With frmPrint.vsPrinter1

64

Printing.bas - PrintAllPatientsSummary

48

```

.FontName = "Arial"
.FontBold = True
.TextAlign = taCenterTop 'center text used in paragraphs
.CurrentY = 1440 * 1 'print name on this line
.FontSize = fFontSize * 1.6 'set font size
.FontItalic = True
frmPrint.vsPrinter1 = "All Patient's Summary" 'print name
.FontSize = fFontSize * 1.2 'set font size
.FontItalic = False
DrawHorizontalLine frmPrint.vsPrinter1, &H40000 'Draw a color line
frmPrint.vsPrinter1 = "" 'skip a line
frmPrint.vsPrinter1 = frmAllPatients.cmbaDataToView.Text + " with " + frmAllPatients.Label1(0).Caption
frmPrint.vsPrinter1 = "Date Range: " + frmAllPatients.cmboDateSelection.Text + " from " + frmAllPatients.txtStartDate.Value + " to "
+ frmAllPatients.txtEndDate.Value

frmPrint.vsPrinter1 = "" 'skip a line
.TextAlign = taCenterTop
.FontBold = False
.TableBorder = tbNone
.FontSize = fFontSize * 1.1 'set font size
.Table = sTable 'send out table

frmPrint.vsPrinter1 = "" 'skip a line
frmPrint.vsPrinter1 = "" 'skip a line
.LineSpacing = 90 '% of current font
.TextAlign = taCenterTop 'center text
End With

'Print the report Data
sTableFormat = "<2400|<1800|>1700|>1700|>1000;"
'Get the column titles from grid
With frmAllPatients.grid
.Row = 0
.Col = 0
sList = .Text
.Col = 1
sList = sList + "|" + .Text
.Col = 2
sList = sList + "|" + .Text
.Col = 3
sList = sList + "|" + .Text
.Col = 4
sList = sList + "|" + .Text
End With
sTable = sTableFormat + sList

With frmPrint.vsPrinter1
.TableBorder = tbBottom
.FontSize = 10
.FontBold = True
.Table = sTable 'send out header
frmPrint.vsPrinter1 = ""
End With

'Print the information from the grid control
With frmAllPatients.grid
iCount = .Rows - 1
sList = ""
For i = 1 To iCount 'number of patients in grid
.Row = i
.Col = 0
sList = sList + .Text
.Col = 1
sList = sList + "|" + .Text
.Col = 2
sList = sList + "|" + .Text
.Col = 3
sList = sList + "|" + .Text

```

Printing.bas - PrintAllPatientsSummr

49

```

.Col = 4
sList = sList + " " + .Text + " "

frmPrint.pnlProgress.FloodPercent = (1 / iCount) * 100
frmPrint.pnlProgressContainer.Refresh
Next i
End With

sTable = sTableFormat + sList
With frmPrint.vsPrinter1
.FontSize = iFontSize * 0.9 'set font size
.LineSpacing = 80 ' % of current font
.TextAlign = taCenterTop
.TableBorder = tbNone
.Table = sTable 'send out table
End With

On Error GoTo 0
frmPrint.pnlProgressContainer.Visible = False 'turn off progress indicator
frmPrint.pnlProgressContainer.Refresh

End Sub

```

Private Sub PrintPatientDosingReport()

```

Dim i As Integer, iErrorCode As Long, sTableFormat As String, sTable As String, sList As String
Dim iFontSize As Single, iCount As Integer, bItemChecked As Boolean

On Error Resume Next

'Print Cover Art to the Print preview control if needed and available
' If frmPrint.lbPictures.ListIndex > 0 Then 'a cover is chosen
' If FileExists("covers") + sgCurrentCoverName, iErrorCode) Then 'look for a bitmap on disk
' Print preview the picture
' gbPageNumberSuspend = True
' LoadPictureToPrinterControl True 'get cover
' InitPageProperties
' frmPrint.vsPrinter1.Action = 4 'start a new page
' gbPageNumberSuspend = False
' End If
' End If

'Prepare progress gauge
With frmPrint
.pnlProgress.FloodPercent = 0
.pnlProgressContainer.Visible = True
.pnlProgressContainer.Refresh
End With
InitPageProperties
iFontSize = 10

'Print the logo on the first page
With frmPrint.vsPrinter1
.X1 = 700
.X2 = frmPrint.vsPrinter1.X1 + frmPrint.picLogo.Width
.Y1 = 500
.Y2 = frmPrint.vsPrinter1.Y1 + frmPrint.picLogo.Height
.Picture = frmPrint.picLogo.Picture
End With

'Print Information Header
gbPageNumberSuspend = False
With frmPrint.vsPrinter1
.FontName = "Arial"

```

Printing.bas - PrintPatientDosingRe

50

```

.FontSize = fFontSize * 1.6 'set font size
.FontBold = True
.FontItalic = True
.TextColor =
.TextAlign = taCenterTop 'center text used in paragraphs
.CurrentY = 1440 * 1 'print name on this line
frmPrint.vsPrinter1 = "Patient Dosing Report" 'print name
.FontSize = fFontSize * 1.2 'set font size
.FontItalic = False
frmPrint.vsPrinter1 = PAT_DATA.sPatientLastName + ", " + PAT_DATA.sPatientFirstName
DrawHorizontalLine frmPrint.vsPrinter1, 8H40000 'Draw a color line
End With

```

```

sTableFormat = "<1400><2800><1400><2800>"
sTable = sTableFormat + gsCustomLblPatientID + ";" + PAT_DATA.sPatientID + ";"
sTable = sTable + gsCustomLblTxCenter + ";" + PAT_DATA.sTxCenter + ";"
sTable = sTable + gsCustomLblDrug + ";" + PAT_DATA.sDrug + ";"
sTable = sTable + gsCustomLblOrgan + ";" + PAT_DATA.sOrgan + ";"
With frmPrint.vsPrinter1
    frmPrint.vsPrinter1 = "" 'skip a line
    .TextAlign = taCenterTop
    .FontBold = True
    .TableBorder = tbNone
    .FontSize = fFontSize * 1.1 'set font size
    .Table = sTable 'send out table

    frmPrint.vsPrinter1 = "" 'skip a line
    frmPrint.vsPrinter1 = "" 'skip a line
    .LineSpacing = 90 '% of current font
    .TextAlign = taCenterTop 'center text
End With

```

```

sList = "Events Types Shown: "
If frmPatientDosingReport.chkDoses.Value Then
    sList = sList + "Doses Taken"
    bItemChecked = True
End If

```

```

If frmPatientDosingReport.chkDoseChanged Then
    If bItemChecked Then sList = sList + " and "
    sList = sList + "Dose Size Changes"
    bItemChecked = True

```

```

If frmPatientDosingReport.chkUserDefined Then
    If bItemChecked Then sList = sList + " and "
    sList = sList + "User Entries"
    bItemChecked = True
End If

```

```

If Not bItemChecked Then
    sList = sList + "None"
End If
frmPrint.vsPrinter1 = "" 'skip a line
.TextAlign = taCenterTop
.FontSize = fFontSize * 0.9 'set font size
frmPrint.vsPrinter1 = sList
End If
End With

```

```

Print the report Date
frmPrint.vsPrinter1 = "" 'skip a line
frmPrint.vsPrinter1.TextAlign = taCenterTop
PrintDosingEventsHeader sTableFormat
Print the information from the grid control
With frmPatientDosingReport.grid
    iCount = .Rows - 1
    sList = ""
    For i = 1 To iCount 'number of patients in grid

```

Printing.bas - PrintPatientDosingRt.

51

```

        .Row = 1
        .Col = 0
        sList = sList + .Text
        .Col = 1
        sList = sList + " " + .Text
        .Col = 2
        sList = sList + " " + .Text
        .Col = 3
        sList = sList + " " + .Text
        .Col = 4
        sList = sList + " " + .Text
        .Col = 5
        sList = sList + " " + .Text + "-"

        frmPrint.pnlProgress.FloodPercent = (1 / iCount) * 100
        frmPrint.pnlProgressContainer.Refresh
    Next i
End With

sTable = sTableFormat + sList
With frmPrint.vsPrinter1
    .FontSize = fFontSize * 0.9 'set font size
    .LineSpacing = 80 ' % of current font
    .TextAlign = taCenterTop
    .Table = sTable 'send out table
End With

On Error GoTo 0
frmPrint.pnlProgressContainer.Visible = False 'turn off progress indicator
frmPrint.pnlProgressContainer.Refresh

End Sub

```

Sub PrintDosingEventsHeader(sTableFormat As String)

```

    Dim sTable As String
    Dim fPrevFont As Single, bPrevBold As Boolean, sList As String

    fPrevFont = frmPrint.vsPrinter1.FontSize

    'Print the information from the grid control
    'Pass table format back to caller
    sTableFormat = "<2100|<1600|<1000|<1700|<1700|<1700;"
    With frmPatientDosingReport.grid
        sList = ""
        .Row = 0
        .Col = 0
        sList = sList + .Text
        .Col = 1
        sList = sList + " " + .Text
        .Col = 2
        sList = sList + " " + .Text
        .Col = 3
        sList = sList + " " + .Text
        .Col = 4
        sList = sList + " " + .Text
        .Col = 5
        sList = sList + " " + .Text + "-"
    End With

    sTable = sTableFormat + sList
    frmPrint.vsPrinter1.TableBorder = tbBottom
    frmPrint.vsPrinter1.FontSize = 10
    frmPrint.vsPrinter1.FontBold = True
    frmPrint.vsPrinter1.Table = sTable 'send out header
    frmPrint.vsPrinter1 = ""

```


Printing.bas - PrintDosingEventsHeader

52

```

'Put setting back to previous ones
frmPrint.vsPrinter1.TableBorder = tbNone
frmPrint.vsPrinter1.FontSize = fPrevFont
frmPrint.vsPrinter1.FontBold = bPrevBold

```

End Sub

Public Sub LoadPictureToPrinterControl(ByVal bCover)

```

'Set the printer control to size a picture and copy
'picture from holding area to the print preview control.
'If the picture to be displayed is a cover.
'then the bCover flag should be set to true by caller.
'otherwise it is assumed to be a border.

Dim iPaperWidth%, iPaperHeight%, iNonPrintWidth%, iNonPrintHeight%
frmPrint.vsPrinter1.PhysicalPage = True 'set physical page to paper dimension
iPaperWidth = frmPrint.vsPrinter1.PageWidth 'etermine size of paper
iPaperHeight = frmPrint.vsPrinter1.PageHeight
frmPrint.vsPrinter1.PhysicalPage = False 'return printer to printable area
iNonPrintWidth = (iPaperWidth - frmPrint.vsPrinter1.PageWidth) / 2
iNonPrintHeight = (iPaperHeight - frmPrint.vsPrinter1.PageHeight) / 2

If iNonPrintWidth < 350 Then iNonPrintWidth = 350 'make a minimum margin
If iNonPrintHeight < 350 Then iNonPrintHeight = 350 'make a minimum margin

frmPrint.vsPrinter1.X1 = iNonPrintWidth
frmPrint.vsPrinter1.X2 = frmPrint.vsPrinter1.PageWidth - iNonPrintWidth
frmPrint.vsPrinter1.Y1 = iNonPrintHeight
frmPrint.vsPrinter1.Y2 = frmPrint.vsPrinter1.PageHeight - iNonPrintWidth
' frmPrint.vsPrinter1.Draw = 2 'picture holder only

frmPrint.vsPrinter1.Picture = LoadPicture("graphics\" & "deco.wmf")

```

End Sub

Private Sub InitPageMargins()

```

'Set margins
'Margins don't seem to set properly until the next page is created.
'That's why they can be set only once before printing begins.

frmPrint.vsPrinter1.MarginTop = 1350 'top margin
frmPrint.vsPrinter1.MarginBottom = 1500 'bottom margin

frmPrint.vsPrinter1.MarginLeft = 1725 'left margin
frmPrint.vsPrinter1.MarginRight = 1700 'right margin (from right edge)

```

End Sub

Printing.bas - InitPageProperties

53

Private Sub InitPageProperties()

```

'Reset margins for text and initialize other items
frmPrint.vsPrinter1.LineSpacing = 100      '100% of current font
'Set the normal attributes here
.
frmPrint.vsPrinter1.TextAlign = 0          'set centering back to normal

```

End Sub

Private Sub PrintPageDate()

```

'Print Date
Dim lTextHeight As Long, lTextWidth As Long, sText$
'print date for the above tabs only
'Rather than using .TextAlign property, text is centered here using this method
'to ensure page centering regardless of margins or paragraph settings
InitPageProperties
frmPrint.vsPrinter1.FontName = "Arial"
frmPrint.vsPrinter1.FontSize = 8
sText = "Printed: " + Date$ + " with " + App.Title + " software."
frmPrint.vsPrinter1.Measure = sText      'set string to measure
lTextHeight = frmPrint.vsPrinter1.TextHeight 'get text height
lTextWidth = frmPrint.vsPrinter1.TextWidth  'get text width
frmPrint.vsPrinter1.CurrentX = (frmPrint.vsPrinter1.PageWidth - lTextWidth) / 2

If frmPrint.vsPrinter1.CurrentY < 13000 Then
    frmPrint.vsPrinter1.CurrentY = frmPrint.vsPrinter1.PageHeight - (frmPrint.vsPrinter1.MarginBottom + (2.5 * lTextHeight)) 'set line to very bottom
Else
    frmPrint.vsPrinter1.CurrentY = frmPrint.vsPrinter1.PageHeight - (frmPrint.vsPrinter1.MarginBottom + (0.1 * lTextHeight)) 'set line to very bottom
End If

frmPrint.vsPrinter1 = sText

sText = "Copyright 1998 by SangStat Medical Corporation"
frmPrint.vsPrinter1.Measure = sText      'set string to measure
lTextWidth = frmPrint.vsPrinter1.TextWidth 'get text width
frmPrint.vsPrinter1.CurrentX = (frmPrint.vsPrinter1.PageWidth - lTextWidth) / 2
frmPrint.vsPrinter1 = sText

```

End Sub

Public Sub PrintPageNumber()

```

'Print page number if check box is active on form
giPrintedPageNumber = giPrintedPageNumber + 1      'increment page number for next time
If gbPageNumberSuspend = False Then
    frmPrint.vsPrinter1.HdrFontSize = 8
    frmPrint.vsPrinter1.Footer = "Dosing Report " + PAT_DATA.sPatientLastName + ", " + PAT_DATA.sPatientFirstName + " " + PAT_DATA.sPatientID + " Page " + CStr(giPrintedPageNumber)
Else
    frmPrint.vsPrinter1.Footer = "" 'must print a blank footer otherwise old page # will show
End If

```

End Sub

Printing.bas - RefreshPreview

54

Public Sub RefreshPreview()

```

Static bRefreshPreviewInProgress As Integer
'prevent recursive calls to here
If bRefreshPreviewInProgress = True Then Exit Sub
If gbPreventPreviewUpdates Then Exit Sub
bRefreshPreviewInProgress = True

frmPrint.MousePointer = vbHourglass
frmPrint.HScroll1.Enabled = False
frmPrint.HScroll1.Refresh
frmPrint.HScroll1.Value = 1
DoEvents
On Error GoTo 0          'reset error processing

frmPrint.btnRefresh.Enabled = False
frmPrint.btnRefresh.Refresh

frmPrint.btnPrintNow.Enabled = False
frmPrint.btnPrintNow.Refresh

frmPrint.btnClose.Enabled = False          'disable buttons until preview build is complete
frmPrint.btnClose.Refresh

' frmPrint.btnFormat.Enabled = False
' frmPrint.btnFormat.Refresh
DoEvents

glTotalPrintPages = 0          'reset the page counter
glPrintedPageNumber = 1

frmPrint.vsPrinter1.Preview = True          'print to screen
frmPrint.vsPrinter1.Footer = ""          'Clear the footer
frmPrint.vsPrinter1.HdrFontName = "font name goes here" 'Controls footer also
frmPrint.vsPrinter1.HdrFontSize = ?? 'Controls footer also

'Send information to the preview screen
Initialize print job
InitPageMargins

If gbPrinterErrorDetected Then GoTo RefreshPreview_Exit

frmPrint.vsPrinter1.PreviewPage = 1          'show 1st page
frmPrint.vsPrinter1.PreviewMode = 0          '0=screen compatible, 1=print compat, 2 = force monochrome
frmPrint.vsPrinter1.PageBorder = 0          'no page border
frmPrint.vsPrinter1.TextAlign = 0          'left align text

' Call LoadPictureToPrinterControl(False)
Select Case gaActiveFormName
    Case "frmPatientSummary"
    Case "frmAllPatients"
        Call PrintAllPatientsSummary
    Case "frmPatientDosingReport"
        Call PrintPatientDosingReport
End Select

PrintPageDate          'print date for last recipe

frmPrint.vsPrinter1.Action = paEndDoc          'END DOC
frmPrint.vsPrinter1.Visible = True
'vb5 says object does not support this method frmPrint.vsPrinter1.Refresh
Call UpdatePageButtons

```

Printing.bas - RefreshPreview

55

```

frmPrint.HScroll1.Max = gTotalPrintPages

RefreshPreview_Exit:
frmPrint.btnClose.Enabled = True      'enable buttons
DoEvents
bRefreshPreviewInProgress = False      'allow future calls to this procedure

frmPrint.btnPrintNow.Enabled = True
frmPrint.MousePointer = vbDefault
DoEvents

End Sub

Public Sub SetPreviewSize()
Dim bHeightLimit%, fTemp As Single

frmPrint.MousePointer = vbHourglass
' frmPrint.Refresh 'a refresh of the form causes controls inside a frame to disappear
' frmPrint.vsPrinter1.Visible = False
' frmPrint.vsViewPort1.Visible = False
DoEvents

If (frmPrint.vsPrinter1.PageHeight / frmPrint.vsPrinter1.PageWidth) > (frmPrint.vsViewPort1.Height / frmPrint.vsViewPort1.Width) Then
bHeightLimit = True

Select Case frmPrint.optZoom(0).Value
Case True      'full page view
If bHeightLimit = True Then 'there is a height restriction in the viewport control for this print orientation
frmPrint.vsPrinter1.Height = frmPrint.vsViewPort1.Height * 0.99
fTemp = frmPrint.vsPrinter1.PageWidth / frmPrint.vsPrinter1.PageHeight
fTemp = frmPrint.vsPrinter1.Height * fTemp
frmPrint.vsPrinter1.Width = fTemp
Else
frmPrint.vsPrinter1.Width = frmPrint.vsViewPort1.Width * 0.99
fTemp = frmPrint.vsPrinter1.PageHeight / frmPrint.vsPrinter1.PageWidth
fTemp = frmPrint.vsPrinter1.Width * fTemp
frmPrint.vsPrinter1.Height = fTemp
End If

'Make viewport virtual screen large enough to show full page of print control
frmPrint.vsViewPort1.VirtualWidth = frmPrint.vsPrinter1.Width * 1
frmPrint.vsViewPort1.VirtualHeight = frmPrint.vsPrinter1.Height * 1
frmPrint.vsViewPort1.BorderStyle = 1 'turn off border

Case Else      'Magnify view
frmPrint.vsPrinter1.Width = frmPrint.vsPrinter1.PageWidth * 1
frmPrint.vsPrinter1.Height = frmPrint.vsPrinter1.PageHeight * 1
frmPrint.vsViewPort1.VirtualWidth = frmPrint.vsPrinter1.Width * 1
frmPrint.vsViewPort1.VirtualHeight = frmPrint.vsPrinter1.Height * 1 'ensure scroll bars will be shown
frmPrint.vsViewPort1.BorderStyle = 0 'turn on border
End Select

frmPrint.vsPrinter1.Visible = True
frmPrint.vsViewPort1.Visible = True
frmPrint.vsViewPort1.Refresh
frmPrint.MousePointer = vbDefault
DoEvents

End Sub

```

Printing.bas - UpdatePageButtons

56

Public Sub UpdatePageButtons()

`frmPrint.LibPageNumber.Caption = "Page " + CStr(frmPrint.HScroll1.Value) + " of " + CStr(giTotalPrintPages)``frmPrint.LibPageNumber.Refresh``If giTotalPrintPages < 2 Then``frmPrint.HScroll1.Enabled = False` *'no scroll bar needed for a single page*`Else``frmPrint.HScroll1.Enabled = True``End If``DoEvents``End Sub`

Fax.bas - File Declarations

57

```

Attribute VB_Name = "modFax"
Option Explicit
Public gcFax As Control
Public gsFaxFileSpec As String
Public gsEditName As String 'a temporary place to hold fax names being edited or created
Public gsEditVoice As String
Public gsEditFax As String

Public gsEditGroupIndexes As String 'holds temporary indexes to all locations associated with a group
Public gsEditGroupName As String

Type FaxDataStructure
    sFaxID As String
    sDialPrefix As String
    iRetries As Integer
    iRetryInterval As Integer
    bFaxResolution As Byte
    sSenderName As String
    sSenderCompany As String
    sSenderFaxNumber As String
    sSenderVoiceNumber As String

    iLocTotal As Integer 'a count of the locations
    sLocPersonName(100) As String 'rgh it may be desirable in the future to make these arrays dynamic
    sLocFaxNumber(100) As String
    sLocVoiceNumber(100) As String

    iGroupsTotal As Integer
    sGroupTitle(50) As String
    sGroupNamesInTitle(50) As String 'indexes to names separated by pipe. (ie J[6]15)

    iGroupLastSelected As Integer
End Type
Public FAX_DATA As FaxDataStructure

```

Public Sub GetFaxLocations()

```
Dim i As Integer, r As Integer, sSection As String
```

```
With FAX_DATA
```

```
    sSection = "Fax Locations"
```

```
    iLocTotal = CInt(GetINISetting(gsFaxFileSpec, sSection, "Total Locations", "0"))
```

```
    For i = 1 To iLocTotal
```

```
        .sLocPersonName(i) = GetINISetting(gsFaxFileSpec, sSection, "Person " + CStr(i), "")
```

```
        .sLocFaxNumber(i) = GetINISetting(gsFaxFileSpec, sSection, "Fax " + CStr(i), "")
```

```
        .sLocVoiceNumber(i) = GetINISetting(gsFaxFileSpec, sSection, "Voice " + CStr(i), "")
```

```
    Next i
```

```
    sSection = "Fax Groups"
```

```
    iGroupsTotal = GetINISetting(gsFaxFileSpec, sSection, "Total Groups", "0")
```

```
    For i = 0 To iGroupsTotal
```

```
        .sGroupTitle(i) = GetINISetting(gsFaxFileSpec, sSection, "Group " + CStr(i), "")
```

```
        .sGroupNamesInTitle(i) = GetINISetting(gsFaxFileSpec, sSection, "Group Locations " + CStr(i), "")
```

```
    Next i
```

```
    sSection = "User Selections"
```

```
    iGroupLastSelected = CInt(GetINISetting(gsFaxFileSpec, sSection, "Last Group Selected", "0"))
```

```
End With
```

```
End Sub
```

Fax.bas - GetIndexToFaxGroupName

58

Public Function GetIndexToFaxGroupName(ByVal sGroup As String) As Integer

*'Find sName in the list of fax names. If found, pass index back to caller.
'otherwise return 0.*

```
Dim I As Integer
sGroup = LCase$(sGroup)
With FAX_DATA
    For I = 1 To JGroupsTotal
        If LCase$(sGroupTitle(I)) = sGroup Then
            GetIndexToFaxGroupName = I
            Exit Function
        End If
    Next I
End With
```

End Function

Public Function GetIndexToFaxLocName(ByVal sName As String) As Integer

*'Find sName in the list of fax names. If found, pass index back to caller.
'otherwise return 0.*

```
Dim I As Integer
sName = LCase$(sName)
With FAX_DATA
    For I = 1 To JLocTotal
        If LCase$(sLocPersonName(I)) = sName Then
            GetIndexToFaxLocName = I
            Exit For
        End If
    Next I
End With
```

End Function

Public Sub RemoveGroupFromFaxList(ByVal sGroup As String)

*'Remove the name from the list and move up all others in the list.
Dim I As Integer, J As Integer, IIndexFound As Integer*

```
With FAX_DATA
    For I = 1 To JGroupsTotal
        If sGroupTitle(I) = sGroup Then
            IIndexFound = I
            Exit For
        End If
    Next I

    For I = IIndexFound To JGroupsTotal - 1
        sGroupTitle(I) = sGroupTitle(I + 1)
        sGroupNamesInTitle(I) = sGroupNamesInTitle(I + 1)
    Next I
```

JGroupsTotal = JGroupsTotal - 1

End With

End Sub

Fax.bas - RemoveNameFromFaxList

59

Public Sub RemoveNameFromFaxList(ByVal sName As String)*'Remove the name from the list and move up all others in the list.**Dim i As Integer, j As Integer, iIndexFound As Integer, r As Integer**Dim sTempList(100) As String, sNewIndexes As String, iTemp As Integer***With FAX_DATA***For i = 1 To iLocTotal* *Took through whole list for name**If .sLocPersonName(i) = sName Then* *'found it here**iIndexFound = i**Exit For**End If**Next i**For i = iIndexFound To iLocTotal - 1**.sLocPersonName(i) = .sLocPersonName(i + 1)**.sLocVoiceNumber(i) = .sLocVoiceNumber(i + 1)**.sLocFaxNumber(i) = .sLocFaxNumber(i + 1)**Next i**iLocTotal = iLocTotal - 1**'Now that the name has been removed, we must look at all of the indexes of**'each fax group to see if an index pointer was in there. If so, it must**'be removed. Additionally, all index greater than the one removed must be**'decremented by one.**If iIndexFound Then**For j = 1 To jGroupsTotal* *Took at each index record in a fax group**'Parse out all of the indexes into a list for easier processing**r = ParseDelimString(.sGroupNamesInTitle(j), " ", sTempList())**sNewIndexes = ""**If r Then* *Indexes were found for this record**For j = 1 To r**'look at each item in the list to see if it equals or great than the one removed**iTemp = Cint(sTempList(j))**If iTemp = iIndexFound Then* *'same index must be removed from list**'nothing to do. Don't add it to new list of indexes**ElseIf iTemp > iIndexFound Then* *'higher indexes must be decremented by one.**iTemp = iTemp - 1**sNewIndexes = sNewIndexes + CStr(iTemp) + " "**Else* *'original value is OK**sNewIndexes = sNewIndexes + CStr(iTemp) + " "**End If**Next j**End If**.sGroupNamesInTitle(j) = sNewIndexes* *'store the new list of indexes back to array**Next j**End If**End With***End Sub**

Fax.bas - SetFaxDeviceLabel

60

Public Sub SetFaxDeviceLabel()

*This label on the options tab displays the status of the fax device.
If a fax device exists, then the label displays the device, otherwise
it shows an appropriate message.*

```
With frmOptions.lblFaxDevice
    If gcFax.DeviceCount > 0 Then      'at least one fax device was found
        For i = 0 To gcFax.DeviceCount - 1
            .Caption = gcFax.Devices(i) 'show name of the device found
            .BackColor = &HFF00&      'green background
            .ForeColor = &H0&
        Next i
    Else                               'no fax devices were found
        .Caption = "A fax device was not found. Please ensure the fax or modem is connected properly."
        .BackColor = &H80&           'red background
        .ForeColor = &HFFFFFF        'white
    End If
End With

End Sub
```

Calendar.bas - File Declarations

61

Attribute VB_Name = "modCalendar"
Option Explicit

Private giCompliedDosesCreated As Integer *'number of Complied Doses to show on the calendar*
Private giNonCompliedDosesCreated As Integer *'number of non-complied Doses to show on the calendar*
Private giDoseSizeChangesCreated As Integer
Private giZoomDosesCreated As Integer *'number of Doses to show in zoom box*
Private giDosesMissedCreated As Integer *'number of objects to show for missed days*
Private gbCalendarUpdateInProgress As Integer *'prevents recursive calls while updating calendar*

Public gsngComplianceTimeRange As Single *'% of hrs on either side of a prescribed dose in which a dose must be taken*

Type CALENDAR_SELECTIONS

chkDosesTaken As Byte
chkDosesNotComplied As Byte
chkDosesMissed As Byte
chkDoseChanged As Byte

End Type

Public CAL_DEFAULTS As CALENDAR_SELECTIONS

Type SUMMARY_SELECTIONS

cmboDataToView As Byte
cmboChartType As Byte

End Type

Public PAT_SUM_DEFAULTS As SUMMARY_SELECTIONS

Public Function CalcDaysInMonth(ByVal iMonth As Integer, ByVal iYear As Integer)

'Calculate the number of days in the month/year that is passed here
Dim i As Integer, iTemp As Long

i = iMonth + 1
If i = 13 Then i = 1
iTemp = CDate(CStr(i) + "/" + "01" + CStr(iYear))
iTemp = iTemp - 1
CalcDaysInMonth = Day(iTemp)
End Function

Public Sub DrawAllDoseSizeChanges()

Dim i As Integer, r As Integer, iDaysInMonth As Integer
Dim sCalendarStartDate As String, dTime As Double
Dim iDateDifference As Long, iCalendarStartDate As Long
Dim bFirstDayAlreadyPlotted As Boolean, bLastDayAlreadyPlotted As Boolean
RemoveDoseSizeChanges *'remove all of the old doses first*

iDaysInMonth = CalcDaysInMonth(iFromDosingCalendar.Calendar.Month, iFromDosingCalendar.Calendar.Year)
sCalendarStartDate = CStr(iFromDosingCalendar.Calendar.Month) + "/" + "01" + CStr(iFromDosingCalendar.Calendar.Year)
iCalendarStartDate = DateValue(sCalendarStartDate)

If iFromDosingCalendar.chkDoseChanged.Value Then
 For i = 1 To PAT_DATA.iEventData(0) *'total number of events*
 If PAT_DATA.byteEventType(i) = giEVENT_DOSE_CHANGED Then *'show only med events (not errors, etc)*
 iDateDifference = Int(PAT_DATA.dEventData(i)) - iCalendarStartDate
 If iDateDifference >= 0 And iDateDifference < iDaysInMonth Then
 dTime = PAT_DATA.dEventData(i) - Int(PAT_DATA.dEventData(i))
 DrawSingleDoseSizeChange CInt(iDateDifference + 1), dTime, i, True
 End If
 End If
 Next i
End If

'This section of code ensures that dosing info is always plotted on the first and last day of the month.

If bFirstDayAlreadyPlotted = False Then

Calendar.bas - DrawAllDoseSizeChanges

62

```

    r = FindPrescribedDoseSizeForSpecificDay(PAT_DATA, iCalendarStartDate)
    DrawSingleDoseSizeChange i, dTime, r, False
End If

If bLastDayAlreadyPlotted = False Then
    r = FindPrescribedDoseSizeForSpecificDay(PAT_DATA, iCalendarStartDate + iDaysInMonth - 1)
    DrawSingleDoseSizeChange iDaysInMonth, dTime, r, False
End If

For i = 1 To giDoseSizeChangesCreated      'show all the Doses
    frmDosingCalendar.shapeDoseSizeChange(i).Visible = True
Next i

End Sub

```

```

Public Sub DrawAllCompliedDosesTaken()
    Dim i As Integer, r As Integer
    Dim sCalendarStartDate As String, dTime As Double
    Dim iDateDifference As Long, iCalendarStartDate As Long
    Dim iDayDoseCount As Integer, iDayNumberBeingPlotted As Integer, iLastDoseDayDrawn As Integer
    Dim iDaysInMonth As Integer, iTemp As Long

    RemoveCompliedDosesTaken      'remove all of the old doses first
    If frmDosingCalendar.chkDosesTaken.Value = False Then Exit Sub

    sCalendarStartDate = CStr(frmDosingCalendar.Calendar.Month) + "/01/" + Str$(frmDosingCalendar.Calendar.Year)
    iCalendarStartDate = DateValue(sCalendarStartDate)

    'Calc the number of days in the month being displayed
    iDaysInMonth = CalcDaysInMonth(frmDosingCalendar.Calendar.Month, frmDosingCalendar.Calendar.Year)

    For i = 1 To PAT_DATA.EventData(0)      'total number of events
        If PAT_DATA.byteEventType(i) = giEVENT_DOSE_TAKEN Then      'show only med events (not errors, etc)
            iDateDifference = Int(PAT_DATA.dEventDate(i)) - iCalendarStartDate
            If iDateDifference >= 0 And iDateDifference < iDaysInMonth Then      'dose occurred during this month
                'Determine if the dose occurred within the compliance parameters
                r = IsDoseWithinPrescribedTimeRange(PAT_DATA, i)      'pass index to event time
                If r Then
                    iDayNumberBeingPlotted = CInt(iDateDifference + 1)
                    If iDayNumberBeingPlotted = iLastDoseDayDrawn Then
                        iDayDoseCount = iDayDoseCount + 1      'plotting same day as last dose
                    Else
                        iDayDoseCount = 1      'this is a new day. Reset counter
                    End If
                    iLastDoseDayDrawn = iDayNumberBeingPlotted      'remember that we are plotting on this day
                    DrawSingleCompliedDoseTaken iDayNumberBeingPlotted, dTime, i, iDayDoseCount
                End If
            End If
        End If
    Next i

    For i = 1 To giCompliedDosesCreated      'show all the Doses
        frmDosingCalendar.shapeDose(i).Visible = True
    Next i

End Sub

```

Calendar.bas - DrawAllDosesMissed

63

Public Sub DrawAllDosesMissed()

```

Dim i As Integer, iDaysInMonth As Integer, i As Long
Dim sCalendarStartDate As String
Dim iDateDifference As Long, iCalendarStartDate As Long
Dim iDayDoseCount As Integer, iDayBeingPlotted As Integer

```

```

RemoveDosesMissed 'remove all of the old doses first
If frmDosingCalendar.chkDosesMissed.Value = False Then Exit Sub

```

```

iDaysInMonth = CalcDaysInMonth(frmDosingCalendar.Calendar.Month, frmDosingCalendar.Calendar.Year)
sCalendarStartDate = CStr(frmDosingCalendar.Calendar.Month) + "/" + CStr(frmDosingCalendar.Calendar.Year)
iCalendarStartDate = DateValue(sCalendarStartDate)

```

```

For i = iCalendarStartDate To iCalendarStartDate + iDaysInMonth - 1 'sequence through all days in month
    If i >= PAT_DATA.dEventData(1) Then 'day being plotted is not earlier than 1st dose in structure
        If i < PAT_DATA.dEventData(PAT_DATA.iEventData(0)) Then 'day being plotted is not later than last dose in structure
            iDayBeingPlotted = i - iCalendarStartDate + 1 'get the current month day to plot
            iDayDoseCount = CalcDosesSumTakenOnSpecificDay(PAT_DATA, i) 'calc missed doses for this day
            For j = 1 To PAT_DATA.iDosesPerDay - iDayDoseCount
                'DrawSingleDoseMissed iDayBeingPlotted, j 'Plot the current day
            Next j
        End If
    End If
Next i

```

```

For i = 1 To giDosesMissedCreated 'show all the Doses
    frmDosingCalendar.shapeDoseMissed(i).Visible = True
Next i

```

End Sub

Public Sub DrawAllNonCompliedDosesTaken()

```

Dim i As Integer, j As Integer, bDoseOutOfCompliance As Boolean
Dim dTimeLimit As Double, dLowLimit As Double, dHighLimit As Double
Dim sCalendarStartDate As String, dTime As Double
Dim iDateDifference As Long, iCalendarStartDate As Long
Dim iDayDoseCount As Integer, iDayNumberBeingPlotted As Integer, iLastDoseDayDrawn As Integer
Dim iDaysInMonth As Integer

```

```

RemoveNonCompliedDosesTaken 'remove all of the old doses first
If frmDosingCalendar.chkDosesNotComplied.Value = False Then Exit Sub

```

```

iDaysInMonth = CalcDaysInMonth(frmDosingCalendar.Calendar.Month, frmDosingCalendar.Calendar.Year)
sCalendarStartDate = CStr(frmDosingCalendar.Calendar.Month) + "/" + CStr(frmDosingCalendar.Calendar.Year)
iCalendarStartDate = DateValue(sCalendarStartDate)

```

```

dTimeLimit = gsngComplianceTimeRange / 24
For i = 1 To PAT_DATA.iEventData(0) 'total number of events
    If PAT_DATA.byteEventType(i) = giEVENT_DOSE_TAKEN Then 'show only med events (not errors, etc)
        iDateDifference = Int(PAT_DATA.dEventData(i)) - iCalendarStartDate
        If iDateDifference >= 0 And iDateDifference <= iDaysInMonth Then 'dose occurred during this month
            dTime = PAT_DATA.dEventData(i) - Int(PAT_DATA.dEventData(i)) 'get time of dose
            'Determine if the dose occurred within the compliance parameters
            'gh: see if we can use our procedure already created
            If gsngComplianceTimeRange Then 'do test if there is a value set in the compliance time range
                bDoseOutOfCompliance = True 'set default to be out of range unless otherwise set below
                For j = 1 To PAT_DATA.iDosesPerDay
                    'compare dose time against all of the alarm times
                    dLowLimit = PAT_DATA.dPrescribedDoseTime(j) - dTimeLimit - 0.0001 'add a factor to prevent rounding error
                    dHighLimit = PAT_DATA.dPrescribedDoseTime(j) + dTimeLimit + 0.0001
                    If dTime >= dLowLimit And dTime <= dHighLimit Then 'this dose is within compliance
                        bDoseOutOfCompliance = False 'set flag to not plot this dose
                    End If
                Next j
            End If
        End If
    End If
Next i

```

Calendar.bas - DrawAllNonCompliedDosesT

64

```

Else 'there is no compliance range
    bDoseOutOfCompliance = False
End If

If bDoseOutOfCompliance Then
    iDayNumberBeingPlotted = CInt(iDateDifference + 1)
    If iDayNumberBeingPlotted = iLastDoseDayDrawn Then
        iDayDoseCount = iDayDoseCount + 1 'plotting same day as last dose
    Else
        iDayDoseCount = 1 'this is a new day. Reset counter
    End If
    iLastDoseDayDrawn = iDayNumberBeingPlotted 'remember that we are plotting on this day

    DrawSingleNonCompliedDoseTaken CInt(iDateDifference + 1), dTime, i, iDayDoseCount
End If
End If
End If
Next i

For i = 1 To glNonCompliedDosesCreated 'show all the Doses
    frmDosingCalendar.shapeDoseNonComply(i).Visible = True
Next i
End Sub

```

Public Sub DrawSingleDoseSizeChange(iDay As Integer, dTime As Double, iEventNumber As Integer, bHighlight

*'Draw dosing Doses for the day of the month and time of day passed in here.
 'Time of time is expressed in decimal places as a portion of a day (VB time format)
 'NOTE: No checks are currently made to determine whether the event is a med event or
 'a non-med event. If both events are kept in the same array, then a test of the med bit
 'must be done before plotting.
 'A Dose is not visible when first created. The caller should display the Doses once they
 'are all created, so as to speed the redraw of the screen.
 'Create another done of the the Dose shape located in the array Dose(0)
 'When this feature is on, a dose size is automatically entered on the first and last day of the month.*

** On Error Resume Next*

Dim i As Integer, iWeekDay As Integer

Dim iDoseLeft As Single, iDoseTop As Single

Dim iTemp As Long, iDayWidth As Long, iDayHeight As Long

glDoseSizeChangesCreated = glDoseSizeChangesCreated + 1 *'Increment counter*

Load frmDosingCalendar.shapeDoseSizeChange(glDoseSizeChangesCreated) *'create a new object*

DoEvents

frmDosingCalendar.shapeDoseSizeChange(glDoseSizeChangesCreated) = "" + CStr(PAT_DATA.IEventData(iEventNumber)) + " mg

If bHighlight Then

frmDosingCalendar.shapeDoseSizeChange(glDoseSizeChangesCreated).BackColor = &HFFFFFFC0 *'blue*

frmDosingCalendar.shapeDoseSizeChange(glDoseSizeChangesCreated).ToolTipText = " Dose Size Was Changed Today "

Else

frmDosingCalendar.shapeDoseSizeChange(glDoseSizeChangesCreated).ToolTipText = " Current Dose Size "

End If

'These lines are a work-around for a bug in the control that causes it to

'return the wrong values for day.left, day.top, etc. When fixed, we can simply use those properties

'and remove these calculations.

iDayWidth = (frmDosingCalendar.Calendar.Width - 50) / 7 *'actual scalewidth of a single day*

*iTemp = (frmDosingCalendar.Calendar.DayLeft(iDay) * 25) / iDayWidth* *'approximate location of the day*

iWeekDay = CInt(iTemp)

*iDoseLeft = (iWeekDay * iDayWidth)* *'get left edge of day to plot*

*iDoseLeft = iDoseLeft - ((iDayWidth / 5) * iPlotPosition - 1) - (iDayWidth / 10)*

*iDoseLeft = iDoseLeft + (iDayWidth * 0.8) - frmDosingCalendar.shapeDoseSizeChange(glDoseSizeChangesCreated).Width*

If Int(iDayWidth / 150) < 7 Then

frmDosingCalendar.shapeDoseSizeChange(glDoseSizeChangesCreated).FontBold = False

Else

Calendar.bas - DrawSingleDoseSizeChange

65

```
frmDosingCalendar.shapeDoseSizeChange(giDoseSizeChangesCreated).FontBold = True
End If
```

```
frmDosingCalendar.shapeDoseSizeChange(giDoseSizeChangesCreated).FontSize = Int(IDayWidth / 150)
frmDosingCalendar.shapeDoseSizeChange(giDoseSizeChangesCreated).Left = IDoseLeft
```

'These lines are a work-around for a bug in the control that causes it to return the wrong values for day left, day top, etc. When fixed, we can simply use those properties and remove these calculations.

'The control is even more stupid than I first suspected. It can not always return the proper vertical location of a day, thus, we have to jump through more hoops to figure out what week that a particular day is in.

```
IDayHeight = (frmDosingCalendar.Calendar.Height - 50) - 625 'an offset is used to compensate for height of title
IDayHeight = IDayHeight / 6
IWeekDay = (frmDosingCalendar.Calendar.DayLeft(1) * 26) / IDayWidth 'the approximate location of the day
ITemp = Int((IDay + IWeekDay - 1) / 7)
```

```
IDoseTop = (ITemp * IDayHeight) + 625 'this number factors in the title bar
IDoseTop = IDoseTop + 50
```

```
frmDosingCalendar.shapeDoseSizeChange(giDoseSizeChangesCreated).Top = IDoseTop
frmDosingCalendar.shapeDoseSizeChange(giDoseSizeChangesCreated).Tag = iEventNumber
```

'Keep event number for updating the zoom box

' On Error GoTo 0

End Sub

Private Sub DrawSingleNonCompliedDoseTaken(iDay As Integer, dTime As Double, iEventNumber As Integer, iPlotPosition As Integer)

'Draw dosing Doses for the day of the month and time of day passed in here.

'Time of time is expressed in decimal places as a portion of a day (VB time format)
NOTE: No checks are currently made to determine whether the event is a med event or a non-med event. If both events are kept in the same array, then a test of the med bit must be done before plotting.

'A Dose is not visible when first created. The caller should display the Doses once they are all created, so as to speed the redraw of the screen.

'Create another clone of the the Dose shape located in the array Dose(0)

On Error Resume Next

Dim I As Integer, IWeekDay As Integer

Dim IDoseLeft As Single, IDoseTop As Single

Dim ITemp As Long, IDayWidth As Long, IDayHeight As Long

```
giNonCompliedDosesCreated = giNonCompliedDosesCreated + 1 'increment Dose counter
Load frmDosingCalendar.shapeDoseNonComply(giNonCompliedDosesCreated) 'create a new Dose
```

'These lines are a work-around for a bug in the control that causes it to return the wrong values for day left, day top, etc. When fixed, we can simply use those properties and remove these calculations.

```
IDayWidth = (frmDosingCalendar.Calendar.Width - 50) / 7 'actual scalewidth of a single day
IWeekDay = (frmDosingCalendar.Calendar.DayLeft(iDay) * 26) / IDayWidth 'approximate location of the day
```

```
IDoseLeft = (IWeekDay * IDayWidth) 'get left edge of day to plot
IDoseLeft = IDoseLeft + ((IDayWidth / 5) * iPlotPosition - 1) - (IDayWidth / 10)
frmDosingCalendar.shapeDoseNonComply(giNonCompliedDosesCreated).Left = IDoseLeft
```

'These lines are a work-around for a bug in the control that causes it to return the wrong values for day left, day top, etc. When fixed, we can simply use those properties and remove these calculations.

'The control is even more stupid than I first suspected. It can not always return the proper vertical location of a day, thus, we have to jump through more hoops to figure out what week that a particular day is in.

```
IDayHeight = (frmDosingCalendar.Calendar.Height - 50) - 625 'an offset is used to compensate for height of title
```

Calendar.bas - DrawSingleNonCompliedDose.

66

```

IDayHeight = IDayHeight / 6
WeekDay = (frmDosingCalendar.Calendar.DayLeft(1) * 26) / IDayWidth 'the approximate location of the day
ITemp = Int((IDay + iWeekDay - 1) / 7)

IDoseTop = (ITemp * IDayHeight) + 625 'this number factors in the title bar
IDoseTop = IDoseTop + IDayHeight - 50 - frmDosingCalendar.shapeDose(0).Height - frmDosingCalendar.shapeDoseNonComply(0).Height 'draw in middle of day
frmDosingCalendar.shapeDoseNonComply(giNonCompliedDosesCreated).Top = IDoseTop

frmDosingCalendar.shapeDoseNonComply(giNonCompliedDosesCreated).Tag = iEventNumber 'keep event number for updating the zoom box.

On Error GoTo 0
End Sub

```

Public Function IsDoseWithinPrescribedTimeRange(DataStruct As DeviceDataStruct, ByVal iIndex As Integer)

'Test to see that the event at the index passed here is a medication event and that it is within the prescribed time range for a daily dose. If yes, then pass TRUE back to the caller.

Dim i As Integer, dTime As Double
Dim dTimeLimit As Double, dLowLimit As Double, dHighLimit As Double

```

dTimeLimit = gsngComplianceTimeRange / 24
dTime = DataStruct.dEventDate(iIndex) - Int(DataStruct.dEventDate(iIndex)) 'get time of dose
If gsngComplianceTimeRange Then
    For i = 1 To DataStruct.iDosesPerDay
        'compare dose time against all of the alarm times
        dLowLimit = DataStruct.dPrescribedDoseTime(i) - dTimeLimit - 0.0001
        dHighLimit = DataStruct.dPrescribedDoseTime(i) + dTimeLimit + 0.0001 'add a factor to prevent rounding error
        If dTime >= dLowLimit And dTime <= dHighLimit Then
            IsDoseWithinPrescribedTimeRange = True
            Exit For 'no need to do any further testing for this dose
        End If
    Next i
Else
    'there is no compliance range, so pass back a success flag
    IsDoseWithinPrescribedTimeRange = True
End If
End Function

```

Private Sub PrintCalendar()

' This routine is called when the user presses the print button on the calendar form

```

Dim sPrintInfo As String
Dim CRLF As String
Dim bcolorCalendar As Long
Dim bcolorpnlZoom As Long
Dim bcolorpnlTime As Long
Dim bcolorForm As Long
Dim fcolorPrescribed As Long
Dim fcolorMissed As Long
Dim fcolorWeek As Long

Dim curTop As Long
Dim curWidth As Integer
Dim curHeight As Integer

Const XOffset = 1926
Const YOffset = 1890

On Error GoTo Error_btnPrint

CRLF = Chr$(13) + Chr$(10)

```

Calendar.bas - PrintCalendar

67

```

*
* curTop = Me.Top
* curWidth = Me.Width
* curHeight = Me.Height
*
* Hide this guy off the screen while we print
* Me.Top = -(curHeight * 2)
*
* Save current background colors
* bcolorCalendar = Calendar.BackColor
* bcolorpnlZoom = pnlZoom.BackColor
* bcolorpnlTime = pnlTime.BackColor
* bcolorForm = Me.BackColor
* fcolorPrescribed = chkDosesTaken.FillColor
* fcolorMissed = chkDosesMissed.FillColor
* fcolorWeek = chkWeekNumbers.FillColor
*
* hide the buttons
* btnClose.Visible = False
* btnPrint.Visible = False
*
* Add date + time info to printed data
* sPrintInfo = "Printed on: " & Format$(Now, "ddddd hh:nn")
* lblPrintInfo.Caption = sPrintInfo
*
* Set titles at top of printed page
* lblTitle.Caption = "Dosing Calendar"
* lblPatient.Caption = "Patient: " & tgDeviceStat.sPatient
* lblDrug.Caption = "Drug: " & tgDeviceStat.sDrug
*
* Set background colors
* Calendar.BackColor = WHITE
* pnlZoom.BackColor = WHITE
* pnlTime.BackColor = WHITE
* Me.BackColor = WHITE
* chkDosesTaken.FillColor = WHITE
* chkDosesMissed.FillColor = WHITE
* chkWeekNumbers.FillColor = WHITE
*
* Let user know we are printing
* Load frm_Status
* frm_Status.lblStatus.Caption = "Preparing to print calendar"
* frm_Status.Show
*
* Move/resize the form and move objects to give space for printing
* Call DeleteAllObjects 'remove extraneous elements from calendar
* Call MoveFormObjects(Me, XOffset, YOffset, True)
* Call UpdateCalendar
* DoEvents
*
* Call UpdateZoomBox
*
* Switch on visibility of titles
* lblPrintInfo.Visible = True
* lblTitle.Visible = True
* lblPatient.Visible = True
* lblDrug.Visible = True
*
* Make Check boxes two dimensional
* chkDosesTaken.CheckBox2d = True
* chkDosesMissed.CheckBox2d = True
* chkWeekNumbers.CheckBox2d = True
*
* Bring the Zoom labels to the front
* lblZoomTime.ZOrder 0
* lblZoomTime.Visible = True
*
* Print the form

```


Calendar.bas - PrintCalendar

68

```
.  
Me.Height = Me.Height - YOffset  
Me.Width = Me.Width - XOffset  
.  
DoEvents  
Me.PrintForm  
DoEvents  
.  
frm_Status.lblStatus.Caption = "Sending calendar to printer"  
.  
' hide the titles and show the buttons  
lblPnnInfo.Visible = False  
lblTitle.Visible = False
```

Calendar.bas - PrintCalendar

69

```

    btpPatient.Visible = False
    btdrug.Visible = False

    * Move everything back
    Call DeleteAllObjects
    Call MoveFormObjects(Me, -XOffset, -YOffset, True)
    Call UpdateCalendar

    * Restore background colors
    Calendar.BackColor = bcolorCalendar
    pnlZoom.BackColor = bcolorpnlZoom
    pnlTime.BackColor = bcolorpnlTime
    Me.BackColor = bcolorForm
    chkDosesTaken.FillColor = fcolorPrescribed
    chkDosesMissed.FillColor = fcolorMissed
    chkWeekNumbers.FillColor = fcolorWeek

    * Restore buttons
    btnClose.Visible = True
    btnPrint.Visible = True

    * Set check boxes back to 3d
    * Make Check boxes two dimensional
    chkDosesTaken.CheckBox2d = False
    chkDosesMissed.CheckBox2d = False
    chkWeekNumbers.CheckBox2d = False

    * Bring box back into view
    Me.Width = curWidth
    Me.Height = curHeight
    Me.Top = curTop      ' bring form back into view

    Unload frm_Status

Exit btnPrint
Exit Sub

Error btnPrint:
Resume Exit btnPrint

```

Calendar.bas - PrintCalendar

70

End Sub

Public Sub RemoveDoseSizeChanges()

Dim i As Integer

On Error Resume Next

```
For i = 1 To giDoseSizeChangesCreated 'remove all previous Doses
```

```
    Unload frmDosingCalendar.shapeDoseSizeChange(i)
```

Next i

giDoseSizeChangesCreated = 0

On Error GoTo 0

End Sub

Public Sub RemoveDosesMissed()

Dim i As Integer

On Error Resume Next

```
For i = 1 To giDosesMissedCreated 'remove all previous objects
```

```
    Unload frmDosingCalendar.shapeDoseMissed(i)
```

Next i

giDosesMissedCreated = 0

On Error GoTo 0

End Sub

Public Sub RemoveCompliedDosesTaken()

Dim i As Integer

On Error Resume Next

```
For i = 1 To giCompliedDosesCreated 'remove all previous Doses
```

```
    Unload frmDosingCalendar.shapeDose(i)
```

Next i

giCompliedDosesCreated = 0

On Error GoTo 0

End Sub

Public Sub RemoveNonCompliedDosesTaken()

Dim i As Integer

On Error Resume Next

```
For i = 1 To giNonCompliedDosesCreated 'remove all previous Doses
```

```
    Unload frmDosingCalendar.shapeDoseNonComply(i)
```

Next i

giNonCompliedDosesCreated = 0

On Error GoTo 0

End Sub

Calendar.bas - UpdateZoomBox

71

Public Sub UpdateZoomBox()

'A different day was clicked on the calendar, so we need to plot the events for the current day into the zoom box.
'This procedure draws doses taken for a given day
NOTE: No checks are currently made to determine whether the event is a med event or a non-med event. If both events are kept in the same array, then a test of the med bit must be done before plotting.

Dim i As Integer, dTime As Double, iDoseDay As Integer, iZoomPanelWidth As Integer
 Static bProcedureInProgress
 If bProcedureInProgress Then Exit Sub
 bProcedureInProgress = True

On Error Resume Next *'prevent error if already unloaded*
 For i = 1 To gZoomDosesCreated *'remove all previous Doses*
 Unload frmDosingCalendar.shapeZoomDose(i)
 Next i
 gZoomDosesCreated = 0

For i = 1 To 4
 Unload frmDosingCalendar.shapeZoomPrescribed(i) *'clear the text box for zoom time*
 Unload frmDosingCalendar.shapeZoomTimeRange(i) *'clear the text box for zoom time*
 Next i

On Error GoTo 0 *'resume normal error status*

iZoomPanelWidth = frmDosingCalendar.pnlZoom.Width *'speed up process by defining width from control*

For i = 1 To 4
 If PAT_DATA.dPrescribedDoseTime(i) >= 0 Then
 dTime = PAT_DATA.dPrescribedDoseTime(i) - Int(PAT_DATA.dPrescribedDoseTime(i))
 Load frmDosingCalendar.shapeZoomTimeRange(i)
 frmDosingCalendar.shapeZoomTimeRange(i).Left = (iZoomPanelWidth * dTime) - (iZoomPanelWidth * (gsngComplianceTimeRange / 24))
 frmDosingCalendar.shapeZoomTimeRange(i).Width = (iZoomPanelWidth * (gsngComplianceTimeRange / 24) * 2)
 frmDosingCalendar.shapeZoomTimeRange(i).ToolTipText = "Compliance Time Range = " & CStr(gsngComplianceTimeRange) & " Hrs."
 frmDosingCalendar.shapeZoomTimeRange(i).Visible = True
 frmDosingCalendar.shapeZoomTimeRange(i).ZOrder
 End If
 Next i

For i = 1 To 4
 If PAT_DATA.dPrescribedDoseTime(i) >= 0 Then
 dTime = PAT_DATA.dPrescribedDoseTime(i) - Int(PAT_DATA.dPrescribedDoseTime(i))
 Load frmDosingCalendar.shapeZoomPrescribed(i)
 frmDosingCalendar.shapeZoomPrescribed(i).Left = (iZoomPanelWidth * dTime) - (frmDosingCalendar.shapeZoomPrescribed(i).Width / 2) + 15
 frmDosingCalendar.shapeZoomPrescribed(i).ToolTipText = Format\$(dTime, gsTimeDisplayFormat)
 frmDosingCalendar.shapeZoomPrescribed(i).Visible = True
 frmDosingCalendar.shapeZoomPrescribed(i).ZOrder
 End If
 Next i

For i = 1 To 4
 frmDosingCalendar.txtZoomTime(i).Caption = "" *'clear the text box for zoom time*
 Next i

Dim iCalendarDate As Long
 iCalendarDate = DateValue(frmDosingCalendar.Calendar.Date)
 For i = 1 To gCompliedDosesCreated
'gh we may later want to use a global array instead of the tag property to prevent flashing and speed things up.
 iDoseDay = frmDosingCalendar.shapeDose(i).Tag *'get the day that the dose was taken on*
 If Int(PAT_DATA.dEventDate(iDoseDay)) = iCalendarDate Then
'Creates another clone of the the Dose shape located in the array Dose(0)
 gZoomDosesCreated = gZoomDosesCreated + 1 *'increment Dose counter*

Calendar.bas - UpdateZoomBox

72

```

Load frmDosingCalendar.shapeZoomDose(giZoomDosesCreated) 'create a new Dose
dTime = PAT_DATA.dEventDate(iDoseDay) - Int(PAT_DATA.dEventDate(iDoseDay))
frmDosingCalendar.shapeZoomDose(giZoomDosesCreated).Left = (frmDosingCalendar.pnZoom.Width * dTime) - (
frmDosingCalendar.shapeZoomDose(giZoomDosesCreated).Width / 2)
frmDosingCalendar.shapeZoomDose(giZoomDosesCreated).ToolTipText = Format$(dTime, gsTimeDisplayFormat)
frmDosingCalendar.shapeZoomDose(giZoomDosesCreated).Visible = True
frmDosingCalendar.shapeZoomDose(giZoomDosesCreated).ZOrder
End If
Next i

For i = 1 To giNonCompliedDosesCreated
iDoseDay = frmDosingCalendar.shapeDoseNonComply(i).Tag 'get the day that the dose was taken on
If Int(PAT_DATA.dEventDate(iDoseDay)) = iCalendarDate Then
'Create another clone of the the Dose shape located in the array Dose(0)
giZoomDosesCreated = giZoomDosesCreated + 1 'increment Dose counter
Load frmDosingCalendar.shapeZoomDose(giZoomDosesCreated) 'create a new Dose
dTime = PAT_DATA.dEventDate(iDoseDay) - Int(PAT_DATA.dEventDate(iDoseDay))
frmDosingCalendar.shapeZoomDose(giZoomDosesCreated).Left = (frmDosingCalendar.pnZoom.Width * dTime) - (
frmDosingCalendar.shapeZoomDose(giZoomDosesCreated).Width / 2)
frmDosingCalendar.shapeZoomDose(giZoomDosesCreated).ToolTipText = Format$(dTime, gsTimeDisplayFormat)
frmDosingCalendar.shapeZoomDose(giZoomDosesCreated).Visible = True
frmDosingCalendar.shapeZoomDose(giZoomDosesCreated).ZOrder
End If
Next i

frmDosingCalendar.pnZoom.Caption = Format$(frmDosingCalendar.Calendar.Date, "General Date") & " Detail View"

'update position of time scale
For i = 2 To 22 Step 2
frmDosingCalendar.lblDetailTime(i).Left = (frmDosingCalendar.pnZoom.Width * (i / 24)) - (frmDosingCalendar.lblDetailTime(i).Width / 2)
Next i
frmDosingCalendar.shapeDayLight(2).Width = frmDosingCalendar.pnZoom.Width * 0.53
frmDosingCalendar.shapeDayLight(1).Width = frmDosingCalendar.pnZoom.Width * 0.03
frmDosingCalendar.shapeDayLight(3).Width = frmDosingCalendar.shapeDayLight(1).Width
frmDosingCalendar.shapeDayLight(2).Left = (frmDosingCalendar.pnZoom.Width - frmDosingCalendar.shapeDayLight(2).Width) / 1.8
frmDosingCalendar.shapeDayLight(1).Left = 20 + frmDosingCalendar.shapeDayLight(2).Left - frmDosingCalendar.shapeDayLight(1).Width
frmDosingCalendar.shapeDayLight(3).Left = frmDosingCalendar.shapeDayLight(2).Left + frmDosingCalendar.shapeDayLight(2).Width - 15

bProcedureInProgress = False
End Sub

```

Private Sub MoveFormObjects(frm As Form, XOffset As Integer, YOffset As Integer, VisibleOnly As Integer)

```

' This routine moves all objects on a form by the specified amount
' Argument Description
' frm Form object
' XOffset offset (in twips) to move in x plane. Positive is to the right
' YOffset offset (in twips) to move in y plane. Positive is down.
' VisibleOnly If true only move visible objects
Dim i As Integer
On Error GoTo Error_MoveFormObjects

```

```

' loop through all the forms on the form
For i = 0 To frm.Controls.Count - 1
' if 1) the processing only visible controls and the controls is visible
' or 2) processing all controls
If frm.Controls(i).Tag <> "contained" Then
If (VisibleOnly And frm.Controls(i).Visible) Or Not VisibleOnly Then
' reset left and top properties
frm.Controls(i).Left = frm.Controls(i).Left + XOffset
frm.Controls(i).Top = frm.Controls(i).Top + YOffset
End If
End If
Next i

```

Calendar.bas - MoveFormObjects

73

Exit_MoveFormObjects:
Exit Sub

Error_MoveFormObjects:
Resume Exit_MoveFormObjects

End Sub

Private Sub DrawSingleCompliedDoseTaken(iDay As Integer, dTime As Double, iEventNumber As Integer, iPlot

Draw dosing Doses for the day of the month and time of day passed in here.

Time is expressed in decimal places as a portion of a day (VB time format)

NOTE: No checks are currently made to determine whether the event is a med event or a non-med event. If both events are kept in the same array, then a test of the med bit must be done before plotting.

A Dose is not visible when first created. The caller should display the Doses once they are all created, so as to speed the redraw of the screen.

Create another clone of the the Dose shape located in the array Dose(0)

On Error Resume Next

Dim i As Integer, iWeekDay As Integer

Dim iDoseLeft As Single, iDoseTop As Single

Dim iTemp As Long, iDayWidth As Long, iDayHeight As Long

giCompliedDosesCreated = giCompliedDosesCreated + 1

Increment Dose counter

Load frmDosingCalendar.shapeDose(giCompliedDosesCreated)

create a new Dose

These lines are a work-around for a bug in the control that causes it to

return the wrong values for day.left, day.top, etc. When fixed, we can simply use those properties

and remove these calculations.

iDayWidth = (frmDosingCalendar.Calendar.Width - 50) / 7

actual scalewidth of a single day

iWeekDay = (frmDosingCalendar.Calendar.DayLeft(iDay) * 25) / iDayWidth

approximate location of the day

iDoseLeft = (iWeekDay * iDayWidth) *'get left edge of day to plot*

iDoseLeft = iDoseLeft + ((iDayWidth / 5) * iPlotPosition - 1) - (iDayWidth / 10)

frmDosingCalendar.shapeDose(giCompliedDosesCreated).Left = iDoseLeft

These lines are a work-around for a bug in the control that causes it to

return the wrong values for day.left, day.top, etc. When fixed, we can simply use those properties

and remove these calculations.

The control is even more stupid than I first suspected. It can not always return the proper vertical location of a day, thus, we have to jump through more hoops to figure out what week that a particular day is in.

iDayHeight = (frmDosingCalendar.Calendar.Height - 50) - 625

'an offset is used to compensate for height of title

iDayHeight = iDayHeight / 6

iWeekDay = (frmDosingCalendar.Calendar.DayLeft(1) * 25) / iDayWidth

'the approximate location of the day

iTemp = Int((iDay + iWeekDay - 1) / 7)

iDoseTop = (iTemp * iDayHeight) + 625 *'this number factors in the title bar*

iDoseTop = iDoseTop + iDayHeight - 25 - frmDosingCalendar.shapeDose(0).Height

'draw in bottom of day

frmDosingCalendar.shapeDose(giCompliedDosesCreated).Top = iDoseTop

frmDosingCalendar.shapeDose(giCompliedDosesCreated).Tag = iEventNumber

'keep event number for updating the zoom box

On Error GoTo 0

End Sub

Calendar.bas - DrawSingleDoseMissed

74

Private Sub DrawSingleDoseMissed(iDay As Integer, iPlotPosition As Integer)

'Draw doses for the day of the month passed in here
'NOTE: No checks are currently made to determine whether the event is a med event or
'a non-med event. A test of the med bit must be done before calling this procedure.
'A Dose is not visible when first created. The caller should display the Doses once they
'are all created, so as to speed the redraw of the screen.

On Error Resume Next
 Dim i As Integer, iWeekDay As Integer
 Dim iDoseLeft As Single, iDoseTop As Single
 Dim iTemp As Long, iDayWidth As Long, iDayHeight As Long

giDosesMissedCreated = giDosesMissedCreated + 1 *'increment Dose counter*
'Create another clone of the the Dose shape located in the array Dose(0)
 Load frmDosingCalendar.shapeDoseMissed(giDosesMissedCreated) *'create a new Dose*

'These lines are a work-around for a bug in the control that causes it to
'return the wrong values for day.left, day.top, etc. When fixed, we can simply use those properties
'and remove these calculations.
 iDayWidth = (frmDosingCalendar.Calendar.Width - 50) / 7 *'actual scalewidth of a single day*
 iWeekDay = (frmDosingCalendar.Calendar.DayLeft(iDay) * 26) / iDayWidth *'approximate location of the day*

iDoseLeft = (iWeekDay * iDayWidth) *'get left edge of day to plot*
 iDoseLeft = iDoseLeft + ((iDayWidth / 5) * iPlotPosition - 1) - (iDayWidth / 10)
 frmDosingCalendar.shapeDoseMissed(giDosesMissedCreated).Left = iDoseLeft

'These lines are a work-around for a bug in the control that causes it to
'return the wrong values for day.left, day.top, etc. When fixed, we can simply use those properties
'and remove these calculations.
'The control is even more stupid than I first suspected. It can not always return the proper vertical
'location of a day, thus, we have to jump through more hoops to figure out what week that a particular
'day is in.
 iDayHeight = (frmDosingCalendar.Calendar.Height - 50) - 625 *'an offset is used to compensate for height of title*
 iDayHeight = iDayHeight / 6
 iWeekDay = (frmDosingCalendar.Calendar.DayLeft(1) * 26) / iDayWidth *'the approximate location of the day*
 iTemp = Int((iDay + iWeekDay - 1) / 7)

iDoseTop = (iTemp * iDayHeight) + 625 *'this number factors in the title bar*
 iDoseTop = iDoseTop + iDayHeight - 75 - frmDosingCalendar.shapeDose(0).Height -
 frmDosingCalendar.shapeDose(0).Height *'draw in bottom of day*
 frmDosingCalendar.shapeDoseMissed(giDosesMissedCreated).Top = iDoseTop

On Error GoTo 0
 End Sub

Public Sub UpdateCalendar()

'The month or year to the calendar has changed, so we need to plot the events for
'the current month and year being shown.

Static bProcedureInProgress As Boolean
 If bProcedureInProgress Then Exit Sub
 bProcedureInProgress = True

Dim iObjectDiameter As Integer

'Show custom labels from config file if there were any
 If Len(gsCustomLbiPatientLastName) > 0 Then frmDosingCalendar.Label1 = gsCustomLbiPatientLastName
 frmDosingCalendar.lbiPatientName = "" + PAT_DATA.sPatientLastName + ", " + PAT_DATA.sPatientFirstName

iObjectDiameter = frmDosingCalendar.Calendar.Width / 45 *'resize the objects drawn on the calendar*
 If iObjectDiameter > frmDosingCalendar.Calendar.Height / 50 Then iObjectDiameter = frmDosingCalendar.Calendar.Height / 50

Calendar.bas - UpdateCalendar

75

```
frmDosingCalendar.shapeDose(0).Width = iObjectDiameter
frmDosingCalendar.shapeDose(0).Height = iObjectDiameter
frmDosingCalendar.shapeDoseNonComply(0).Width = iObjectDiameter
frmDosingCalendar.shapeDoseNonComply(0).Height = iObjectDiameter
frmDosingCalendar.shapeDoseMissed(0).Width = iObjectDiameter
frmDosingCalendar.shapeDoseMissed(0).Height = iObjectDiameter
frmDosingCalendar.shapeDoseSizeChange(0).Width = iObjectDiameter
frmDosingCalendar.shapeDoseSizeChange(0).Height = iObjectDiameter
```

```
DrawAllCompliedDosesTaken
DrawAllNonCompliedDosesTaken
DrawAllDosesMissed
DrawAllDoseSizeChanges
```

```
UpdateZoomBox
```

```
bProcedureInProgress = False
End Sub
```

Public Sub RemoveAllObjects()

```
'Remove all objects from calendar
RemoveCompliedDosesTaken
RemoveDosesMissed
DoEvents
End Sub
```


frmMain.frm - File Declarations

76

```

Attribute VB_Name = "frmMain"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit

```

```

Private Sub CommTimer_Timer()
    gbCommTimerExpired = True
    CommTimer.Enabled = False
End Sub

```

```

Private Sub FaxMan1_ConfigurationDone()
    Dim i As Integer
    frmOptions.MousePointer = vbDefault
    frmOptions.btnConfigureFax.Enabled = True
    SetFaxDeviceLabel
End Sub

```

```

Private Sub FaxMan1_FaxStatus(Device As Integer, Status As Integer)
    Beep
    If gcFax.Status(Device) = "Initializing Modem" Or gcFax.Status(Device) = "Answering" Then
        frmFaxStatus.Show
    ElseIf gcFax.Status(Device) = "Port Closed" Then
        Unload frmFaxStatus
    End If

    frmFaxStatus.lblRemoteID = gcFax.StatusRemoteID(Device)

    If gcFax.StatusConnectSpeed(Device) > 0 Then
        frmFaxStatus.lblSpeed = gcFax.StatusConnectSpeed(Device)
    Else
        frmFaxStatus.lblSpeed = ""
    End If

    If gcFax.StatusPages(Device) Then
        frmFaxStatus.lblPage = CStr(gcFax.StatusPagesSent(Device)) + " of " + Str$(gcFax.StatusPages(Device))
    Else
        frmFaxStatus.lblPage = CStr(gcFax.StatusPagesSent(Device))
    End If

    If gcFax.StatusPercentage(Device) > 0 Then
        frmFaxStatus.lblPercent = CStr(gcFax.StatusPercentage(Device)) + "% Complete"
    Else
        frmFaxStatus.lblPercent = ""
    End If

    frmFaxStatus.lblStatus = gcFax.Status(Device)
    frmFaxStatus.lblDestination = gcFax.StatusDestination(Device)
    frmFaxStatus.lblFaxNumber = gcFax.StatusNumber(Device)
End Sub

```

frmMain.frm - MDIForm_Load

77

Private Sub MDIForm_Load()

```

On Error Resume Next
Me.Left = CLng(GetINISetting(gsAppIniFileSpec, "Windows", "Main Left", "1000"))
Me.Top = CLng(GetINISetting(gsAppIniFileSpec, "Windows", "Main Top", "1000"))
Me.Width = CLng(GetINISetting(gsAppIniFileSpec, "Windows", "Main Width", "6500"))
Me.Height = CLng(GetINISetting(gsAppIniFileSpec, "Windows", "Main Height", "6500"))
Me.WindowState = CLng(GetINISetting(gsAppIniFileSpec, "Windows", "Main WindowState", "0"))
On Error GoTo 0
End Sub

```

Private Sub MDIForm_Unload(Cancel As Integer)

```

Dim r As Integer

r = ValidatePatientDataSaved 'make sure any device data has first been saved
'Save Window positions
If Me.WindowState <> vbMinimized Then
    SaveINISetting gsAppIniFileSpec, "Windows", "Main Left", CStr(Me.Left)
    SaveINISetting gsAppIniFileSpec, "Windows", "Main Top", CStr(Me.Top)
    SaveINISetting gsAppIniFileSpec, "Windows", "Main Width", CStr(Me.Width)
    SaveINISetting gsAppIniFileSpec, "Windows", "Main Height", CStr(Me.Height)
    SaveINISetting gsAppIniFileSpec, "Windows", "Main WindowState", CStr(Me.WindowState)
End If

SaveProgramPreferences
End Sub

```

Private Sub mnuAccessWebSite_Click()

```

If the form is minimized then set it back to normal
Call LogonToWebSite
If frmBrowser.WindowState = vbMinimized Then
    frmBrowser.WindowState = vbNormal
End If
frmBrowser.ZOrder
End Sub

```

Private Sub mnuFaxConfigure_Click()

```

giLatestOptionsTabSelected = 2 'display the fax tab once the dialog is opened
frmOptions.Show vbModal
End Sub

```

Private Sub mnuFaxSend_Click()

```

frmFaxSend.Show
End Sub

```

Private Sub mnuFaxViewLogs_Click()

```

frmFaxLog.Show
End Sub

```

frmMain.frm - mnuFileProperties_Click

78

```
Private Sub mnuFileProperties_Click()
    frmOptions.Show vbModal
End Sub
```

```
Private Sub mnuFileSave_Click()
    Dim r As Integer
```

```
    If PAT_DATA.sPatientDataFileName = "" Then
        r = SaveDataToNewFile
    Else
        r = SavePatientData(PAT_DATA.sPatientDataFileName)
    End If
```

```
    If r = False Then
        Beep
        MsgBox "An error occurred while attempting to save the data file. It was not saved.", vbCritical, "File Not Saved"
    End If
End Sub
```

```
Private Sub mnuGenError_Click()
```

```
    MsgBox "This is a temporary test error handler. When you click OK, a synthetic error (Divide by 0) will be generated. The same dialog will be shown when any error is generated. It generates a log file that provides valuable information for the developer. This will be removed from the next build.", vbInformation, "Test Error Handler"
    Error 11
End Sub
```

```
Private Sub mnuHelpDeviceDiag_Click()
```

```
    Dim sMSG, sReply As String
    sMSG = "Performing a device diagnostics test could cause loss of vital device information and should be done only with the assistance of technical support."
    sMSG = sMSG + vbCrLf + vbCrLf + "Please contact our technical support department at 1-800-777-7777 for a password and assistance."
    ' Display message, title, and default value.
    sReply = InputBox(sMSG, "Password Required")
    If LCase$(sReply) = "h2o" Then frmDeviceDiagnostics.Show
End Sub
```

```
Private Sub mnuHelpTips_Click()
```

```
    frmTip.Show
    ' If the form is minimized then set it back to normal
    If frmTip.WindowState = vbMinimized Then frmTip.WindowState = vbNormal
    frmTip.ZOrder
End Sub
```

frmMain.frm - mnuReadDeviceData_Click

79

```
Private Sub mnuReadDeviceData_Click()  
    frmReadDeviceData.Show  
    'If the form is minimized then set it back to normal  
    If frmReadDeviceData.WindowState = vbMinimized Then frmReadDeviceData.WindowState = vbNormal  
    frmReadDeviceData.ZOrder  
End Sub
```

```
Private Sub mnuSendDeviceData_Click()  
    frmDeviceInitialize.Show  
    'If the form is minimized then set it back to normal  
    If frmDeviceInitialize.WindowState = vbMinimized Then frmDeviceInitialize.WindowState = vbNormal  
    frmDeviceInitialize.ZOrder  
End Sub
```

```
Private Sub mnuViewAllPatients_Click()  
    frmAllPatients.Show  
    'If the form is minimized then set it back to normal  
    If frmAllPatients.WindowState = vbMinimized Then frmAllPatients.WindowState = vbNormal  
    frmAllPatients.ZOrder  
End Sub
```

```
Private Sub mnuHelpAbout_Click()  
    frmAbout.Show vbModal, Me  
End Sub
```

```
Private Sub mnuViewCalendar_Click()  
    frmDosingCalendar.Show  
    'If the form is minimized then set it back to normal  
    If frmDosingCalendar.WindowState = vbMinimized Then frmDosingCalendar.WindowState = vbNormal  
    frmDosingCalendar.ZOrder  
End Sub
```

```
Private Sub mnuViewExplorer_Click()  
    mnuViewExplorer.Checked = Not mnuViewExplorer.Checked 'toggle the state of the check box  
    SSListBar1.Visible = mnuViewExplorer.Checked  
End Sub
```

```
Private Sub mnuViewOptions_Click()  
    frmOptions.Show vbModal, Me  
End Sub
```

frmMain.frm - mnuViewPatientDosingReport

80

Private Sub mnuViewPatientDosingReport_Click()

```

frmPatientDosingReport.Show
'if the form is minimized then set it back to normal
If frmPatientDosingReport.WindowState = vbMinimized Then frmPatientDosingReport.WindowState = vbNormal
frmPatientDosingReport.ZOrder
End Sub

```

Private Sub mnuViewPatientSummary_Click()

```

frmPatientSummary.Show
'if the form is minimized then set it back to normal
If frmPatientSummary.WindowState = vbMinimized Then frmPatientSummary.WindowState = vbNormal
frmPatientSummary.ZOrder
End Sub

```

Private Sub mnuViewStatusBar_Click()

```

If mnuViewStatusBar.Checked Then
    sbStatusBar.Visible = False
    mnuViewStatusBar.Checked = False
Else
    sbStatusBar.Visible = True
    mnuViewStatusBar.Checked = True
End If
End Sub

```

Private Sub mnuViewToolBar_Click()

```

If mnuViewToolBar.Checked Then
    tbToolBar.Visible = False
    mnuViewToolBar.Checked = False
Else
    tbToolBar.Visible = True
    mnuViewToolBar.Checked = True
End If
End Sub

```

Private Sub SSListBar1_ListItemClick(ByVal ItemClicked As Listbar.SSListItem)

```

Select Case SSListBar1.CurrentGroupKey
Case "Patient Data" 'patient data
    Select Case ItemClicked.Key
    Case "Event Calendar" 'calendar
        mnuViewCalendar_Click
    Case "Summary" 'summary
        mnuViewPatientSummary_Click
    Case "Dosing Information" 'grid
        mnuViewPatientDosingReport_Click
    Case "All Patients" 'all patients
        mnuViewAllPatients_Click
    End Select

```

```

Case "Device" 'device data
    Select Case ItemClicked.Key

```

frmMain.frm - SSListBar1_ListItemClick

81

```

Case "Retrieve Data" 'read device data
    mnuReadDeviceData_Click
Case "Program Device" 'send to DosPro Device
    mnuSendDeviceData_Click
End Select
End Sub

```

Private Sub tbToolBar_ButtonClick(ByVal Button As ComctlLib.Button)

Select Case Button.Key

Case "Open"
mnuFileOpen_Click

Case "Save"
mnuFileSave_Click

Case "Print"
mnuFilePrint_Click

Case "Cut"
mnuEditCut_Click

Case "Copy"
mnuEditCopy_Click

```

Clipboard.Clear
If TypeOf ActiveForm.ActiveControl Is TextBox Then
    Select Case Index
        Case 0 ' Cut
            Copy selected text to Clipboard.

```

```

Clipboard.SetText ActiveForm.ActiveControl.Text
Delete selected text.
ActiveForm.ActiveControl.Text = ""
Case 1 ' Copy
    Copy selected text to Clipboard.

Clipboard.SetText ActiveForm.ActiveControl.Text
Case 2 ' Paste
    Put Clipboard text in text box.
    ActiveForm.ActiveControl.Text = Clipboard.GetText()
Case 3 ' Delete
    Delete selected text.
    ActiveForm.ActiveControl.Text = ""

```

frmMain.frm - tbToolBar_ButtonClick

82

```

        End Select
    End If

    Case "Paste"
        mnuEditPaste_Click

    Case "Bold"

    Case "Italic"

    Case "Underline"

    Case "Left"

    Case "Center"

    Case "Right"

    End Select
End Sub

```

Private Sub mnuHelpContents_Click()

Dim nRet As Integer

*If there is no helpfile for this project display a message to the user
you can set the HelpFile for your application in the
Project Properties dialog*

If Len(App.HelpFile) = 0 Then

MsgBox "Unable to display Help Contents. There is no Help associated with this project.", vbInformation, Me.Caption

Else

On Error Resume Next

nRet = OSWinHelp(Me.hWnd, App.HelpFile, 3, 0)

If Err Then MsgBox Err.Description

End If

End Sub

Private Sub mnuHelpSearch_Click()

Dim nRet As Integer

*If there is no helpfile for this project display a message to the user
you can set the HelpFile for your application in the
Project Properties dialog*

If Len(App.HelpFile) = 0 Then

MsgBox "Unable to display Help Contents. There is no Help associated with this project.", vbInformation, Me.Caption

Else

On Error Resume Next

nRet = OSWinHelp(Me.hWnd, App.HelpFile, 261, 0)

If Err Then MsgBox Err.Description

End If

End Sub

frmMain.frm - mnuWindowArrangeIcons_

83

```
Private Sub mnuWindowArrangeIcons_Click()
    Me.Arrange vbArrangeIcons
End Sub
```

```
Private Sub mnuWindowCascade_Click()
    Me.Arrange vbCascade
End Sub
```

```
Private Sub mnuWindowTileHorizontal_Click()
    Me.Arrange vbTileHorizontal
End Sub
```

```
Private Sub mnuWindowTileVertical_Click()
    Me.Arrange vbTileVertical
End Sub
```

```
Private Sub mnuFileOpen_Click()
    Dim r As Integer
    r = OpenPatientData()
```

*If any of these forms are open at the time a new file is loaded,
then refresh them.*

```
For r = 0 To Forms.Count - 1
```

```
    Select Case Forms(r).Name
```

```
        Case "frmPatientDosingReport"
```

```
            frmPatientDosingReport.UpdatePatientGridDisplay
```

```
        Case "frmDosingCalendar"
```

```
            UpdateCalendar
```

```
        Case "frmPrint"
```

```
            RefreshPreview
```

```
    End Select
```

```
Next r
```

```
End Sub
```

```
Private Sub mnuFileSaveAs_Click()
    SaveDataToNewFile
End Sub
```


frmMain.frm - mnuFilePageSetup_67

84

```
Private Sub mnuFilePageSetup_Click()  
On Error GoTo mnuFilePageSetup_Click_Error
```

```
    dlgCommonDialog.ShowPrinter
```

```
mnuFilePageSetup_Click_Exit:  
On Error GoTo 0  
Exit Sub
```

```
mnuFilePageSetup_Click_Error:  
Resume mnuFilePageSetup_Click_Exit 'any error message would have already been sent by the common dialog  
End Sub
```

```
Private Sub mnuFilePrint_Click()
```

```
    frmPrint.Show  
    'if the form is minimized then set it back to normal  
    If frmPrint.WindowState = vbMinimized Then frmPrint.WindowState = vbNormal  
    frmPrint.ZOrder  
End Sub
```

```
Private Sub mnuFileSend_Click()
```

```
    'To Do  
    MsgBox "Ability to send a file will be active in a future release"  
End Sub
```

```
Private Sub mnuFileMRU_Click(Index As Integer)
```

```
    Dim r As Integer  
    r = OpenPatientData(mnuFileMRU(Index).Caption)  
End Sub
```

```
Private Sub mnuFileExit_Click()
```

```
    'unload the form  
    Unload Me  
End Sub
```

frmSplash.frm - File Declarations

```
Attribute VB_Name = "frmSplash"  
Attribute VB_GlobalNameSpace = False  
Attribute VB_Creatable = False  
Attribute VB_PredeclaredId = True  
Attribute VB_Exposed = False  
Option Explicit
```

85

```
Private Sub Form_Load()  
    lblVersion.Caption = "Version " & App.Major & "." & App.Minor & " " & App.Revision  
    lblProductName.Caption = App.Title  
End Sub
```

frmLogin.frm - File Declaration

```
Attribute VB_Name = "frmLogin"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit
```

```
Private Declare Function GetUserName Lib "advapi32.dll" Alias "GetUserNameA" (ByVal lpBuffer As String, nSize As Long) As Long
Public OK As Boolean
```

Private Sub Form_Load()

```
Dim sBuffer As String
Dim lSize As Long

Me.Move frmSplash.Left + 4000, frmSplash.Top + 3500

sBuffer = Space$(255)
lSize = Len(sBuffer)
Call GetUserName(sBuffer, lSize)
If lSize > 0 Then
    txtUserName = Left$(sBuffer, lSize)
Else
    txtUserName = vbNullString
End If
End Sub
```

Private Sub cmdCancel_Click()

```
OK = False
Me.Hide
End Sub
```

Private Sub cmdOK_Click()

```
'To Do - create test for correct password
'check for correct password
Me.MousePointer = vbHourglass
If txtPassword = "" Then
    OK = True
    imgLocked.Visible = False
    imgUnlocked.Visible = True
    Wait 1.5
    Me.MousePointer = vbDefault
    Me.Hide
Else
    imgLocked.Visible = False
    Wait 0.05
    imgLocked.Visible = True
    Wait 0.05
    imgLocked.Visible = False
    Wait 0.05
    imgLocked.Visible = True
    Wait 0.05
    imgLocked.Visible = False
    Wait 0.05
    imgLocked.Visible = True
    Me.MousePointer = vbDefault
    Beep
    MsgBox "Invalid Password, try again.", "Login"
    txtPassword.SetFocus
    txtPassword.SelStart = 0
```

frmLogin.frm - cmdOK_Click

```
txtPassword.SetLength = Len(txtPassword)  
End If  
End Sub
```

87

frmOptions.frm - File Declarations

88

```

Attribute VB_Name = "frmOptions"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit

```

```

Private Sub btnConfigureFax_Click()
    btnConfigureFax.Enabled = False
    Me.MousePointer = vbHourglass
    gcFax.AutoDetect
    With lblFaxDevice
        .Caption = "Searching for a Fax Device. Please wait a few seconds."
        .BackColor = &HCOFFFF highlight background
        .ForeColor = &H0&
        .Refresh
    End With
End Sub

```

```

Private Sub cmdApply_Click()
    Dim sSection As String

    glLatestOptionsTabSelected = sstab1.Tab

    'set the global value to the user's selection
    gsDateDisplayFormat = Choose(cmbDates.ListIndex + 1, "Short Date", "Medium Date", "Long Date")
    gsTimeDisplayFormat = Choose(cmbTimes.ListIndex + 1, "Short Time", "Medium Time", "Long Time")

```

```

Select Case cmbComplianceTimeRange.ListIndex
Case 0
    gsngComplianceTimeRange = 0.5
Case 1
    gsngComplianceTimeRange = 1
Case 2
    gsngComplianceTimeRange = 1.5
Case 3
    gsngComplianceTimeRange = 2
Case 4
    gsngComplianceTimeRange = 2.5
Case 5
    gsngComplianceTimeRange = 3
Case 6
    gsngComplianceTimeRange = 3.5
Case 7
    gsngComplianceTimeRange = 4
Case 8
    gsngComplianceTimeRange = 4.5
Case 9
    gsngComplianceTimeRange = 5
Case 10
    gsngComplianceTimeRange = 5.5
Case 11
    gsngComplianceTimeRange = 6

```

```
End Select
```

```
RefreshAllOpenForms
```

```

'Save the Fax Information
With FAX_DATA
    .sSenderName = txtName
    .sSenderCompany = txtCompany
    .sSenderVoiceNumber = txtVoiceNumber
    .sSenderFaxNumber = txtFaxNumber

```

frmOptions.frm - cmdApply_Click

89

```

.sFaxID = txtFaxID
.sDialPrefix = txtDialPrefix
.iRetries = Val(txtRetries)
.iRetryInterval = Val(txtRetryInterval)
.bFaxResolution = chkResolution.Value '0=low, 1 = high
End With

```

'Save values for the Fax control that was last set by user
sSection = "User Selections"

With FAX_DATA

```

SaveINISetting gsFaxFileSpec, sSection, "Sender Name", .sSenderName
SaveINISetting gsFaxFileSpec, sSection, "Sender Company", .sSenderCompany
SaveINISetting gsFaxFileSpec, sSection, "Sender Voice Number", .sSenderVoiceNumber
SaveINISetting gsFaxFileSpec, sSection, "Sender Fax Number", .sSenderFaxNumber
SaveINISetting gsFaxFileSpec, sSection, "Fax ID", .sFaxID
SaveINISetting gsFaxFileSpec, sSection, "Dial Prefix", .sDialPrefix
SaveINISetting gsFaxFileSpec, sSection, "Retries", CStr(.iRetries)
SaveINISetting gsFaxFileSpec, sSection, "Retry Interval", CStr(.iRetryInterval)
SaveINISetting gsFaxFileSpec, sSection, "Resolution", CStr(.bFaxResolution)
End With
End Sub

```

```

Private Sub cmdCancel_Click()
Unload Me
End Sub

```

```

Private Sub cmdOK_Click()
'Code goes here to set options and close dialog.
cmdApply_Click
Unload Me
End Sub

```

```

Private Sub Form_Activate()
SetPrinterIcon False, ""
End Sub

```

```

Private Sub Form_Load()
'Define the mask for the telephone and fax numbers text box
Load the available choices into the list boxes
cmboDates.AddItem Format$(Now, "Short Date")
cmboDates.AddItem Format$(Now, "Medium Date")
cmboDates.AddItem Format$(Now, "Long Date")

cmboTimes.AddItem Format$(Now, "Short Time") + " (24 hour)"
cmboTimes.AddItem Format$(Now, "Medium Time") + " (12 hour)"
cmboTimes.AddItem Format$(Now, "Long Time")

```

```

Select Case gsngComplianceTimeRange
Case 0.5
cmboComplianceTimeRange.ListIndex = 0
Case 1
cmboComplianceTimeRange.ListIndex = 1
Case 1.5
cmboComplianceTimeRange.ListIndex = 2
Case 2
cmboComplianceTimeRange.ListIndex = 3
Case 2.5
cmboComplianceTimeRange.ListIndex = 4

```

frmOptions.frm - Form_Load

90

```

Case 3
  cmboComplianceTimeRange.ListIndex = 5
Case 3.5
  cmboComplianceTimeRange.ListIndex = 6
Case 4
  cmboComplianceTimeRange.ListIndex = 7
Case 4.5
  cmboComplianceTimeRange.ListIndex = 8
Case 5
  cmboComplianceTimeRange.ListIndex = 9
Case 5.5
  cmboComplianceTimeRange.ListIndex = 10
Case 6
  cmboComplianceTimeRange.ListIndex = 11
End Select

```

```

' set the list box to the last selected user state/
Select Case gsDateDisplayFormat
Case ""
  cmboDates.ListIndex = 0
Case "Short Date"
  cmboDates.ListIndex = 0
Case "Medium Date"
  cmboDates.ListIndex = 1
Case "Long Date"
  cmboDates.ListIndex = 2
End Select

```

```

Select Case gsTimeDisplayFormat
Case ""
  cmboTimes.ListIndex = 0
Case "Short Time"
  cmboTimes.ListIndex = 0
Case "Medium Time"
  cmboTimes.ListIndex = 1
Case "Long Time"
  cmboTimes.ListIndex = 2
End Select

```

```

' Get Fax Info from global settings to the Fax Tab
With FAX_DATA
  txtName = .sSenderName
  txtCompany = .sSenderCompany
  txtVoiceNumber = .sSenderVoiceNumber
  txtFaxNumber = .sSenderFaxNumber
  txtFaxID = .sFaxID
  txtDialPrefix = .sDialPrefix
  txtRetries = CStr(.iRetries)
  txtRetryInterval = CStr(.iRetryInterval)
  chkResolution.Value = .bFaxResolution '0=low, 1=high
End With

```

```

ssTab1.Tab = giLatestOptionsTabSelected
SetFaxDeviceLabel 'update the devices label
End Sub

```

frmOptions.frm - Form_Load

91

```
Private Sub sstab1_Click(PreviousTab As Integer)
    Tax tab
    If sstab1.Tab = 2 Then SetFaxDeviceLabel    update the devices label
End Sub
```

```
Private Sub txtFaxNumber_GotFocus()
    txtFaxNumber.SetStart = 1
End Sub
```

```
Private Sub txtRetries_GotFocus()
    txtRetries.SetStart = 1
End Sub
```

```
Private Sub txtRetryInterval_GotFocus()
    txtRetryInterval.SetStart = 1
End Sub
```

```
Private Sub txtVoiceNumber_GotFocus()
    txtVoiceNumber.SetStart = 1
End Sub
```

frmAbout.frm - File Declarations

92

```
Attribute VB_Name = "frmAbout"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit
```

```
* Reg Key Security Options...
Const KEY_ALL_ACCESS = &H2003F
```

```
* Reg Key ROOT Types...
Const HKEY_LOCAL_MACHINE = &H80000002
Const ERROR_SUCCESS = 0
Const REG_SZ = 1           'Unicode nul terminated string
Const REG_DWORD = 4        '32-bit number
```

```
Const gREGKEYSYSINFOLOC = "SOFTWARE\Microsoft\Shared Tools Location"
Const gREGVALSYSINFOLOC = "MSINFO"
Const gREGKEYSYSINFO = "SOFTWARE\Microsoft\Shared Tools\MSINFO"
Const gREGVALSYSINFO = "PATH"
```

```
Private Declare Function RegOpenKeyEx Lib "advapi32" Alias "RegOpenKeyExA" (ByVal hKey As Long, ByVal lpSubKey As String, ByVal  
  ulOptions As Long, ByVal samDesired As Long, ByRef phkResult As Long) As Long
Private Declare Function RegQueryValueEx Lib "advapi32" Alias "RegQueryValueExA" (ByVal hKey As Long, ByVal lpValueName As String,  
  ByVal lpReserved As Long, ByRef lpType As Long, ByVal lpData As String, ByRef lpcbData As Long) As Long
Private Declare Function RegCloseKey Lib "advapi32" (ByVal hKey As Long) As Long
```

```
Private Sub Form_Load()
  lblVersion.Caption = "Version " & App.Major & "." & App.Minor & "." & App.Revision
  lblTitle.Caption = App.Title
End Sub
```

```
Private Sub cmdSysInfo_Click()
  Call StartSysInfo
End Sub
```

```
Private Sub cmdOK_Click()
  Unload Me
End Sub
```

```
Public Sub StartSysInfo()
  On Error GoTo SysInfoErr
```

```
Dim rc As Long
Dim SysInfoPath As String
```

```
* Try To Get System Info Program PathName From Registry...
If GetKeyValue(HKEY_LOCAL_MACHINE, gREGKEYSYSINFO, gREGVALSYSINFO, SysInfoPath) Then
  * Try To Get System Info Program Path Only From Registry...
ElseIf GetKeyValue(HKEY_LOCAL_MACHINE, gREGKEYSYSINFOLOC, gREGVALSYSINFOLOC, SysInfoPath) Then
  * Validate Existence Of Known 32 Bit File Version
```

frmAbout.frm - StartSysInfo

93

```

If (Dir(SysInfoPath + "\MSINFO32.EXE") <> "") Then
    SysInfoPath = SysInfoPath + "\MSINFO32.EXE"
    ' Error - File Can Not Be Found...
Else
    GoTo SysInfoErr
End If
' Error - Registry Entry Can Not Be Found...
Else
    GoTo SysInfoErr
End If

```

```

Call Shell(SysInfoPath, vbNormalFocus)
Exit Sub

```

```

SysInfoErr:
    MsgBox "System Information Is Unavailable At This Time", vbOKOnly
End Sub

```

Public Function GetKeyValue(KeyRoot As Long, KeyName As String, SubKeyRef As String, ByRef KeyVal As S

```

Dim i As Long          ' Loop Counter
Dim rc As Long          ' Return Code
Dim hKey As Long        ' Handle To An Open Registry Key
Dim hDepth As Long      '
Dim KeyValType As Long  ' Data Type Of A Registry Key
Dim tmpVal As String    ' Temporary Storage For A Registry Key Value
Dim KeyValSize As Long  ' Size Of Registry Key Variable

' Open RegKey Under KeyRoot (HKEY_LOCAL_MACHINE...)
rc = RegOpenKeyEx(KeyRoot, KeyName, 0, KEY_ALL_ACCESS, hKey) ' Open Registry Key
If (rc <> ERROR_SUCCESS) Then GoTo GetKeyError ' Handle Error...
tmpVal = String$(1024, 0) ' Allocate Variable Space
KeyValSize = 1024 ' Mark Variable Size

' Retrieve Registry Key Value...
rc = RegQueryValueEx(hKey, SubKeyRef, 0, KeyValType, tmpVal, KeyValSize) ' Get/Create Key Value
If (rc <> ERROR_SUCCESS) Then GoTo GetKeyError ' Handle Errors

If (Asc(Mid(tmpVal, KeyValSize, 1)) = 0) Then ' Win95 Adds Null Terminated String...
    tmpVal = Left(tmpVal, KeyValSize - 1) ' Null Found. Extract From String
Else ' WinNT Does NOT Null Terminate String...
    tmpVal = Left(tmpVal, KeyValSize) ' Null Not Found. Extract String Only
End If

' Determine Key Value Type For Conversion...
Select Case KeyValType
    Case REG_SZ ' Search Data Types...
        KeyVal = tmpVal ' String Registry Key Data Type
        ' Copy String Value
    Case REG_DWORD ' Double Word Registry Key Data Type
        For i = Len(tmpVal) To 1 Step -1 ' Convert Each Bit
            KeyVal = KeyVal + Hex(Asc(Mid(tmpVal, i, 1))) ' Build Value Char. By Char.
        Next
        KeyVal = Format$("&h" + KeyVal) ' Convert Double Word To String
End Select

GetKeyValue = True ' Return Success
rc = RegCloseKey(hKey) ' Close Registry Key
Exit Function ' Exit

```

```

GetKeyError:
    KeyVal = "" ' Cleanup After An Error Has Occured...
                ' Set Return Val To Empty String

```

frmAbout.frm - GetKeyValue

GetKeyValue = False
rc = RegCloseKey(hKey)
End Function

' Return Failure
' Close Registry key

94

///

frmBrowser.frm - File Declarations

95

```
Attribute VB_Name = "frmBrowser"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit
```

```
Public StartingAddress As String
Dim mbDontNavigateNow As Boolean
```

```
Private Sub Form_Load()
    Dim i As Integer
```

```
    'On Error Resume Next
    Me.Show
    lbToolBar.Refresh
    Form_Resize
```

```
    If Len(StartingAddress) > 0 Then
        cboAddress = StartingAddress
        cboAddress.AddItem cboAddress
        'try to navigate to the starting address
        tmTimer.Enabled = True
        brwWebBrowser.Navigate StartingAddress
        Me.MousePointer = vbHourglass
    End If
End Sub
```

```
Private Sub brwWebBrowser_DownloadComplete()
```

```
    On Error Resume Next
    Me.Caption = brwWebBrowser.LocationName
    Me.MousePointer = vbDefault
End Sub
```

```
Private Sub brwWebBrowser_NavigateComplete(ByVal URL As String)
```

```
    Dim i As Integer, r As Integer
```

```
    Dim bFound As Boolean
    On Error Resume Next
    Me.Caption = brwWebBrowser.LocationName
    For i = 0 To cboAddress.ListCount - 1
        If cboAddress.List(i) = brwWebBrowser.LocationURL Then
            bFound = True
            Exit For
        End If
    Next i
```

```
    mbDontNavigateNow = True
    If bFound Then cboAddress.RemoveItem i
    cboAddress.AddItem brwWebBrowser.LocationURL, 0
    cboAddress.ListIndex = 0
    mbDontNavigateNow = False
    On Error GoTo 0
    Me.MousePointer = vbDefault
```

```
    'Last time to visit the Internet
    'Save new date in INI file that an attempt (or success) was made to visit
    'the Internet web site on this date
    r = GetINISetting(gsAppIniFileSpec, "Web Data", "Connection Reminder Days", 100)
    SaveINISetting gsAppIniFileSpec, "Web Data", "Next Web Visit Reminder Date", Format$(Now + r, "Medium Date")
```

frmBrowser.frm - brwWebBrowser_NavigateCol :10

96

```
SaveINISetting gsAppIniFileSpec, "Web Data", "Last Web Visit Date", Format$(Now, "Medium Date")
End Sub
```

```
Private Sub cboAddress_Click()
    If mbDontNavigateNow Then Exit Sub
    tmTimer.Enabled = True
    brwWebBrowser.Navigate cboAddress.Text
    Me.MousePointer = vbHourglass
End Sub
```

```
Private Sub cboAddress_KeyPress(KeyAscii As Integer)
    On Error Resume Next
    If KeyAscii = vbKeyReturn Then cboAddress_Click
End Sub
```

```
Private Sub Form_Resize()
    Me.Refresh
    If Me.WindowState = vbMinimized Then Exit Sub
    brwWebBrowser.Move brwWebBrowser.Left, brwWebBrowser.Top, Me.ScaleWidth - 100, Me.ScaleHeight - (pnlAddress.Top +
    pnlAddress.Height) - 100
    brwWebBrowser.Width = Me.ScaleWidth - 100
    brwWebBrowser.Height = Me.ScaleHeight - (pnlAddress.Top + pnlAddress.Height) - 100
    cboAddress.Move cboAddress.Left, cboAddress.Top, pnlAddress.Width - cboAddress.Left - 100
End Sub
```

```
Private Sub tmTimer_Timer()
    If brwWebBrowser.Busy = False Then
        tmTimer.Enabled = False
        Me.Caption = brwWebBrowser.LocationName
    Else
        Me.Caption = "Locating Web Site..."
    End If
End Sub
```

```
Private Sub tbToolBar_ButtonClick(ByVal Button As Button)
    On Error Resume Next
    tmTimer.Enabled = True
    Select Case Button.Key
        Case "Back"
            brwWebBrowser.GoBack
        Case "Forward"
            brwWebBrowser.GoForward
        Case "Refresh"
            brwWebBrowser.Refresh
        Case "Home"
            brwWebBrowser.GoHome 'normally takes browser to the registered home page
            cboAddress = StartingAddress
            try to navigate to the starting address
            tmTimer.Enabled = True
            brwWebBrowser.Navigate StartingAddress
            Me.MousePointer = vbHourglass
        Case "Search"
            brwWebBrowser.GoSearch
        Case "Stop"
    End Select
End Sub
```

frmBrowser.frm - lbToolBar_ButtonClick

97

```
timTimer.Enabled = False  
brwWebBrowser.Stop  
Me.Caption = brwWebBrowser.LocationName  
End Select  
End Sub
```

frmTip.frm - File Declarations

98

```

Attribute VB_Name = "frmTip"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit

```

```

' The in-memory database of tips.
Dim Tips As New Collection

```

```

' Name of tips file
Const TIP_FILE = "TIPOFDAY.TXT"

```

Private Sub DoNextTip()

```

' Index in collection of tip currently being displayed.
' Cycle through the tips in order

```

```

giCurrentTip = giCurrentTip + 1
If Tips.Count < giCurrentTip Then giCurrentTip = 1

```

```

' Show it.
frmTip.DisplayCurrentTip
End Sub

```

Function LoadTips(sFile As String) As Boolean

```

Dim NextTip As String ' Each tip read in from file.
Dim InFile As Integer ' Descriptor for file.

```

```

' Obtain the next free file descriptor.
InFile = FreeFile

```

```

' Make sure a file is specified.

```

```

If sFile = "" Then
    LoadTips = False
    Exit Function
End If

```

```

' Make sure the file exists before trying to open it.

```

```

If Dir(sFile) = "" Then
    LoadTips = False
    Exit Function
End If

```

```

' Read the collection from a text file.

```

```

Open sFile For Input As InFile
While Not EOF(InFile)
    Line Input #InFile, NextTip
    Tips.Add NextTip
Wend
Close InFile

```

```

' Display a tip at random.
DoNextTip

```

```

LoadTips = True

```

```

End Function

```

frmAllPatients.frm - File Declarations

100

```

Attribute VB_Name = "frmAllPatients"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit
Dim xbAllPatientsFormLoading As Boolean
Dim xsPatientFileSpecs() As String 'a dynamic array holding the names of patient files on disk

```

Private Function CalculateSinglePatientCompliance(DataStruct As DeviceDataStruct) As Double

*'Calculate the compliance for the patient in memory and pass result back to caller
'Use the settings of the dialog to determine calculations and date ranges.*

```

Dim iDateBegin As Long, iDateEnd As Long, i As Long, iScoreSum As Long
Dim iPlotValue As Integer 'keeps a tally of the value to be plotted for each day

```

```

iDateBegin = txtStartDate 'Int(DataStruct.dEventDate(1))
iDateEnd = txtEndDate 'Int(DataStruct.dEventDate(DataStruct.iEventData(0)))

```

```

If iDateBegin Then 'there is at least a value in there
    Select Case cmbDataToView.ListIndex
    Case 0 'Doses per day score (all doses on this day regardless of time taken)
        For i = iDateBegin To iDateEnd
            iPlotValue = CalcDayDoseScore_AllDoses(DataStruct, i)
            iScoreSum = iScoreSum + iPlotValue
        Next i
        CalculateSinglePatientCompliance = iScoreSum / (iDateEnd - iDateBegin + 1)

    Case 1 'Compliance Doses per day (on time doses per day)
        For i = iDateBegin To iDateEnd
            iPlotValue = CalcDayDoseScore_OnTime(DataStruct, i)
            iScoreSum = iScoreSum + iPlotValue
        Next i
        CalculateSinglePatientCompliance = iScoreSum / (iDateEnd - iDateBegin + 1)

    Case 2 'Doses Taken
        For i = iDateBegin To iDateEnd
            iPlotValue = CalcDosesSumTakenOnSpecificDay(DataStruct, i)
            iScoreSum = iScoreSum + iPlotValue
        Next i
        CalculateSinglePatientCompliance = iScoreSum / (iDateEnd - iDateBegin + 1)
    End Select
End If
End Function

```

Private Sub btnClose_Click()

```

Unload Me
End Sub

```


frmAllPatients.frm - cmbaDataToView_Click

101

```

Private Sub cmbaDataToView_Click()
    CalculateAllPatientsComplianceOnDisk
    Slider1_SlideChange
End Sub

```

```

Private Sub cmbaDateSelection_Click()

```

```

    Select Case cmbaDateSelection.ListIndex
    Case 0 'recent 7 days
        txtEndDate = CDate(Int(Now))
        txtStartDate = CDate(Int(Now) - 7)

    Case 1 'recent 14 days
        txtEndDate = CDate(Int(Now))
        txtStartDate = CDate(Int(Now) - 14)

    Case 2 'recent 30 days
        txtEndDate = CDate(Int(Now))
        txtStartDate = CDate(Int(Now) - 30)

    Case 3 'recent 6 months
        txtEndDate = CDate(Int(Now))
        txtStartDate = CDate(Int(Now) - 180)

    Case 4 'all data
        If PAT_DATA.IEventDate(0) Then 'there are some events in array
            txtStartDate = CDate(Int(PAT_DATA.dEventDate(1)))
            txtEndDate = CDate(Int(PAT_DATA.dEventDate(PAT_DATA.IEventDate(0))))
        End If

    Case 5 'custom dates
        If gsLastStartDateChosen = "" Then
            txtStartDate = CDate(Int(PAT_DATA.dEventDate(1)))
        Else
            txtStartDate = gsLastStartDateChosen
        End If

        If gsLastEndDateChosen = "" Then
            txtEndDate = CDate(Int(PAT_DATA.dEventDate(PAT_DATA.IEventDate(0))))
        Else
            txtEndDate = gsLastEndDateChosen
        End If

    End Select

    CalculateAllPatientsComplianceOnDisk
    Slider1_SlideChange
End Sub

```

frmAllPatients.frm - Form_Activate

102

Private Sub Form_Activate()

```

Me.Refresh
grid.Refresh
Slider1_SlideChange
SetPrinterIcon True, "&Print All Patient's Summary..."
End Sub

```

Public Sub CalculateAllPatientsComplianceOnDisk()

*This procedure is called when it is necessary to update the display
due to some element or feature being changed.*

*Look at all device data files in the specified directory.
Retrieve appropriate data from each file and put into a global
structure holding all patients.*

```
If xbAllPatientsFormLoading = True Then Exit Sub
```

```

Dim r As Integer, l As Integer, dCompliance As Double
Dim sPath As String, sFileName As String, sFileSpec As String
Dim sTab As String, sTemp As String, lErrorCode As Long

```

```

On Error GoTo CalculateAllPatientsComplianceOnDisk_error
Me.MousePointer = vbHourglass
Me.Refresh

```

```

sTab = Chr$(9)
ReDim xsPatientFileSpecs(1) 'clear out the old array
grid.Clear
grid.Rows = 1

```

```

If cmbaDataToView.ListIndex = 2 Then
    'List item number 3 (index 2) was requested to be taken out. The code is still in the program
    'in case any iteration of it is needed later on.

```

```

    grid.FormatString = "< Patient Name |< Patient ID |< Start Date |< Last Dose |> Doses "
Else
    grid.FormatString = "< Patient Name |< Patient ID |< Start Date |< Last Dose |> Score "
End If

```

```
Form_Resize
```

```

grid.Col = 1 'set to column 1
grid.Redraw = False 'turn off redraw to speed up processing

```

```

sPath = App.Path + "\Patient Data\"
sFileName = LCase$(Dir$(sPath + "*.cpd")) 'get all filenames
Do While sFileName <> "" 'read all strings from directory
    sFileSpec = sPath + sFileName
    l = l + 1
    'Load the data for this patient into global array
    r = GetPatientDataFromDisk(sFileSpec, TEMP_DATA, lErrorCode)
    'ugh xxx it
    'If a checksum error or other error occurred on in the above function.
    'don't include the file in the summary and warn user.

```

```

'Call routine to calculate compliance based on dialog settings
dCompliance = CalculateSinglePatientCompliance(TEMP_DATA)

```

```
'Put results into grid
```

```
sTemp = TEMP_DATA.sPatientLastName + ", " + TEMP_DATA.sPatientFirstName + sTab + TEMP_DATA.sPatientID + sTab 'get name and ID
```

```
sTemp = sTemp + Format$(TEMP_DATA.dEventDate(1), "Short Date") + sTab 'get first dose date
```

```
sTemp = sTemp + Format$(TEMP_DATA.dEventDate(TEMP_DATA.lEventDate(0)), "Short Date") 'get last dose date
```

```
sTemp = sTemp + sTab + Format$(CStr(dCompliance), "#0") 'get compliance
```

```
If cmbaDataToView.ListIndex <> 2 Then sTemp = sTemp + " %"
```

```
grid.AddItem sTemp
```

\\AllPatients.frm) - CalculateAllPatientsComplianceOnDisk

103

```

grid.RowData() = 1 - 1
If i >= UBound(xsPatientFileSpecs) Then ReDim Preserve xsPatientFileSpecs(i + 10)

xsPatientFileSpecs(i + 1) = sFileSpec      'keep the name of the file here for when user clicks on cell
sFileName = LCase$(Dir)                  'get next file (one by one)
grid.Redraw = True
Loop                                     'process next filename

grid.Redraw = True

CalculateAllPatientsComplianceOnDisk_Exit:
Me.MousePointer = vbDefault
Exit Sub

CalculateAllPatientsComplianceOnDisk_error:
' Resume 0 'testing only
Resume CalculateAllPatientsComplianceOnDisk_Exit
End Sub

```

```

Private Sub Form_Load()
    frmMain.MousePointer = vbHourglass
    DoEvents

    xbAllPatientsFormLoading = True
    ReDim xsPatientFileSpecs(2)
    If cmbDateSelection.ListIndex < 0 Then cmbDateSelection.ListIndex = 2 'set a default

    cmbDataToView.ListIndex = 0

    xbAllPatientsFormLoading = False
    CalculateAllPatientsComplianceOnDisk
    grid.Col = 0
    grid.Sort = 1 'generic ascending
    frmMain.MousePointer = vbDefault
    RefreshAllOpenForms
End Sub

```

```

Private Sub Form_Resize()
    Dim iWidthRemaining As Integer
    Static bProcedureInProgress As Boolean
    If bProcedureInProgress Then Exit Sub
    If Me.WindowState = vbMinimized Then Exit Sub
    bProcedureInProgress = True

    If Me.Width < 5000 Then
        Me.Width = 5000
        bProcedureInProgress = False
    End If

    If Me.Height < 5000 Then
        Me.Height = 5000
        bProcedureInProgress = False
    End If

    SSPanel1.Left = Me.Width - SSPanel1.Width - 100
    grid.Width = SSPanel1.Left - grid.Left - 150
    grid.Height = Me.Height - grid.Top - 425

    grid.ColWidth(4) = 625 'Score
    grid.ColWidth(0) = (grid.Width - grid.ColWidth(4)) / (grid.Cols - 1) 'Name
    grid.ColWidth(1) = grid.ColWidth(0) 'ID
    iWidthRemaining = grid.Width - grid.ColWidth(0) - grid.ColWidth(1) - grid.ColWidth(4) - 120
    grid.ColWidth(2) = iWidthRemaining / 2 'Start Date
    grid.ColWidth(3) = grid.ColWidth(2) 'Last Dose Date

```

frmAllPatients.frm - Form_Resize

104

```

bProcedureInProgress = False
End Sub

```

Private Sub grid_Click()

```

Dim iRow As Integer

'Find out which column was clicked
'Sort the array only if the header was clicked

If grid.Rows < 2 Then Exit Sub

iRow = grid.MouseRow
If iRow = 0 Then
    grid.Col = grid.MouseCol
    grid.Sort = 1 'genenc ascending
    Exit Sub
End If

End Sub

```

Private Sub grid_DbClick()

```

Dim sFileName As String, r As Integer, iRow As Integer
If grid.Rows < 2 Then Exit Sub

iRow = grid.MouseRow
sFileName = xsPatientFileSpecs(grid.RowData(iRow))

'open the document that was double-clicked
r = OpenPatientData(sFileName)
End Sub

```

Private Sub Slider1_SlideChange()

```

Dim i As Integer, j As Integer

Label1(0).Caption = "Compliance Threshold = < " + CStr(Slider1.Value) + "%"
grid.Redraw = False

For i = 1 To grid.Rows - 1
    If gtAllPatients.sScore(grid.RowData(i)) < Slider1.Value Then
        grid.Row = i
        grid.Col = 4
        If grid.Value < Slider1.Value Then
            For j = 0 To grid.Cols - 1
                grid.Col = j
                grid.CellBackColor = &HC0FFFF '3HC0C0FF
            Next j
        Else
            For j = 0 To grid.Cols - 1
                grid.Col = j
                grid.CellBackColor = 0
            Next j
        End If
    End If
Next i

grid.Row = 0
grid.Col = 0
grid.Redraw = True
End Sub

```

frmAllPatients.frm - Slider1_SlideChange

105

```
Private Sub txtEndDate_HideDropDown()  
    DoEvents  
    gsLastEndDateChosen = txtEndDate  
    cmboDateSelection.ListIndex = 5      'select the custom setting  
    CalculateAllPatientsComplianceOnDisk  
    Slider1_SlideChange  
End Sub
```

```
Private Sub txtStartDate_HideDropDown()  
    DoEvents  
    gsLastStartDateChosen = txtStartDate  
    cmboDateSelection.ListIndex = 5      'select the custom setting  
    CalculateAllPatientsComplianceOnDisk  
    Slider1_SlideChange  
End Sub
```

frmRecentDosingGraph.frm - File Declaration

106

```
Attribute VB_Name = "frmPatientSummary"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit
```

Public Sub UpDatefrmPatientSummaryHeader()

```
Dim DataStruct As DeviceDataStruct
DataStruct = PAT_DATA
```

```
txtPatientLastName = "" + DataStruct.sPatientLastName
txtPatientFirstName = "" + DataStruct.sPatientFirstName
txtPatientID = "" + DataStruct.sPatientID
txtDrug = "" + DataStruct.sDrug
txtTxCenter = "" + DataStruct.sTxCenter
txtOrgan = "" + DataStruct.sOrgan
```

```
If DataStruct.dLastDownloadDate Then
```

```
txtRetrievalDate = "" + Format$(DataStruct.dLastDownloadDate, gsDateDisplayFormat)
```

```
Else
txtRetrievalDate = ""
```

```
End If
```

```
txtSerialNumber = "" + DataStruct.sSerialNumber
```

```
MSChart1.Visible = True
```

```
If cmbDateSelection.ListIndex < 0 Then cmbDateSelection.ListIndex = 1 'pick a default range
```

```
txtStartDate = Format$(CDate(DataStruct.dEventDate(1)), gsDateDisplayFormat)
```

```
txtEndDate = Format$(CDate(DataStruct.dEventDate(DataStruct.iEventData(0))), gsDateDisplayFormat)
```

```
Me.Refresh
```

```
End Sub
```

Public Sub UpdatePatientDosingGraph()

```
'Update the graph due to a check box being changed.
```

```
Dim sTab As String, I As Long
```

```
Dim IDateBegin As Long, IDateEnd As Long, IScoreSum As Long
```

```
Dim IPlotValue As Integer 'keeps a tally of the value to be plotted for each day
```

```
Dim dPlotDate As Double, IDayEventsFound As Integer, IDateIndex As Integer
```

```
Me.MousePointer = vbHourglass
```

```
DoEvents
```

```
If IsDate(txtStartDate) Then
```

```
IDateBegin = Int(CDate(txtStartDate))
```

```
Else
```

```
IDateBegin = Int(PAT_DATA.dEventDate(1))
```

```
End If
```

```
If IsDate(txtEndDate) Then
```

```
IDateEnd = Int(CDate(txtEndDate))
```

```
Else
```

```
IDateEnd = Int(PAT_DATA.dEventDate(PAT_DATA.iEventData(0)))
```

```
End If
```

```
MSChart1.RowCount = 0
```

```
dPlotDate = Int(PAT_DATA.dEventDate(1)) 'get the first event date
```

```
If IDateBegin Then 'there is at least a value in there
```

```
Select Case cmbDataToView.ListIndex
```

```
Case 0 'Doses per day score (all doses on this day regardless of time taken
```

```
For I = IDateBegin To IDateEnd 'the number of events is stored here
```

```
IPlotValue = CalcDayDoseScore_AllDoses(PAT_DATA, I)
```

```
MSChart1.RowCount = MSChart1.RowCount + 1 'increment the row count
```

```
MSChart1.Row = MSChart1.RowCount
```

```
MSChart1.Data = IPlotValue 'plot in last row
```

```
IScoreSum = IScoreSum + IPlotValue
```

nRecentDosingGraph.frm - UpdatePatientDosem aph

107

```

Next I
MSChart1.Plot.Axes(VtChAxisIdY).AxisScale.Type = VtChScaleTypePercent * 100 'set scale to percent
MSChart1.Plot.Axes(VtChAxisIdY).AxisScale.PercentBasis = VtChPercentAxisBasisMaxChart
MSChart1.Plot.Axes(VtChAxisIdY).AxisScale.Type = VtChScaleTypeLinear
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.Auto = False
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.Maximum = 100
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.MajorDivision = 10
MSChart1.TitleText = "Doses/Day Score (regardless of prescribed time)"
MSChart1.Plot.Axes(VtChAxisIdY).AxisTitle.Text = "Percent"
txtScore = Format$(CStr((ScoreSum / (IDateEnd - IDateBegin + 1)), "##") + " %" 'put the most recent score in the text box

Case 1 'Compliance Doses per day (on time doses per day)
For I = IDateBegin To IDateEnd 'the number of events is stored here
    iPlotValue = CalcDayDoseScore_OnTime(PAT_DATA, I)
    MSChart1.RowCount = MSChart1.RowCount + 1 'increment the row count
    MSChart1.Row = MSChart1.RowCount 'plot in last row
    MSChart1.Data = iPlotValue
    iScoreSum = iScoreSum + iPlotValue
Next I
MSChart1.Plot.Axes(VtChAxisIdY).AxisScale.Type = VtChScaleTypeLinear
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.Auto = False
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.Maximum = 100
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.MajorDivision = 10
MSChart1.TitleText = "On-Time Doses/Day Score (within prescribed time)"
MSChart1.Plot.Axes(VtChAxisIdY).AxisTitle.Text = "Percent"
txtScore = Format$(CStr((iScoreSum / (IDateEnd - IDateBegin + 1)), "##") + " %" 'put the most recent score in the text box

Case 2 'Doses Taken per day
For I = IDateBegin To IDateEnd 'the number of events is stored here
    iPlotValue = CalcDosesSumTakenOnSpecificDay(PAT_DATA, I)
    MSChart1.RowCount = MSChart1.RowCount + 1 'increment the row count
    MSChart1.Row = MSChart1.RowCount 'plot in last row
    MSChart1.Data = iPlotValue
Next I
MSChart1.Plot.Axes(VtChAxisIdY).AxisScale.Type = VtChScaleTypeLinear
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.Auto = False
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.Maximum = 10 'client wants to hard code this at 10
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.MajorDivision = 10
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.MinorDivision = 1
MSChart1.TitleText = "Total Doses Taken Per Day"
MSChart1.Plot.Axes(VtChAxisIdY).AxisTitle.Text = "Doses"
txtScore = "" 'no score appears for doses sum

Case 3 'Doses Missed
'note, this section is not used currently. Client decided to remove it from the display
For I = IDateBegin To IDateEnd 'the number of events is stored here
    iPlotValue = PAT_DATA.iDosesPerDay - CalcDosesSumTakenOnSpecificDay(PAT_DATA, I)
    MSChart1.RowCount = MSChart1.RowCount + 1 'increment the row count
    MSChart1.Row = MSChart1.RowCount 'plot in last row
    MSChart1.Data = iPlotValue
Next I
MSChart1.Plot.Axes(VtChAxisIdY).AxisScale.Type = VtChScaleTypeLinear
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.Auto = False
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.Maximum = 10 'client wants to hard code this at 10
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.MajorDivision = 10
MSChart1.Plot.Axes(VtChAxisIdY).ValueScale.MinorDivision = 1
MSChart1.Plot.Axes(VtChAxisIdY).AxisScale.Hide = False
MSChart1.TitleText = "Doses Missed Per Day"
MSChart1.Plot.Axes(VtChAxisIdY).AxisTitle.Text = "Doses"
txtScore = "" 'no score appears for doses missed
End Select
End If

MSChart1.xaxisid = "ies"
Me.MousePointer = vbDefault

End Sub

```

frmRecentDosingGraph.frm - btnClose_Click

108

```
Private Sub btnClose_Click()
    Unload Me
End Sub
```

```
Private Sub cmboAverageDays_Click()
    UpdatePatientDosingGraph
End Sub
```

```
Private Sub cmboChartType_Click()
    Select Case cmboChartType.Text
        Case "Line"
            MSChart1.chartType = VtChChartType2dLine
        Case "Area"
            MSChart1.chartType = VtChChartType2dArea
        Case "Bar"
            MSChart1.chartType = VtChChartType2dCombination
        Case "Step"
            MSChart1.chartType = VtChChartType2dStep
    End Select
End Sub
```

```
Private Sub cmbaDataToView_Click()
    UpdatePatientDosingGraph
End Sub
```

```
Public Sub cmboDateSelection_Click()
```

```
    If PAT_DATA.iEventData(0) = 0 Then 'no data appears to be loaded
        txtStartDate = CDate(Now)
        txtEndDate = CDate(Now)
        Exit Sub
    End If
```

```
    Select Case cmboDateSelection.ListIndex
        Case 0 'recent 7 days
            txtEndDate = CDate(Int(PAT_DATA.dEventDate(PAT_DATA.iEventData(0))))
            txtStartDate = CDate(Int(PAT_DATA.dEventDate(PAT_DATA.iEventData(0))) - 7)
        Case 1 'recent 14 days
            txtEndDate = CDate(Int(PAT_DATA.dEventDate(PAT_DATA.iEventData(0))))
            txtStartDate = CDate(Int(PAT_DATA.dEventDate(PAT_DATA.iEventData(0))) - 14)
        Case 2 'recent 30 days
            txtEndDate = CDate(Int(PAT_DATA.dEventDate(PAT_DATA.iEventData(0))))
            txtStartDate = CDate(Int(PAT_DATA.dEventDate(PAT_DATA.iEventData(0))) - 30)
        Case 3 'recent 6 months
            txtEndDate = CDate(Int(PAT_DATA.dEventDate(PAT_DATA.iEventData(0))))
            txtStartDate = CDate(Int(PAT_DATA.dEventDate(PAT_DATA.iEventData(0))) - 180)
```


frmRecentDosingGraph.frm - Form_Resiz

110

```

End If

pnlControls.Left = Me.Width - pnlControls.Width - 200
pnlChart.Width = pnlControls.Left - 100
MSChart1.Width = pnlChart.Width - 100
pnlChart.Height = Me.Height - pnlChart.Top - 500
MSChart1.Height = pnlChart.Height - 100

bProcedureInProgress = False

End Sub

Private Sub Form_Unload(Cancel As Integer)
    'Save last settings selected by user
    PAT_SUM_DEFAULTS.cmboDataToView = cmboDataToView.ListIndex
    PAT_SUM_DEFAULTS.cmboChartType = cmboChartType.ListIndex
End Sub

```

```

Private Sub txtEndDate_Change()
    DoEvents
    gsLastEndDateChosen = txtEndDate
    'In case the user chose a date far removed from the date of the first
    'dose, then notify the user and set the date to the first dose.
    If DateValue(gsLastEndDateChosen) > Int(PAT_DATA.dEventData(PAT_DATA.iEventData(0))) Then
        MsgBox "The ending date you chose is later than the last dose taken by this patient. The starting date is being set to the first dose.",
            vbInformation, "Ending Date Too Late"
        txtEndDate = CDate(PAT_DATA.dEventData(PAT_DATA.iEventData(0)))
        gsLastEndDateChosen = txtEndDate
    End If

    cmboDataSelection.ListIndex = 5 'select the user setting
    UpdatePatientDosingGraph
End Sub

```

```

Private Sub txtStartDate_Change()
    DoEvents
    gsLastStartDateChosen = txtStartDate
    'In case the user chose a date far removed from the date of the last
    'dose, then notify the user and set the date to the last dose.
    If DateValue(gsLastStartDateChosen) < Int(PAT_DATA.dEventData(1)) Then
        MsgBox "The starting date you chose is sooner than the first dose taken by this patient. The starting date is being set to the first dose",
            vbInformation, "Start Date Too Early"
        txtStartDate = CDate(PAT_DATA.dEventData(1))
        gsLastStartDateChosen = txtStartDate
    End If

    cmboDataSelection.ListIndex = 5 'select the user setting
    UpdatePatientDosingGraph
End Sub

```

125

frmDosingCalendar.frm - File Declarations

111

```

Attribute VB_Name = "frmDosingCalendar"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit

```

```

Private bgResizedCalendar As Boolean

```

```

Private Sub btnChangeCompliance_Click()
    glLatestOptionsTabSelected = 1 'display the proper tab once the dialog is opened
    frmOptions.Show vbModal

```

```

End Sub

```

```

Private Sub btnClose_Click()
    Unload Me
End Sub

```

```

Private Sub Calendar_DayChange()
    Static bProcedureInProgress As Boolean

    If bProcedureInProgress Then Exit Sub 'prevent recursive calls
    bProcedureInProgress = True
    frmDosingCalendar.MousePointer = vbHourglass 'hour glass
    DoEvents
    UpdateZoomBox
    frmDosingCalendar.MousePointer = vbDefault 'default glass
    bProcedureInProgress = False 'allow another call to this sub
End Sub

```

```

Private Sub Calendar_MonthChange()
    Static bProcedureInProgress

    If bProcedureInProgress Then Exit Sub 'prevent recursive calls
    bProcedureInProgress = True
    frmDosingCalendar.MousePointer = vbHourglass 'hour glass
    DoEvents
    UpdateCalendar
    frmDosingCalendar.MousePointer = vbDefault 'default
    bProcedureInProgress = False
End Sub

```

frmDosingCalendar.frm - Calendar_Mouse

112

```
Private Sub Calendar_MouseMove(Button As Integer, Shift As Integer, X As Single, Y As Single)
    If bgResizedCalendar Then
        UpdateCalendar
        bgResizedCalendar = False
    End If
End Sub
```

```
Private Sub Calendar_YearChange()
    Calendar_MonthChange
End Sub
```

```
Private Sub chkDoseChanged_Click()
    DrawAllDoseSizeChanges
    UpdateZoomBox
End Sub
```

```
Private Sub chkDosesMissed_Click()
    DrawAllDosesMissed
    UpdateZoomBox
End Sub
```

```
Private Sub chkDosesNotComplied_Click()
    DrawAllNonCompliedDosesTaken
    UpdateZoomBox
End Sub
```

```
Private Sub chkDosesTaken_Click()
    DrawAllCompliedDosesTaken
    UpdateZoomBox
End Sub
```

```
Private Sub chkWeekNumbers_Click()
    RemoveAllObjects
    frmDosingCalendar.Calendar.WeekNumbers = chkWeekNumbers
    DoEvents
    UpdateCalendar
End Sub
```

frmDosingCalendar.frm - Form_Activate

113

```
Private Sub Form_Activate()
    SetPrinterIcon False, -
End Sub
```

```
Private Sub Form_Load()
    Dim i As Integer
    'start calendar with date of latest dose
    If PAT_DATA.dEventData(PAT_DATA.iEventData(0)) > 0 Then
        frmDosingCalendar.Calendar.Date = CDate(PAT_DATA.dEventData(PAT_DATA.iEventData(0)))
    Else
        frmDosingCalendar.Calendar.Date = Now
    End If
    Me.Show

    Load lblDetailTime(2)
    Load lblDetailTime(4)
    Load lblDetailTime(6)
    Load lblDetailTime(8)
    Load lblDetailTime(10)
    Load lblDetailTime(12)
    Load lblDetailTime(14)
    Load lblDetailTime(16)
    Load lblDetailTime(18)
    Load lblDetailTime(20)
    Load lblDetailTime(22)
    CreateCalendarTimeScale
    DoEvents
    frmDosingCalendar.Calendar.MouseExpand = 5    'expand the hot spot around date arrows
    Me.Show

    'Set the dialog controls to the settings last set by user
    chkDosesMissed = CAL_DEFAULTS.chkDosesMissed
    chkDosesNotComplied = CAL_DEFAULTS.chkDosesNotComplied
    chkDosesTaken = CAL_DEFAULTS.chkDosesTaken
    chkDoseChanged = CAL_DEFAULTS.chkDoseChanged

    UpdateCalendar
End Sub
```

```
Private Sub CreateCalendarTimeScale()
    On Error Resume Next
    'Create the time scale on detail area
    Dim sAM As String, sPM As String, i As Integer

    If frmDosingCalendar.Width < 5000 Then
        sAM = ""
        sPM = ""
    Else
        sAM = "am"
        sPM = "pm"
    End If

    lblDetailTime(2).Caption = "2" + sAM
    With lblDetailTime(2)
        .Left = (Me.priZoom.Width * (1 / 24)) - (Me.lblDetailTime(i).Width / 2)
        .ForeColor = &HFFFFFF
        .Visible = True
        .ZOrder
    End With

    lblDetailTime(4).Caption = "4" + sAM
    With lblDetailTime(4)
        .Left = (Me.priZoom.Width * (1 / 24)) - (Me.lblDetailTime(i).Width / 2)
        .ForeColor = &HFFFFFF
        .Visible = True
    End With
```

frmDosingCalendar.frm - CreateCalendarTimeS

114

```

.ZOrder
End With

lblDetailTime(6).Caption = "6" + sAM
With lblDetailTime(6)
.Left = (Me.pnlZoom.Width * (1 / 24)) - (Me.lblDetailTime(i).Width / 2)
.ForeColor = &HFFFFFF
.Visible = True
.ZOrder
End With

lblDetailTime(8).Caption = "8" + sAM
With lblDetailTime(8)
.Left = (Me.pnlZoom.Width * (1 / 24)) - (Me.lblDetailTime(i).Width / 2)
.ForeColor = &HFFFFFF
.Visible = True
.ZOrder
End With

lblDetailTime(10).Caption = "10" + sAM
With lblDetailTime(10)
.Left = (Me.pnlZoom.Width * (1 / 24)) - (Me.lblDetailTime(i).Width / 2)
.ForeColor = &HFFFFFF
.Visible = True
.ZOrder
End With

lblDetailTime(12).Caption = "12" + sPM
With lblDetailTime(12)
.Left = (Me.pnlZoom.Width * (1 / 24)) - (Me.lblDetailTime(i).Width / 2)
.ForeColor = &HFFFFFF
.Visible = True
.ZOrder
End With

lblDetailTime(14).Caption = "2" + sPM
With lblDetailTime(14)
.Left = (Me.pnlZoom.Width * (1 / 24)) - (Me.lblDetailTime(i).Width / 2)
.ForeColor = &HFFFFFF
.Visible = True
.ZOrder
End With

lblDetailTime(16).Caption = "4" + sPM
With lblDetailTime(16)
.Left = (Me.pnlZoom.Width * (1 / 24)) - (Me.lblDetailTime(i).Width / 2)
.ForeColor = &HFFFFFF
.Visible = True
.ZOrder
End With

lblDetailTime(18).Caption = "6" + sPM
With lblDetailTime(18)
.Left = (Me.pnlZoom.Width * (1 / 24)) - (Me.lblDetailTime(i).Width / 2)
.ForeColor = &HFFFFFF
.Visible = True
.ZOrder
End With

lblDetailTime(20).Caption = "8" + sPM
With lblDetailTime(20)
.Left = (Me.pnlZoom.Width * (1 / 24)) - (Me.lblDetailTime(i).Width / 2)
.ForeColor = &HFFFFFF
.Visible = True
.ZOrder
End With

```

frmDosingCalendar.frm - CreateCalendarTimeScale

115

```

ibiDetailTime(22).Caption = "10" + sPM
With ibiDetailTime(22)
    .Left = (Me.pniZoom.Width * (i / 24)) - (Me.ibiDetailTime(i).Width / 2)
    .ForeColor = &HFFFFFF
    .Visible = True
    .ZOrder
End With

On Error GoTo 0
End Sub

```

```

Private Sub Form_MouseMove(Button As Integer, Shift As Integer, X As Single, Y As Single)
    If bgResizedCalendar Then
        UpdateCalendar
        bgResizedCalendar = False
    End If
End Sub

```

```

Private Sub Form_Resize()

```

```

    Me.Refresh
    Static bProcedureInProgress As Boolean

    If bProcedureInProgress Then Exit Sub
    If Me.WindowState = vbMinimized Then Exit Sub

```

```

    bProcedureInProgress = True

```

```

    If Me.Width < 6000 Then
        Me.Width = 6000
        bProcedureInProgress = False
    End If

```

```

    If Me.Height < 5000 Then
        Me.Height = 5000
        bProcedureInProgress = False
    End If

```

```

    DeleteAllObjects

```

```

    CreateCalendarTimeScale
    pniControls.Left = Me.Width - pniControls.Width - 200
    Calendar.Width = pniControls.Left - Calendar.Left - 150

```

```

    pniZoom.Top = Me.Height - pniZoom.Height - 450
    pniTime.Top = pniZoom.Top

```

```

    Calendar.Height = pniZoom.Top - Calendar.Top - 100
    bgResizedCalendar = True 'med events on form can not be updated until resize is done
    pniZoom.Width = Calendar.Width

```

```

    bProcedureInProgress = False
End Sub

```

frmDosingCalendar.frm - Form_Unload

116

```
Private Sub Form_Unload(Cancel As Integer)
    CAL_DEFAULTS.chkDosesMissed = chkDosesMissed
    CAL_DEFAULTS.chkDosesNotComplied = chkDosesNotComplied
    CAL_DEFAULTS.chkDosesTaken = chkDosesTaken
    CAL_DEFAULTS.chkDoseChanged = chkDoseChanged
End Sub
```

```
Private Sub frameView_MouseMove(Button As Integer, Shift As Integer, X As Single, Y As Single)
    If bgResizedCalendar Then
        UpdateCalendar
        bgResizedCalendar = False
    End If
End Sub
```

```
Private Sub pnlControls_MouseMove(Button As Integer, Shift As Integer, X As Single, Y As Single)
    If bgResizedCalendar Then
        UpdateCalendar
        bgResizedCalendar = False
    End If
End Sub
```

frmPatientDosingRpt.frm - File Declaration

117

```
Attribute VB_Name = "frmPatientDosingReport"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit
```

Private Sub RescaleGrid()

```
Dim iRemainder As Integer
'put any fixed width columns first
grid.ColWidth(6) = 0 'don't show this column
grid.ColWidth(1) = 1300
grid.ColWidth(2) = 900
grid.ColWidth(3) = 1000
grid.ColWidth(4) = 950
```

```
iRemainder = grid.Width - grid.ColWidth(4) - grid.ColWidth(3) - grid.ColWidth(2) - grid.ColWidth(1)
grid.ColWidth(0) = iRemainder * 0.25
grid.ColWidth(5) = iRemainder * 0.75 - 370
```

```
End Sub
```

Public Sub UpdatefrmPatientDosingReportHeader()

```
'Show custom labels from config file if there were any
```

```
Label1(3) = gsCustomLblPatientLastName
Label1(1) = gsCustomLblPatientFirstName
Label1(2) = gsCustomLblOrgan
Label1(0) = gsCustomLblPatientID
Label1(6) = gsCustomLblTxCenter
Label1(7) = gsCustomLblDrug
```

```
txtPatientLastName = "" + PAT_DATA.sPatientLastName
txtPatientFirstName = "" + PAT_DATA.sPatientFirstName
txtPatientID = "" + PAT_DATA.sPatientID
txtDrug = "" + PAT_DATA.sDrug
txtTxCenter = "" + PAT_DATA.sTxCenter
txtOrgan = "" + PAT_DATA.sOrgan
```

```
End Sub
```

Public Sub UpdatePatientGridDisplay()

```
'This proc is called when it is necessary to update the display
'due to some element or feature being changed.
```

```
Dim sTab As String, i As Integer, sTime As String, sDate As String, sTemp As String
```

```
Dim bShowDosesTaken As Boolean, bShowDoseChanges As Boolean, bShowUserEvents As Boolean
```

```
sTab = Chr$(9)
```

```
grid.Clear 'the erase the grid
```

```
grid.Rows = 1
```

```
grid.FormatString = "<Date & Time | Event Type | Dose Size |" + gsLabelGridColumnCustom1 + " |" + gsLabelGridColumnCustom2  
+ " |" + gsLabelGridColumnCustom3
```

```
RescaleGrid
```

```
'grid.Row = 0
```

```
grid.Col = 1 'set to column 1
```

```
grid.Redraw = False
```

```
bShowDosesTaken = chkDoses.Value 'speed up the display in loop by assigning control value to a var
bShowDoseChanges = chkDoseChanged.Value
bShowUserEvents = chkUserDefined.Value
```

frmPatientDosingRpt.frm - UpdatePatientGridDi.

118

```

For i = 1 To CInt(PAT_DATA.iEventData(0))      'the number of events is stored here
If PAT_DATA.byteEventType(i) = giEVENT_DOSE_TAKEN Then
If bShowDosesTaken Then
sTemp = Format$(PAT_DATA.dEventData(i), gsDateDisplayFormat) + " " + Format$(PAT_DATA.dEventData(i),
gsTimeDisplayFormat)
sTemp = sTemp + sTab + "Dose Taken" + sTab + CStr(PAT_DATA.iEventData(i)) + " mg"
sTemp = sTemp + sTab + PAT_DATA.sUserData1(i) + sTab + PAT_DATA.sUserData2(i) + sTab + PAT_DATA.sUserData3(i)
grid.AddItem sTemp
'xxx grid.ListApplyTo = 12 'LC_LISTAPPLYTO_SINGLE_ITEM
'xxx grid.ForeColor = &H0 'black
grid.RowData(grid.Rows - 1) = CStr(i)
End If

ElseIf PAT_DATA.byteEventType(i) = giEVENT_DOSE_CHANGED Then
If bShowDoseChanges Then
sTemp = Format$(PAT_DATA.dEventData(i), gsDateDisplayFormat) + " " + Format$(PAT_DATA.dEventData(i),
gsTimeDisplayFormat)
sTemp = sTemp + sTab + "Dose Change" + sTab + CStr(PAT_DATA.iEventData(i)) + " mg"
sTemp = sTemp + sTab + PAT_DATA.sUserData1(i) + sTab + PAT_DATA.sUserData2(i) + sTab + PAT_DATA.sUserData3(i)
grid.AddItem sTemp
'xxx grid.ListApplyTo = 12 'LC_LISTAPPLYTO_SINGLE_ITEM
'xxx grid.ForeColor = &HC0FFFF 'yellow
grid.RowData(grid.Rows - 1) = i
End If

ElseIf PAT_DATA.byteEventType(i) = giEVENT_USER_DEFINED Then
If bShowUserEvents Then
sTemp = Format$(PAT_DATA.dEventData(i), gsDateDisplayFormat) + " " + Format$(PAT_DATA.dEventData(i),
gsTimeDisplayFormat)
sTemp = sTemp + sTab + "Custom Event" + sTab      'no doses data to be saved with custom events + sTab + CStr(PAT_DATA.
iEventData(i)) + " mg"
sTemp = sTemp + sTab + PAT_DATA.sUserData1(i) + sTab + PAT_DATA.sUserData2(i) + sTab + PAT_DATA.sUserData3(i)
grid.AddItem sTemp
'xxx grid.ListApplyTo = 12 'LC_LISTAPPLYTO_SINGLE_ITEM
'xxx grid.ForeColor = &HC0FFFF 'yellow
grid.RowData(grid.Rows - 1) = i
End If
End If

Next i
grid.Redraw = True
grid.Row = 0
grid.Col = 0
btnDeleteUserEvent.Enabled = False

End Sub

Private Sub btnClose_Click()
Unload Me
End Sub

```

frmPatientDosingRpt.frm - btnDeleteUserEvent

119

```
Private Sub btnDeleteUserEvent_Click()
    Dim r As Integer, lindex As Integer, sMSG As String
```

```
    sMSG = "The selected event will be permanently removed from the file." + vbCrLf + vbCrLf + "Do you want to delete the event?"
    r = MsgBox(sMSG, vbYesNo + vbQuestion, "Verify Event Deletion")
```

```
    If r = vbYes Then
```

```
        lindex = CInt(grid.RowData(grid.Row))
```

```
        If PAT_DATA.byteEventType(lindex) = giEVENT_USER_DEFINED Then EventDelete PAT_DATA, lindex
```

```
        event
```

```
        Call UpdatePatientGridDisplay
```

```
        'delete this
```

```
    End If
```

```
End Sub
```

```
Private Sub btnNewUserEvent_Click()
```

```
    'Add an event to the grid for a time and date defined by the user
    Dim lindex As Integer, lDate As Long, i As Integer
```

```
    frmGetDateTime.Show vbModal
```

```
    'Get the date of event from user
```

```
    If gdTempDateTime Then
```

```
        'don't add a date if user cancelled out of entry dialog
```

```
        'Find the date in the structure
```

```
        lindex = FindClosestDateInArray(PAT_DATA, gdTempDateTime)
```

```
        EventInsert PAT_DATA, lindex, gdTempDateTime
```

```
        'insert a new custom event
```

```
        If chkUserDefined.Value = 0 Then
```

```
            chkUserDefined.Value = vbChecked
```

```
        Else
```

```
            Call UpdatePatientGridDisplay
```

```
        End If
```

```
        grid.SetFocus
```

```
        grid.Col = 3
```

```
        For i = 1 To grid.Rows - 1
```

```
            If grid.RowData(i) = lindex Then
```

```
                'this is the index we just added
```

```
                grid.Row = i
```

```
                'highlight this row
```

```
                grid.TopRow = i
```

```
            End If
```

```
        Next i
```

```
    End If
```

```
End Sub
```

```
Private Sub chkDoseChanged_Click()
```

```
    Call UpdatePatientGridDisplay
```

```
End Sub
```

```
Private Sub chkDoses_Click()
```

```
    Call UpdatePatientGridDisplay
```

```
End Sub
```

frmPatientDosingRpt frm - chkUserDefined_

120

```
Private Sub chkUserDefined_Click()
    Call UpdatePatientGridDisplay
End Sub
```

```
Private Sub Form_Activate()
    Me.Refresh
    Form_Resize
    grid.Refresh
    Call UpdatePatientGridDisplay
    SetPrinterIcon True, "&Print Dosing Report..."
End Sub
```

```
Private Sub Form_Load()
    UpdatefrmPatientDosingReportHeader
End Sub
```

```
Private Sub Form_Resize()
    Static bProcedureInProgress As Boolean
    If bProcedureInProgress Then Exit Sub
    If Me.WindowState = vbMinimized Then Exit Sub
    bProcedureInProgress = True

    If Me.Width < 8100 Then
        Me.Width = 8100
        bProcedureInProgress = False
    End If

    If Me.Height < 5000 Then
        Me.Height = 5000
        bProcedureInProgress = False
    End If

    frameView.Left = Me.Width - frameView.Width - 250
    btnClose.Left = Me.Width - btnClose.Width - 250

    grid.Width = btnClose.Left + btnClose.Width - grid.Left
    grid.Height = Me.Height - grid.Top - 425
    RescaleGrid

    bProcedureInProgress = False
End Sub
```

```
Private Sub grid_AfterEdit(ByVal Row As Long, ByVal Col As Long)
```

```
    Dim Index As Integer
```

```
    Select Case Col
```

```
        Case 3 'user column 1
```

```
            PAT_DATA.sUserData1(grid.RowData(Row)) = grid.Text 'put the change into the structure
```

```
        Case 4 'user column 2
```

```
            PAT_DATA.sUserData2(grid.RowData(Row)) = grid.Text 'put the change into the structure
```

```
        Case 5 'user column 3
```

```
            PAT_DATA.sUserData3(grid.RowData(Row)) = grid.Text 'put the change into the structure
```

```
    End Select
```

```
    gbPatientDataNotSaved = True 'set flag to indicate that the file has changed but not yet been saved
```

```
    If grid.Col = 3 Then
```

```
        grid.Col = 4 'go to next cell
```

```
    ElseIf grid.Col = 4 Then
```

```
        grid.Col = 5 'go to next cell
```

frmPatientDosingRpt.frm - grid_AfterEdit

121

```

ElseIf grid.Col = 5 Then
    If grid.Rows > grid.Row + 1 Then 'prevent run time error by looking for last row
        grid.Row = grid.Row + 1
        grid.Col = 3 'gc to next cell
    End If
End If
End Sub

```

```

Private Sub grid_KeyDown(KeyCode As Integer, Shift As Integer)

```

```

    If KeyCode = 46 And grid.Col > 4 Then
        'delete key was pressed and column is editable
        grid.Text = ""
        gbPatientDataNotSaved = True 'set flag to indicate that the file has changed but not yet been saved
        grid_AfterEdit grid.Row, grid.Col
    End If
End Sub

```

```

Private Sub grid_RowColChange()
    Dim iIndex As Integer

```

```

    If grid.Col > 4 Then 'allow cell to be edited if it is a custom column
        grid.Editable = True
    Else
        grid.Editable = False
    End If

    iIndex = CInt(grid.RowData(grid.Row))
    If PAT_DATA.byteEventType(iIndex) = giEVENT_USER_DEFINED Then
        btnDeleteUserEvent.Enabled = True
    Else
        btnDeleteUserEvent.Enabled = False
    End If
End Sub

```

frmReadDeviceData.frm - File Declarations

122

```
Attribute VB_Name = "frmReadDeviceData"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit
```

Private Sub btnClose_Click()

```
Unload Me
End Sub
```

Private Sub btnReadEntireContents_Click()

```
Dim r As Integer, IErrorCode As Long, i As Integer, sMSG As String
```

```
r = ValidatePatientDataSaved 'ensure that previous patient data was saved before proceeding
If r = vbCancel Then Exit Sub
```

```
btnReadEntireContents.Enabled = False 'prevent recursive calls to device
gbKeepPollingDevice = False 'stop polling for now
Wait 0.25
```

```
txtPatientLastName = "" 'clear out the text boxes before reading data
txtPatientFirstName = "" 'clear out the text boxes before reading data
txtDrug = ""
txtPatientID = ""
txtTxCenter = ""
txtOrgan = ""
txtSerialNumber = ""
txtDoseSize = ""
txtDoseTime(1) = ""
txtDoseTime(2) = ""
txtDoseTime(3) = ""
txtDoseTime(4) = ""
txtDosesPerDay = ""
txtDoseLockoutHours = ""
txtMedicationRemaining = ""
txtEventCount = ""
txtLastRetrievalDate = ""
txtDeviceStarted = ""
```

```
r = Comm_ReadEntireMemoryContents(PAT_DATA, IErrorCode)
```

```
If r Then
```

```
PopulateDeviceCommDialog PAT_DATA, Me
PAT_DATA.sPatientDataFileName = ""
frmMain.mnuFileSave.Enabled = False
```

```
Else
```

```
EraseDataInMemory PAT_DATA
DisplayErrorMessage IErrorCode
```

```
End If
```

```
btnReadEntireContents.Enabled = True 're-enable button
RefreshAllOpenForms
```

```
'Compare battery time to value retrieved from ini file to determine if a
'reminder should be given to the user to change the batteries.
```

```
i = CInt(GetINISetting(gsAppIniFileSpec, "Options", "Battery Change Days", 180))
```

```
If i And Val(PAT_DATA.sBatteryChangeTimer) >= i Then
```

```
sMSG = "The battery in this device needs to be changed." + vbCrLf + vbCrLf
```

```
'Also look at the error flag returned from the device to see if the brownout
'bit was set. If so, append a different notice to the message than the normal one.
```

```
If PAT_DATA.bErrorBrownOut Then
```

```
sMSG = sMSG + "The device indicates that power was briefly lost due to low voltage."
```

```
Else
```

```
sMSG = sMSG + "They have been in place for over " + CStr(i / 30) + " months."
```

nReadDeviceData.frm - btnReadEntireContents_ rck

123

```

End If

sMSG = sMSG + vbCrLf + vbCrLf + "Do you want to change the battery now?"
Beep
r = MsgBox(sMSG, vbExclamation + vbYesNo + vbDefaultButton2, "Battery Change Needed")
If r = vbYes Then
    Call ChangeBatteriesRequest
End If

End If

End Sub

```

Private Sub Form_Activate()

```

PopulateDeviceCommDialog PAT_DATA, Me
Comm_InitializeCommPort 'initialize the comm port from ini file settings

gbKeepPollingDevice = True 'continue polling device
PollDeviceContinually Me
SetPrinterIcon False, ""
End Sub

```

Private Sub Form_Load()

```

Me.Left = 0
Me.Top = 0
Unload frmDeviceInitialize 'don't need this form
gbCommOK = 99 'reset flag that will give an indication as to the communication status.
End Sub

```

Private Sub Form_Unload(Cancel As Integer)

```

gbKeepPollingDevice = False 'stop polling the device
Wait 0.1
End Sub

```

frmPrint.frm - File Declarations

124

```

Attribute VB_Name = "frmPrint"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit

```

```

Private Sub btnClose_Click()
    Unload Me
End Sub

```

```

Private Sub btnPrintNow_Click()
    btnPrintNow.Enabled = False
    btnPrintNow.Refresh
    vsPrinter1.Action = paPrintAll
    btnPrintNow.Enabled = False
End Sub

```

```

Private Sub btnPrintPage_Click()
    btnPrintPage.Enabled = False
    btnPrintPage.Refresh
    vsPrinter1.Action = paPrintPage
    btnPrintPage.Enabled = True
End Sub

```

```

Private Sub btnRefresh_Click()
    RefreshPreview
End Sub

```

```

Private Sub Form_Load()
    Dim i As Integer

    gbPrintFormLoading = True
    frmPrint.vsViewPort1.BorderStyle = 1    'turn off til needed

    gbPrinterErrorDetected = False

    Me.Width = 7500
    Me.Move (Screen.Width - Me.Width) / 2, (Screen.Height - Me.Height) / 2    'center form on screen
    datPrintingredient.DatabaseName = sgDataBaseName

    lblActivePrinter.Caption = "" + vsPrinter1.Device
    gbPreventPreviewUpdates = False    'allow controls to update when called

    Me.Show
    DoEvents
    RefreshPreview
    SetPrinterIcon False, ""
End Sub

```

frmPrint.frm - Form_QueryUnload

125

Private Sub Form_QueryUnload(Cancel As Integer, UnloadMode As Integer)

Dim r As Integer

If gbPrintSpoolingInProgress Then *'user tried to exit while print spooling*

Beep

r = MsgBox("Items are still waiting to be printed. If you continue, the print job may be lost." + vbCrLf + vbCrLf + "Do you still want to close this form?", vbQuestion + vbYesNo, "Waiting For Printer")

If r = vbNo Then Cancel = True *'prevent crash error*

End If

gbPrintSpoolingInProgress = False

End Sub

Private Sub Form_Resize()
If frmPrint.WindowState <> vbMinimized Then *'not maximized*

Panel3D1.Left = frmPrint.Width - Panel3D1.Width - 100

vsViewPort1.Height = (frmPrint.Height - 400)

vsViewPort1.Width = (Panel3D1.Left - 100)

End If

SetPreviewSize

End Sub

Private Sub HScroll1_Change()

Static bProcedureActive

If bProcedureActive Then

HScroll1.Refresh

Exit Sub

End If

bProcedureActive = True *'prevent recursive calls to this procedure*

HScroll1.Enabled = False

frmPrint.vsPrinter1.PreviewPage = HScroll1.Value

UpdatePageButtons

bProcedureActive = False *'prevent recursive calls to this procedure*

End Sub

Private Sub HScroll1_Scroll()

Static bProcedureActive

If bProcedureActive Then Exit Sub

bProcedureActive = True *'prevent recursive calls to this procedure*

lblPageNumber = HScroll1.Value

UpdatePageButtons

bProcedureActive = False *'prevent recursive calls to this procedure*

End Sub

frmPrint.frm - lblActivePrinter_Click

126

Private Sub lblActivePrinter_Click()

```

Static bProcedureActive
If bProcedureActive Then Exit Sub
bProcedureActive = True 'prevent recursive calls to this procedure

On Error GoTo btnChangePrinter_Click_Error
CommonDialog1.Min = 1 'set lowest page number to print
CommonDialog1.Max = giTotalPrintPages 'set highest page number to print
CommonDialog1.FromPage = 1 'set lowest page number to print
CommonDialog1.ToPage = giTotalPrintPages 'set highest page number to print

Set flags
'PD_HIDEPRINTTOFILE &H100000 The Print to File check box is not displayed
'PD_NOPAGENUMS &H30 Disables the Pages option button and the associated edit control
'PD_PRINTSETUP &H40 Causes the system to display the Print Setup dialog box rather than the Print dialog box
CommonDialog1.Flags = &H408
CommonDialog1.CancelError = True
CommonDialog1.Action = 5 'call printer common dialog
'Update caption to most current printer selection
lblActivePrinter.Caption = "" + vsPrinter1.Device
frmPrint.MousePointer = vbHourglass
DoEvents
SetPreviewSize 'this is mainly for layout if portrait/landscape is changed
RefreshPreview
frmPrint.MousePointer = vbDefault
DoEvents

btnChangePrinter_Click_Exit:
bProcedureActive = False 'prevent recursive calls to this procedure
Exit Sub

btnChangePrinter_Click_Error:
Resume btnChangePrinter_Click_Exit

End Sub

```

Private Sub optZoom_Click(Index As Integer)

```

SetPreviewSize
End Sub

```

Private Sub vsPrinter1_EndPage()

```

Call PrintPageNumber
End Sub

```

Private Sub vsPrinter1_Error()

```

gbPrinterErrorReceived = True 'tells other procs that error occurred. Proc must reset flag
If vsPrinter1.Error = 5 Then 'a cancel was received from the print options dialog
    vsPrinter1.Action = paStartDoc 'start doc
    vsPrinter1.Action = paEndDoc 'end doc
ElseIf vsPrinter1.Error = 3 Or vsPrinter1.Error = 4 Then 'can't access printer, or can't start job
    'an error code of 3 is generated when user presses the 'CANCEL' button from options dialog
    vsPrinter1.Action = paStartDoc 'start doc
    vsPrinter1.Action = paEndDoc 'end doc
    If gbPrinterErrorDetected = False Then 'warning has not yet been issued
        Beep
        MsgBox "The printer is not available. Please ensure it is powered on and is on-line.", , "Can't Print"
        gbPrinterErrorDetected = True
    End If
ElseIf vsPrinter1.Error = 6 Then 'already printing

```

141

frmPrint.frm - vsPrinter1_Error

127

```

If gbPrinterErrorDetected = False Then 'warning has not yet been issued
    Beep
    MsgBox "The printer is not available. Please ensure it is powered on and is on-line.", , "Can't Print"
    gbPrinterErrorDetected = True
End If
End Sub

```

```

Private Sub vsPrinter1_NewPage()

```

```

    Dim fCurrentFontSize As Single, bCurrentFontItalic As Boolean, sCurrentFontName As String
    Dim iCurrentTextAlign As Integer, iCurrentY As Long

```

```

    With frmPrint.vsPrinter1

```

```

        fCurrentFontSize = .FontSize 'remember the existing settings, so they can be changed back
        bCurrentFontItalic = .FontItalic
        iCurrentTextAlign = .TextAlign
        iCurrentY = .CurrentY
        sCurrentFontName = .FontName

```

```

        .FontName = "Arial"
        .FontItalic = False
        .TextAlign = taRightTop
        .CurrentY = 1440 * 0.5 'print name on of program
        .FontSize = 9 'set font size
        .FontItalic = False
        frmPrint.vsPrinter1 = App.Title

```

```

        .FontSize = fCurrentFontSize
        .FontItalic = bCurrentFontItalic
        .TextAlign = iCurrentTextAlign
        .CurrentY = iCurrentY
        .FontName = sCurrentFontName
    End With

```

```

    Select Case gsActiveFormName
        Case "tmPatientDosingReport"
            If giTotalPrintPages Then PrintDosingEventsHeader
    End Select

```

```

    giTotalPrintPages = giTotalPrintPages + 1
End Sub

```

frmDeviceDiagnostics.frm - File Declaration:

128

```

Attribute VB_Name = "frmDeviceDiagnostics"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit

```

```

Private Sub btnChangeBatteries_Click()
    Call ChangeBatteriesRequest
End Sub

```

```

Private Sub btnClose_Click()
    Unload Me
End Sub

```

```

Private Sub btnReadEntireContents_Click()
    Dim r As Integer, IErrorCode As Long, i As Integer

    r = ValidatePatientDataSaved 'ensure that previous patient data was saved before proceeding
    If r = vbCancel Then Exit Sub

    btnReadEntireContents.Enabled = False 'prevent recursive calls to device
    btnSendData.Enabled = False 'prevent recursive calls to device
    gbKeepPollingDevice = False 'stop polling for now
    Wait 0.25

    txtPatientLastName = "" 'clear out the text boxes before reading data
    txtPatientFirstName = "" 'clear out the text boxes before reading data
    txtDrug.Clear
    txtPatientID = ""
    txtTxCenter = ""
    txtOrgan.Clear
    txtSerialNumber = ""
    txtDoseSize = ""
    txtDoseTime(1) = ""
    txtDoseTime(2) = ""
    txtDoseTime(3) = ""
    txtDoseTime(4) = ""
    txtDosesPerDay = ""
    txtDoseLockoutHours = ""
    txtDeviceStarted = ""
    txtMedicationRemaining = ""
    txtBatteryChangeTimer = ""
    txtEventCount = ""
    txtFirmwareVer = ""

    r = Comm_ReadEntireMemoryContents(PAT_DATA, IErrorCode)
    If r Then
        'PopulateDeviceDiagDialog PAT_DATA, Me
        PAT_DATA.sPatientDataFileName = ""
        frmMain.mnuFileSave.Enabled = False
    Else
        EraseDataInMemory PAT_DATA
    End If

    gbKeepPollingDevice = True 'start polling again
    btnReadEntireContents.Enabled = True 're-enable button
    btnSendData.Enabled = True
    RefreshAllOpenForms
End Sub

```

frmDeviceDiagnostics.frm - btnSendData_Click

129

Private Sub btnSendData_Click()

Dim i As Integer, r As Integer, IErrorCode As Long

r = ValidateDoseNumbers(Me)

If r = False Then Exit Sub

Beep

r = MsgBox("Patient Information and Dosing Information currently in the CycloTech device will be changed if you continue. Medication data will be preserved." + vbCrLf + vbCrLf + "Do you want to continue?", vbYesNo + vbQuestion, "Device Data being changed")

If r = vbNo Then Exit Sub

btnSendData.Enabled = False

'prevent recursive calls to device

btnReadEntireContents.Enabled = False

gbKeepPollingDevice = False

'stop polling for now

Wait 0.25

On Error GoTo btnSendData_Click_Error

r = Comm_SendCustomData(PAT_DATA, DATA_BEGIN_CUSTOM1, IErrorCode)

If IErrorCode Then Error IErrorCode

'error number

r = Comm_SendCustomData(PAT_DATA, DATA_BEGIN_CUSTOM2, IErrorCode)

'send to device

If IErrorCode Then Error IErrorCode

'error number

r = Comm_SendCustomData(PAT_DATA, DATA_BEGIN_CUSTOM3, IErrorCode)

'send to device

If IErrorCode Then Error IErrorCode

'error number

r = Comm_SendCustomData(PAT_DATA, DATA_BEGIN_CUSTOM4, IErrorCode)

'send to device

If IErrorCode Then Error IErrorCode

'error number

'ensure that the values in the text boxes are converted into the global structure

For i = 1 To 4

If IsDate(txtDoseTime(i)) Then

PAT_DATA.dPrescribedDoseTime(i) = TimeValue(txtDoseTime(i))

'save Dose Interval

Else

PAT_DATA.dPrescribedDoseTime(i) = -1

'Indicate that no time was set

End If

Next i

r = Comm_SendDosingParams(PAT_DATA, IErrorCode)

If IErrorCode Then Error IErrorCode

'error number

btnSendData_Click_Exit:

btnSendData.Enabled = True

're-enable button

btnReadEntireContents.Enabled = True

gbKeepPollingDevice = True

'continue polling device

Exit Sub

btnSendData_Click_Error:

DisplayErrorMessage IErrorCode

Resume 0 'temp test

Resume btnSendData_Click_Exit

End Sub

frmDeviceDiagnostics.frm - Form_Activate

130

Private Sub Form_Activate()

```

PopulateDeviceDialog PAT_DATA, Me
Comm_InitializeCommPort 'initialize the comm port from INI file settings

lblCommPort.Caption = "" + CStr(giCommPort)
lblSettings.Caption = "" + gsCommDeviceSettings
lblDeviceWaitTime = "" + CStr(giDeviceResponseWait)

gbCommSusy = False 'reset flag
gbCommReplyPending = False 'reset flag
gbKeepPollingDevice = True 'continue polling device
PollDeviceContinually Me
SetPrinterIcon False, ""
End Sub

```

Private Sub Form_Initialize()

```

Me.Left = 0
Me.Top = 0
End Sub

```

Private Sub Form_Load()

```

Unload frmDeviceInitialize
Unload frmReadDeviceData
gbCommOK = 99 'reset flag that will give an indication as to the communication status.
End Sub

```

Private Sub Form_Unload(Cancel As Integer)

```

Dim r As Integer

r = ValidateDoseNumbers(Me)
If r = False Then Cancel = True

gbKeepPollingDevice = False 'stop polling the device
Wait 0.1
End Sub

```

Private Sub txtDoseTime_Change(Index As Integer)

```

If IsDate(txtDoseTime(Index)) Then
    PAT_DATA.dPrescribedDoseTime(Index) = TimeValue(txtDoseTime(Index)) 'save Dose Interval
Else
    PAT_DATA.dPrescribedDoseTime(Index) = -1 'indicate that no time was set
End If
End Sub

```

DeviceDiagnostics.frm - txtDoseLockoutHours_Change

131

```
Private Sub txtDoseLockoutHours_Change()  
    PAT_DATA.sDoseLockoutHours = txtDoseLockoutHours    'save Dose Lockout Hours  
End Sub
```

```
Private Sub txtDoseSize_Change()  
    PAT_DATA.sDoseSize = txtDoseSize    'save Dose Size  
End Sub
```

```
Private Sub txtDosesPerDay_Change()  
    PAT_DATA.sDosesPerDay = Val(txtDosesPerDay)    'save Doses per day  
End Sub
```

```
Private Sub txtDrug_Click()  
    PAT_DATA.sDrug = txtDrug    'save field  
End Sub
```

```
Private Sub txtOrgan_Click()  
    PAT_DATA.sOrgan = txtOrgan    'save field  
End Sub
```

```
Private Sub txtPatientFirstName_Change()  
    PAT_DATA.sPatientFirstName = txtPatientFirstName    'save field  
End Sub
```

```
Private Sub txtPatientID_Change()  
    PAT_DATA.sPatientID = txtPatientID    'save Patient ID  
End Sub
```

```
Private Sub txtPatientLastName_Change()  
    PAT_DATA.sPatientLastName = txtPatientLastName    'save field  
End Sub
```

```
Private Sub txtSerialNumber_Change()  
    PAT_DATA.sSerialNumber = txtSerialNumber    'save serial number  
End Sub
```

frmDeviceDiagnostics.frm - txtTxCenter_Change

132

```

Private Sub txtTxCenter_Change()
    PAT_DATA.sTxCenter = txtTxCenter    'save field
End Sub

```

```

Private Sub UpDownDoseTime_DownClick(Index As Integer)
    Dim fDalyIncrement As Single, iIndex As Integer

```

```

    If IsDate(txtDoseTime(Index)) Then
        iIndex = TimeValue(txtDoseTime(Index)) * 24
        iIndex = iIndex - 1
        If iIndex < 0 Then
            txtDoseTime(Index) = ""
            Exit Sub
        End If
    Else
        iIndex = 23
        End If

    fDalyIncrement = (iIndex / 24)
    txtDoseTime(Index) = "" & Format$(TimeValue(CDate(fDalyIncrement)), gsTimeDisplayFormat)
End Sub

```

```

Private Sub UpDownDoseTime_UpClick(Index As Integer)
    Dim fDalyIncrement As Single, iIndex As Integer

```

```

    If IsDate(txtDoseTime(Index)) Then
        iIndex = TimeValue(txtDoseTime(Index)) * 24
        iIndex = iIndex + 1
        If iIndex > 23 Then
            txtDoseTime(Index) = ""
            Exit Sub
        End If
    Else
        iIndex = 0
        End If

    fDalyIncrement = (iIndex / 24)
    txtDoseTime(Index) = "" & Format$(TimeValue(CDate(fDalyIncrement)), gsTimeDisplayFormat)
End Sub

```

frmFaxStatus.frm - File Declarations

133

```
Attribute VB_Name = "frmFaxStatus"  
Attribute VB_GlobalNameSpace = False  
Attribute VB_Creatable = False  
Attribute VB_PredeclaredId = True  
Attribute VB_Exposed = False  
Option Explicit
```

```
Private Sub cmdCancel_Click()  
    On Error Resume Next  
    gcFax.CancelFax gcFax.FaxLogID  
    Unload Me  
End Sub
```

```
Private Sub Form_Activate()  
    SetPrinterIcon False, ""  
End Sub
```

```
Private Sub Form_Load()  
    lblDestination = ""  
    lblPage = ""  
    lblSpeed = ""  
End Sub
```


frmFaxSend.frm - File Declarations

134

```

Attribute VB_Name = "frmFaxSend"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit

```

Sub ReloadGroupsList()

```
Dim i As Integer
```

```
With cmbGroups
```

```
.Clear
```

```
.AddItem "Select Recipients Manually"
```

```
For i = 1 To FAX_DATA.iGroupsTotal
```

```
If FAX_DATA.sGroupTitle(i) <> "" Then
```

```
.AddItem FAX_DATA.sGroupTitle(i)
```

```
.ItemData(.NewIndex) = i
```

```
End If
```

```
Next i
```

```
End With
```

```
End Sub
```

Sub ReloadLocationsList()

```
Dim i As Integer
```

```
With lstLocations
```

```
.Clear
```

```
For i = 1 To FAX_DATA.iLocTotal
```

```
.AddItem FAX_DATA.sLocPersonName(i)
```

```
.ItemData(.NewIndex) = i 'keep index
```

```
Next i
```

```
End With
```

```
End Sub
```

Private Sub btnAddGroup_Click()

```
Dim i As Integer, r As Integer, sSection As String
```

```
gsEditGroupName = ""
```

```
gsEditGroupIndexes = ""
```

```
Load frmFaxEditGroups
```

```
frmFaxEditGroups.Show vbModal
```

```
If gsEditGroupName = "" Then Exit Sub 'no name was entered
```

```
i = GetIndexToFaxGroupName(gsEditGroupName) 'see if name is already in the list
```

```
If i = 0 Then 'name is not in list yet
```

```
sSection = "Fax Groups"
```

```
With FAX_DATA
```

```
'Add the new name to the fax data structure
```

```
'Before saving new data, clear out the old strings
```

```
r = WritePrivateProfileString(sSection, ByVal 0&, ByVal 0&, gsFaxFileSpec)
```

```
.iGroupsTotal = .iGroupsTotal + 1 'increment count by one
```

```
.sGroupTitle(.iGroupsTotal) = gsEditGroupName
```

```
.sGroupNamesInTitle(.iGroupsTotal) = gsEditGroupIndexes
```

```
SaveINISetting gsFaxFileSpec, sSection, "Total Groups", CStr(.iGroupsTotal)
```

```
For i = 0 To .iGroupsTotal
```

```
SaveINISetting gsFaxFileSpec, sSection, "Group " + CStr(i), .sGroupTitle(i)
```

```
SaveINISetting gsFaxFileSpec, sSection, "Group Locations " + CStr(i), .sGroupNamesInTitle(i)
```

```
Next i
```

```
cmbGroups.AddItem gsEditGroupName
```

```
cmbGroups.ItemData(cmbGroups.NewIndex) = .iGroupsTotal 'save index
```

```
cmbGroups.ListIndex = cmbGroups.NewIndex
```

frmFaxSend.frm - btnAddGroup_Click

135

```

End With
Else this name already exist
MsgBox "The name " + gsEditGroupName + " is already entered."
End If

```

End Sub

Private Sub btnClose_Click()

```

Unload Me
End Sub

```

Private Sub btnDeleteGroup_Click()

Dim i As Integer, r As Integer, sMSG As String

gsEditGroupName = cmboGroups.List(cmboGroups.ListIndex)

sMSG = "The following name and related information will be permanently deleted." + vbCrLf + gsEditGroupName + vbCrLf + vbCrLf + "Do you want to delete it?"

r = MsgBox(sMSG, vbYesNo + vbDefaultButton2 + vbQuestion, "Confirm Name Deletion")

If r = vbNo Then Exit Sub

RemoveGroupFromFaxList gsEditGroupName

ReloadGroupsList

cmboGroups.ListIndex = 0 *default to manual selections*

End Sub

Private Sub btnDeleteName_Click()

Dim i As Integer, r As Integer, sMSG As String

gsEditName = lstLocations.List(lstLocations.ListIndex)

sMSG = "The following name and related information will be permanently deleted." + vbCrLf + gsEditName + vbCrLf + vbCrLf + "Do you want to delete it?"

r = MsgBox(sMSG, vbYesNo + vbDefaultButton2 + vbQuestion, "Confirm Name Deletion")

If r = vbNo Then Exit Sub

RemoveNameFromFaxList gsEditName

ReloadLocationsList

If lstLocations.ListCount < 1 Then

btnDeleteName.Enabled = False

End If

cmboGroups_Click *cause appropriate boxes to be reselected*

End Sub

frmFaxSend.frm - btnEditGroup_Click

136

Private Sub btnEditGroup_Click()

Dim i As Integer, r As Integer, sSection As String

With cmboGroups

gsEditGroupName = .List(.ListIndex)

i = GetIndexToFaxGroupName(gsEditGroupName) *'get index from structure*

gsEditGroupIndexes = FAX_DATA.sGroupNamesInTitle(i)

frmFaxEditGroups.Show vbModal

.List(.ListIndex) = gsEditGroupName

End With

With FAX_DATA

.sGroupTitle(i) = gsEditGroupName

.sGroupNamesInTitle(i) = gsEditGroupIndexes

r = WritePrivateProfileString(sSection, ByVal 0&, ByVal 0&, gsFaxFileSpec)

SaveINISetting gsFaxFileSpec, sSection, "Total Groups", CStr(.iGroupsTotal)

For i = 0 To .iGroupsTotal

SaveINISetting gsFaxFileSpec, sSection, "Group " + CStr(i), .sGroupTitle(i)

SaveINISetting gsFaxFileSpec, sSection, "Group Locations " + CStr(i), .sGroupNamesInTitle(i)

Next i

End With

UpdateListBoxSelections gsEditGroupIndexes

End Sub

Private Sub btnEditName_Click()

Dim i As Integer

gsEditName = lstLocations.List(lstLocations.ListIndex)

i = GetIndexToFaxLocName(gsEditName) *'get index from structure*

gsEditVoice = FAX_DATA.sLocVoiceNumber(i)

gsEditFax = FAX_DATA.sLocFaxNumber(i)

frmFaxEditLocations.Show vbModal

lstLocations.List(lstLocations.ListIndex) = gsEditName

FAX_DATA.sLocPersonName(i) = gsEditName

FAX_DATA.sLocVoiceNumber(i) = gsEditVoice

FAX_DATA.sLocFaxNumber(i) = gsEditFax

End Sub

Private Sub btnNew_Click()

End Sub

frmFaxSend.frm - btnNewName_G

12

Private Sub btnNewName_Click()

Dim i As Integer

gsEditName = ""

gsEditVoice = ""

gsEditFax = ""

Load frmFaxEditLocations

frmFaxEditLocations.Caption = "Enter New Name"

frmFaxEditLocations.Show vbModal

If gsEditName = "" Then Exit Sub 'no name was entered

i = GetIndexToFaxLocName(gsEditName)

If i = 0 Then 'name is not in list yet 'see if name is already in the list

With FAX_DATA

'Add the new name to the fax data structure

.iLocTotal = .iLocTotal + 1 'increment count by one

.sLocPersonName(.iLocTotal) = gsEditName

.sLocVoiceNumber(.iLocTotal) = gsEditVoice

.sLocFaxNumber(.iLocTotal) = gsEditFax

lstLocations.AddItem gsEditName

lstLocations.ItemData(lstLocations.NewIndex) = .iLocTotal 'save index

End With

Else 'this name already exist

MsgBox "The name " + gsEditName + " is already entered."

End If

End Sub

Private Sub btnSendFax_Click()

Dim i As Integer, r As Integer, sFileSpec As String, iErrorCode As Long

Dim sSourceFileSpec As String, sDestFileSpec As String

On Error GoTo btnSendFax_Error

If Len(txtFileToSend) < 2 Then

MsgBox "The is no information to fax. Please open a patient file.", vbExclamation, "No File Selected"

Exit Sub

End If

CreateTxtSummaryFile

'Copy the report information to a text file for conversion to a fax document.

sSourceFileSpec = App.Path + "Vaxes\" + PAT_DATA.sPatientLastName + " " + PAT_DATA.sPatientFirstName + " " + PAT_DATA.sPatientID + ".txt"

sDestFileSpec = App.Path + "Vaxes\" + PAT_DATA.sPatientID + " " + PAT_DATA.sPatientID + ".frm"

r = FileExists(sDestFileSpec, iErrorCode)

If r Then

On Error Resume Next 'this won't be needed if we can determine that the file is already open

Kill sFileSpec

On Error GoTo btnSendFax_Error

End If

If lstLocations.ListCount Then

With gcFax

frmMain.FaxMan1.ImportFiles "c:\VaxesTemp.frm", "c:\sample if c:\cover.tif"

ImportFiles App.Path + "VaxesTemp.frm", App.Path + "VaxesTemp.txt"

.ImportFiles sDestFileSpec, sSourceFileSpec

.FaxFiles = sDestFileSpec

.FaxResolution = FAX_DATA.bFaxResolution

.UserCompany = FAX_DATA.sSenderCompany

frmFaxSend.frm - btnSendFax_On

```

13
.UserName = FAX_DATA.sSenderName
.UserVoiceNumber = FAX_DATA.sSenderVoiceNumber
.UserFaxNumber = FAX_DATA.sSenderFaxNumber
.FaxID = FAX_DATA.sFaxID
.FaxRetries = FAX_DATA.iRetries
.FaxRetryInterval = FAX_DATA.iRetryInterval
End With
End If

'Loop through all check boxes in the list to see which ones to send faxes to
For i = 0 To lstLocations.ListCount - 1
    With gcFax
        If lstLocations.Selected(i) Then 'this location is selected
            .FaxSubject = App.Title + " Report"
            .FaxName = FAX_DATA.sLocPersonName(lstLocations.ItemData(i) + 1)
            .FaxCompany = ""
            .FaxNumber = FAX_DATA.sDialPrefix + FAX_DATA.sLocFaxNumber(lstLocations.ItemData(i) + 1)
            .SendFax
        End If
    End With
Next i

btnSendFax_Exit:
'Unload Me 'must be unloaded because the status form is a nonmodal child form
'and can not be displayed while a modal form is being displayed.
Exit Sub

btnSendFax_Error:
MsgBox "An uncorrectable error occurred while trying to fax the document. Please try again.", vbExclamation, "Fax Error - " & Err$
Resume btnSendFax_Exit

End Sub

Private Sub cmbGroups_Click()
    Dim i As Integer, r As Integer, j As Integer, sName As String

    With cmbGroups
        i = GetIndexToFaxGroupName(.List(.ListIndex))
        UpdateListBoxSelections FAX_DATA.sGroupNamesInTitle(i)

        If .ListIndex > 0 Then 'the manual selection was made
            btnEditGroup.Enabled = True
            btnDeleteGroup.Enabled = True
        Else
            btnEditGroup.Enabled = False
            btnDeleteGroup.Enabled = False
        End If
    End With

    End With

End Sub

```

frmFaxSend.frm - UpDateListBoxSel .jns

11

```
Private Sub UpDateListBoxSelections(ByVal sGroup As String)
    Dim i As Integer, j As Integer, r As Integer, sTempList(100) As String
```

```
    With lstLocations
```

```
        .Clear
```

```
        For i = 0 To FAX_DATA.iLocTotal - 1
```

```
        .AddItem FAX_DATA.sLocPersonName(i + 1)
```

```
        .ItemData(.NewIndex) = i 'keep index
```

```
        Next i
```

```
    'Parse the attached locations and check the appropriate boxes
```

```
    r = ParseDelimString(sGroup, "T", sTempList())
```

```
    If r Then 'some locations are attached
```

```
        For j = 1 To r
```

```
            For i = 0 To .ListCount - 1 'step through the names and check appropriate ones
```

```
            If .ItemData(i) = sTempList(j) Then
```

```
                .Selected(i) = True
```

```
            Exit For
```

```
        End If
```

```
        Next i
```

```
    Next j
```

```
    End If
```

```
End With
```

```
End Sub
```

```
Private Sub Form_Activate()
```

```
    SetPrinterIcon False, "
```

```
End Sub
```

```
Private Sub Form_Load()
```

```
    brtFileToSend = PAT_DATA.sPatientDataFileName
```

```
    GetFaxLocations
```

```
    ReloadLocationsList
```

```
    ReloadGroupsList
```

```
    'make sure it is in range
```

```
    If cmbGroups.ListCount >= FAX_DATA.iGroupLastSelected Then cmbGroups.ListIndex = FAX_DATA.iGroupLastSelected
```

```
End Sub
```

```
Private Sub Form_Unload(Cancel As Integer)
```

```
    Dim i As Integer, r As Integer, sSection As String
```

```
    With FAX_DATA
```

```
        sSection = "User Selections"
```

```
        iGroupLastSelected = cmbGroups.ListIndex
```

```
        SaveINISetting gsFaxFileSpec, sSection, "Last Group Selected", CStr(iGroupLastSelected)
```

```
    End With
```

```
    With FAX_DATA
```

```
        sSection = "Fax Locations"
```

```
        'Before saving new data, clear out the old strings
```

```
        r = WritePrivateProfileString(sSection, ByVal 0&, ByVal 0&, gsFaxFileSpec)
```

```
        SaveINISetting gsFaxFileSpec, sSection, "Total Locations", CStr(iLocTotal)
```

```
        For i = 1 To iLocTotal
```

```
            SaveINISetting gsFaxFileSpec, sSection, "Person " + CStr(i), sLocPersonName(i)
```

```
            SaveINISetting gsFaxFileSpec, sSection, "Fax " + CStr(i), sLocFaxNumber(i)
```

```
            SaveINISetting gsFaxFileSpec, sSection, "Voice " + CStr(i), sLocVoiceNumber(i)
```

```
        Next i
```

```
        sSection = "Fax Groups"
```

frmFaxSend.frm - Form_Unload

14

```

'Before saving new data, clear out the old strings
r = WritePrivateProfileString(sSection, ByVal 0&, ByVal 0&, gsFaxFileSpec)
SaveINISetting gsFaxFileSpec, sSection, "Total Groups", CStr(.iGroupsTotal)
For i = 0 To .iGroupsTotal
    SaveINISetting gsFaxFileSpec, sSection, "Group " + CStr(i), .sGroupTitle(i)
    SaveINISetting gsFaxFileSpec, sSection, "Group Locations " + CStr(i), .sGroupNamesInTitle(i)
Next i
End With
End Sub

```

```

Private Sub lstLocations_Click()
    Dim i As Integer, sTemp As String

```

```

    btnEditName.Enabled = True
    btnDeleteName.Enabled = True

```

```

    If cmbGroups.ListIndex = 0 Then 'manual selection of groups is enabled
        With lstLocations

```

```

            For i = 0 To .ListCount - 1

```

```

                If .Selected(i) Then sTemp = sTemp + CStr(.ItemData(i)) + " "

```

```

            Next i

```

```

        End With

```

```

        FAX_DATA.sGroupNamesInTitle(0) = sTemp

```

```

    End If 'position 0 holds manual selections
End Sub

```

```

Private Sub lstLocations_DblClick()
    btnEditName_Click
End Sub

```

frmFaxLog.frm - File Declaration

```

Attribute VB_Name = "frmFaxLog"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit

```

```

Private Sub btnClose_Click()
    Me.Hide 'stay loaded because it contains the fax control
End Sub

```

```

Private Sub Form_Activate()
    SetPrinterIcon False, -
End Sub

```

```

Private Sub Form_Load()
    optViewFaxes_Click 1 'cause the "Sent" button to be clicked
End Sub

```

```

Private Sub Form_Resize()
    btnClose.Left = Me.Width - btnClose.Width - 250
    FaxMan1.Width = Me.Width - FaxMan1.Left - 250
    FaxMan1.Height = Me.Height - FaxMan1.Top - 500
End Sub

```

```

Private Sub optViewFaxes_Click(Index As Integer)
    Select Case Index
        Case 0
            FaxMan1.Log = Pending
        Case 1
            FaxMan1.Log = Completed
        Case 2
            FaxMan1.Log = Failed
    End Select
End Sub

```


frmFaxEditGroups.frm - File Declarations

```
Attribute VB_Name = "frmFaxEditGroups"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit
```

```
Private Sub btnClose_Click()
    Unload Me
```

```
End Sub
```

```
Private Sub Form_Activate()
    SetPrinterIcon False, ""
```

```
End Sub
```

```
Private Sub Form_Load()
```

```
    Dim i As Integer, j As Integer, r As Integer, sTempList(100) As String
```

```
    txtName = gsEditGroupName
```

```
    With lstLocations
```

```
        .Clear
```

```
        For i = 0 To FAX_DATA.JLocTotal - 1
```

```
            .AddItem FAX_DATA.sLocPersonName(i + 1)
```

```
            .ItemData(.NewIndex) = i 'keep index
```

```
        Next i
```

```
    'Parse the attached locations and check the appropriate boxes
```

```
    r = ParseDelimString(gsEditGroupIndexes, T, sTempList())
```

```
    If r Then 'some locations are attached
```

```
        For j = 1 To r
```

```
            For i = 0 To .ListCount - 1 'step through the names and check appropriate ones
```

```
                If .ItemData(i) = sTempList(j) Then
```

```
                    .Selected(i) = True
```

```
                    Exit For
```

```
                End If
```

```
            Next i
```

```
        Next j
```

```
    End If
```

```
    End With
```

```
End Sub
```

```
Private Sub Form_Unload(Cancel As Integer)
```

```
    Dim i As Integer, r As Integer, sSection As String
```

```
    gsEditGroupIndexes = ""
```

```
    With lstLocations
```

```
        For i = 0 To .ListCount - 1
```

```
            'add this index to the list
```

```
            If .Selected(i) Then gsEditGroupIndexes = gsEditGroupIndexes + CStr(.ItemData(i)) + T
```

```
        Next i
```

```
    End With
```

```
    gsEditGroupName = txtName
```

```
End Sub
```

frmFaxEditGroups.frm - lstLocations_KeyPress

```
Private Sub lstLocations_KeyPress(KeyAscii As Integer)
    If KeyAscii = 13 Then 'the "Enter" key was pressed
        btnClose.SetFocus
        KeyAscii = 0 'change it to a tab key
    End If
End Sub
```

```
Private Sub txtName_KeyPress(KeyAscii As Integer)
    If KeyAscii = 13 Then 'the "Enter" key was pressed
        lstLocations.SetFocus
        KeyAscii = 0 'change it to a tab key
    End If
End Sub
```

frmFaxEditLocations.frm - File Declarations

```

Attribute VB_Name = "frmFaxEditLocations"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit
Dim bxSaveData As Boolean

```

```

Private Sub btnCancel_Click()
    Unload Me
End Sub

```

```

Private Sub btnClose_Click()
    bxSaveData = True
    Unload Me
End Sub

```

```

Private Sub Form_Activate()
    txtName.SetFocus
    SetPrinterIcon False, =
End Sub

```

```

Private Sub Form_Load()
    txtName = gsEditName
    txtVoiceNumber = gsEditVoice
    txtFaxNumber = gsEditFax
End Sub

```

```

Private Sub Form_Unload(Cancel As Integer)
    If bxSaveData Then
        gsEditName = Trim$(txtName)
        gsEditVoice = Trim$(txtVoiceNumber)
        gsEditFax = Trim$(txtFaxNumber)
    End If
End Sub

```

```

Private Sub txtFax_KeyPress(KeyAscii As Integer)
    If KeyAscii = 13 Then 'the "Enter" key was pressed
        btnClose.SetFocus
        KeyAscii = 0 'change it to a tab key
    End If
End Sub

```

frmFaxEditLocations.frm - txtFaxNumber_.yPress

Private Sub txtFaxNumber_KeyPress(KeyAscii As Integer)
If KeyAscii = 13 Then *the "Enter" key was pressed*
txtVoiceNumber.SetFocus
KeyAscii = 0 *change it to a tab key*
End If
End Sub

Private Sub txtName_KeyPress(KeyAscii As Integer)
If KeyAscii = 13 Then *the "Enter" key was pressed*
txtFaxNumber.SetFocus
KeyAscii = 0 *change it to a tab key*
End If
End Sub

Private Sub txtTelephone_KeyPress(KeyAscii As Integer)
If KeyAscii = 13 Then *the "Enter" key was pressed*
txtFaxNumber.SetFocus
KeyAscii = 0 *change it to a tab key*
End If
End Sub

Private Sub txtVoiceNumber_KeyPress(KeyAscii As Integer)
If KeyAscii = 13 Then *the "Enter" key was pressed*
btnClose.SetFocus
KeyAscii = 0 *change it to a tab key*
End If
End Sub

frmDeviceInitialize.frm - File Declarations

14

```
Attribute VB_Name = "frmDeviceInitialize"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit
```

```
Private Sub btnChangeBatteries_Click()
    btnChangeBatteries.Enabled = False 'prevent recursive calls to device
    Call ChangeBatteriesRequest
    btnChangeBatteries.Enabled = True 'enable button again
End Sub
```

```
Private Sub btnClose_Click()
    Unload Me
```

```
End Sub
```

```
Private Sub btnReadEntireContents_Click()
    Dim r As Integer, IErrorCode As Long, i As Integer
```

```
    r = ValidatePatientDataSaved 'ensure that previous patient data was saved before proceeding
    If r = vbCancel Then Exit Sub
```

```
    btnReadEntireContents.Enabled = False 'prevent recursive calls to device
    btnSendData.Enabled = False 'prevent recursive calls to device
    gbKeepPollingDevice = False 'stop polling for now
    Wait 0.25
```

```
    txtPatientLastName = "" 'clear out the text boxes before reading data
    txtPatientFirstName = "" 'clear out the text boxes before reading data
    txtPatientID = ""
    txtTxCenter = ""
    txtDrug.Clear
    txtOrgan.Clear
    txtSerialNumber = ""
    txtDoseSize = ""
    txtDoseTime(1) = ""
    txtDoseTime(2) = ""
    txtDoseTime(3) = ""
    txtDoseTime(4) = ""
    txtDosesPerDay = ""
    txtDoseLockoutHours = ""
    txtDeviceStarted = ""
    txtMedicationRemaining = ""
    txtBatteryChangeTimer = ""
    txtEventCount = ""
    txtPatientLastName.SetFocus 'take focus away from listbox
    Me.Refresh
```

```
r = Comm_ReadEntireMemoryContents(PAT_DATA, IErrorCode)
```

```
If r Then
```

```
    PopulateDeviceCommDialog PAT_DATA, Me
```

```
    PAT_DATA.sPatientDataFileName = ""
```

```
    frmMain.mnuFileSave.Enabled = False
```

```
Else
```

```
    EraseDataInMemory PAT_DATA
```

```
    gbPatientDataNotSaved = False
```

```
    DisplayErrorMessage IErrorCode
```

```
End If
```

```
gbKeepPollingDevice = True 'start polling again
btnReadEntireContents.Enabled = True 're-enable button
```

161

frmDeviceInitialize.frm - btnReadEntireCon, is_Click

```

btnSendData.Enabled = True
RefreshAllOpenForms

```

```

End Sub

```

```

Private Sub btnSendData_Click()

```

```

    Dim i As Integer, r As Integer, IErrorCode As Long

```

```

    r = ValidateDoseNumbers(Me)
    If r = False Then Exit Sub

```

```

    Beep

```

```

    r = MsgBox("Patient Information and Dosing Information currently in the CycloTech device will be changed if you continue. Medication data will be preserved." + vbCrLf + vbCrLf + "Do you want to continue?", vbYesNo + vbQuestion, "Device Data being changed")
    If r = vbNo Then Exit Sub

```

```

    btnSendData.Enabled = False           'prevent recursive calls to device
    btnReadEntireContents.Enabled = False
    gbKeepPollingDevice = False          'stop polling for new
    Wait 0.25

```

```

    On Error GoTo btnSendData_Click_Error

```

```

    r = Comm_SendCustomData(PAT_DATA, DATA_BEGIN_CUSTOM1, IErrorCode)
    If IErrorCode Then Error IErrorCode    'error number

```

```

    r = Comm_SendCustomData(PAT_DATA, DATA_BEGIN_CUSTOM2, IErrorCode)
    If IErrorCode Then Error IErrorCode    'error number           'send to device

```

```

    r = Comm_SendCustomData(PAT_DATA, DATA_BEGIN_CUSTOM3, IErrorCode)
    If IErrorCode Then Error IErrorCode    'error number           'send to device

```

```

    r = Comm_SendCustomData(PAT_DATA, DATA_BEGIN_CUSTOM4, IErrorCode)
    If IErrorCode Then Error IErrorCode    'error number           'send to device

```

```

    'ensure that the values in the text boxes are converted into the global structure
    For i = 1 To 4

```

```

        If IsDate(txtDoseTime(i)) Then
            PAT_DATA.dPrescribedDoseTime(i) = TimeValue(txtDoseTime(i))    'save Dose Interval
        Else

```

```

            PAT_DATA.dPrescribedDoseTime(i) = -1    'Indicate that no time was set
        End If
    Next i

```

```

    r = Comm_SendDosingParams(PAT_DATA, IErrorCode)
    If IErrorCode Then Error IErrorCode    'error number

```

```

btnSendData_Click_Exit:

```

```

    btnSendData.Enabled = True    're-enable button
    btnReadEntireContents.Enabled = True
    gbKeepPollingDevice = True    'continue polling device
    Exit Sub

```

```

btnSendData_Click_Error:

```

```

    DisplayErrorMessage IErrorCode
    Resume 0 'temp test
    Resume btnSendData_Click_Exit
End Sub

```

frmDeviceInitialize.frm - Form_Activ...

14

```

Private Sub Form_Activate()
    PopulateDeviceCommDialog PAT_DATA, Me
    Comm_InitializeCommPort 'initialize the comm port from INI file settings

    gbKeepPollingDevice = True 'continue polling device
    PollDeviceContinually Me
    SetPrinterIcon False, ""
End Sub

```

```

Private Sub Form_Load()
    Me.Left = 0
    Me.Top = 0
    Unload frmReadDeviceData 'don't need this form
    gbCommOK = 99 'reset flag that will give an indication as to the communication status.
End Sub

```

```

Private Sub Form_QueryUnload(Cancel As Integer, UnloadMode As Integer)
    Dim r As Integer

    r = ValidateDoseNumbers(Me)
    If r = False Then Cancel = True
End Sub

```

```

Private Sub Form_Unload(Cancel As Integer)
    gbKeepPollingDevice = False 'stop polling the device
    Wait 0.1
End Sub

```

```

Private Sub txtDoseLockoutHours_Change()
    PAT_DATA.sDoseLockoutHours = txtDoseLockoutHours 'save Dose Lockout Hours
End Sub

```

```

Private Sub txtDoseSize_Change()
    PAT_DATA.sDoseSize = txtDoseSize 'save Dose Size
End Sub

```

```

Private Sub txtDosesPerDay_Change()
    PAT_DATA.iDosesPerDay = Val(txtDosesPerDay) 'save Doses per day
End Sub

```

frmDeviceInitialize.frm - txtDrug_Click...

14

```
Private Sub txtDrug_Click()
    PAT_DATA.sDrug = txtDrug
End Sub
```

```
Private Sub txtOrgan_Click()
    PAT_DATA.sOrgan = txtOrgan
End Sub
```

```
Private Sub txtPatientFirstName_Change()
    PAT_DATA.sPatientFirstName = txtPatientFirstName
End Sub
```

'save Patient name'

```
Private Sub txtPatientID_Change()
    PAT_DATA.sPatientID = txtPatientID
End Sub
```

```
Private Sub txtPatientLastName_Change()
    PAT_DATA.sPatientLastName = txtPatientLastName
End Sub
```

'save Patient name'

```
Private Sub txtSerialNumber_Change()
    PAT_DATA.sSerialNumber = txtSerialNumber
End Sub
```

'save serial number'

```
Private Sub txtTxCenter_Change()
    PAT_DATA.sTxCenter = txtTxCenter
End Sub
```

```
Private Sub UpDownDoseTime_DownClick(Index As Integer)
    Dim fDalyIncrement As Single, lIndex As Integer
```

```
    If IsDate(txtDoseTime(Index)) Then
        lIndex = TimeValue(txtDoseTime(Index)) * 24
        lIndex = lIndex - 1
        If lIndex < 0 Then
            txtDoseTime(Index) = ""
            Exit Sub
        End If
    Else
        lIndex = 23
    End If
```

```
    fDalyIncrement = (lIndex / 24)
    txtDoseTime(Index) = "" + Format$(TimeValue(CDate(fDalyIncrement)), "gsTimeDisplayFormat")
End Sub
```

164

frmDeviceInitialize.frm - UpDownDoseTime_Click

15

Private Sub UpDownDoseTime_Click(Index As Integer)
Dim fDalyIncrement As Single, iIndex As Integer

If IsDate(txtDoseTime(iIndex)) Then
iIndex = TimeValue(txtDoseTime(iIndex)) * 24
iIndex = iIndex + 1

If iIndex > 23 Then
txtDoseTime(iIndex) = ""
Exit Sub

End If

Else
iIndex = 0
End If

fDalyIncrement = (iIndex / 24)
txtDoseTime(iIndex) = "" + Format\$(TimeValue(CDate(fDalyIncrement)), "gsTimeDisplayFormat")
End Sub

frmGetDateTime: frm - File Declaration

15

```

Attribute VB_Name = "frmGetDateTime"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = False
Option Explicit

```

Private Sub btnDateCancel_Click()

```

    gdTempDateTime = 0      'indicate to user that a cancel occurred
    giTempCya = 0
    giTempCreatinine = 0
    gsTempCustomInfo = ""
    Unload frmGetDateTime
End Sub

```

Private Sub cmdDateOK_Click()

'validate the text boxes before exit

```

If Val(txtCya.Text) = 0 Then
    MsgBox "Please enter a CYA level", vbExclamation, "Value Required"
    txtCya.SetFocus
Exit Sub
End If

```

```

If Val(txtCreatinine.Text) = 0 Then
    MsgBox "Please enter a Creatinine level", vbExclamation, "Value Required"
    txtCreatinine.SetFocus
Exit Sub
End If

```

'On Error Resume Next

```

gdTempDateTime = CDate(txtDateEntry.Value)      'get date from control
gdTempDateTime = gdTempDateTime + CDate(txtTimeEntry.Time)
giTempCya = txtCya.Text
giTempCreatinine = txtCreatinine.Text
gsTempCustomInfo = txtCustomInfo

```

```

Unload frmGetDateTime
'On Error GoTo 0
End Sub

```

Private Sub Form_Activate()

```

SetPrinterIcon False, ""
End Sub

```

Private Sub Form_Load()

```

Me.Width = prlGetDate.Width + 90
Me.Height = prlGetDate.Height + 90
txtTimeEntry.Time = CStr(Time)
End Sub

```

Function Index

B

brwWebBrowser_DownloadComplete (frmBrowser.frm) , 95
brwWebBrowser_NavigateComplete (frmBrowser.frm) , 95
btnAddGroup_Click (frmFaxSend.frm) , 134
btnCancel_Click (frmFaxEditLocations.frm) , 144
btnChangeBatteries_Click (frmDeviceDiagnostics.frm) , 128
btnChangeBatteries_Click (frmDeviceInitialize.frm) , 146
btnChangeCompliance_Click (frmDosingCalendar.frm) , 111
btnClose_Click (frmAllPatients.frm) , 100
btnClose_Click (frmDeviceDiagnostics.frm) , 128
btnClose_Click (frmDeviceInitialize.frm) , 146
btnClose_Click (frmDosingCalendar.frm) , 111
btnClose_Click (frmFaxEditGroups.frm) , 142
btnClose_Click (frmFaxEditLocations.frm) , 144
btnClose_Click (frmFaxLog.frm) , 141
btnClose_Click (frmFaxSend.frm) , 135
btnClose_Click (frmPatientDosingRpt.frm) , 118
btnClose_Click (frmPrint.frm) , 124
btnClose_Click (frmReadDeviceData.frm) , 122
btnClose_Click (frmRecentDosingGraph.frm) , 108
btnConfigureFax_Click (frmOptions.frm) , 88
btnDateCancel_Click (frmGetDateTime.frm) , 151
btnDeleteGroup_Click (frmFaxSend.frm) , 135
btnDeleteName_Click (frmFaxSend.frm) , 135
btnDeleteUserEvent_Click (frmPatientDosingRpt.frm) , 119
btnEditGroup_Click (frmFaxSend.frm) , 136
btnEditName_Click (frmFaxSend.frm) , 136
btnNew_Click (frmFaxSend.frm) , 136
btnNewName_Click (frmFaxSend.frm) , 137
btnNewUserEvent_Click (frmPatientDosingRpt.frm) , 119
btnPrinter_Preview_Click_Proc (Printing.bas) , 46
btnPrintNow_Click (frmPrint.frm) , 124
btnPrintPage_Click (frmPrint.frm) , 124
btnReadEntireContents_Click (frmDeviceDiagnostics.frm) , 128

Function Index

btnReadEntireContents_Click (frmDeviceInitialize.frm) , 146
btnReadEntireContents_Click (frmReadDeviceData.frm) , 122
btnRefresh_Click (frmPrint.frm) , 124
btnSendData_Click (frmDeviceDiagnostics.frm) , 129
btnSendData_Click (frmDeviceInitialize.frm) , 147
btnSendFax_Click (frmFaxSend.frm) , 137

C

CalcDayDoseScore_AllDoses (General.bas) , 4
CalcDayDoseScore_OnTime (General.bas) , 4
CalcDaysInMonth (Calendar.bas) , 61
CalcDosesSumTakenOnSpecificDay (General.bas) , 5
CalculateAllPatientsComplianceOnDisk (frmAllPatients.frm) , 102
CalculateSinglePatientCompliance (frmAllPatients.frm) , 100
Calendar_DayChange (frmDosingCalendar.frm) , 111
Calendar_MonthChange (frmDosingCalendar.frm) , 111
Calendar_MouseMove (frmDosingCalendar.frm) , 112
Calendar_YearChange (frmDosingCalendar.frm) , 112
cboAddress_Click (frmBrowser.frm) , 96
cboAddress_KeyPress (frmBrowser.frm) , 96
ChangeBatteriesRequest (Comm.bas) , 28
chkDoseChanged_Click (frmDosingCalendar.frm) , 112
chkDoseChanged_Click (frmPatientDosingRpt.frm) , 119
chkDoses_Click (frmPatientDosingRpt.frm) , 119
chkDosesMissed_Click (frmDosingCalendar.frm) , 112
chkDosesNotComplied_Click (frmDosingCalendar.frm) , 112
chkDosesTaken_Click (frmDosingCalendar.frm) , 112
chkLoadTipsAtStartup_Click (frmTip.frm) , 99
chkUserDefined_Click (frmPatientDosingRpt.frm) , 120
chkWeekNumbers_Click (frmDosingCalendar.frm) , 112
cmboAverageDays_Click (frmRecentDosingGraph.frm) , 108
cmboChartType_Click (frmRecentDosingGraph.frm) , 108
cmboDataToView_Click (frmAllPatients.frm) , 101
cmboDataToView_Click (frmRecentDosingGraph.frm) , 108
cmboDateSelection_Click (frmAllPatients.frm) , 101

Function Index

cmboDateSelection_Click (frmRecentDosingGraph.frm) , 108
cmboGroups_Click (frmFaxSend.frm) , 138
cmdApply_Click (frmOptions.frm) , 88
cmdCancel_Click (frmFaxStatus.frm) , 133
cmdCancel_Click (frmLogin.frm) , 86
cmdCancel_Click (frmOptions.frm) , 89
cmdDateOK_Click (frmGetDateTime.frm) , 151
cmdNextTip_Click (frmTip.frm) , 99
cmdOK_Click (frmAbout.frm) , 92
cmdOK_Click (frmLogin.frm) , 86
cmdOK_Click (frmOptions.frm) , 89
cmdOK_Click (frmTip.frm) , 99
cmdSysInfo_Click (frmAbout.frm) , 92
Comm_CheckComm (Comm.bas) , 28
Comm_GetDeviceReply (Comm.bas) , 29
Comm_InitializeCommPort (Comm.bas) , 39
Comm_ReadEntireMemoryContents (Comm.bas) , 34
Comm_ReadFirmwareVersion (Comm.bas) , 30
Comm_SendCustomData (Comm.bas) , 35
Comm_SendDataToDevice (Comm.bas) , 45
Comm_SendDosingParams (Comm.bas) , 37
Comm_SendResetClockAndBattery (Comm.bas) , 35
Comm_SendSerialNumber (Comm.bas) , 38
CommTimer_Timer (frmMain.frm) , 76
ComputeIniSectionChecksum (General.bas) , 2
ConvertHexStringToAscii (Comm.bas) , 38
CreateCalendarTimeScale (frmDosingCalendar.frm) , 113
CreateChecksum (Comm.bas) , 39
CreateTxtSummaryFile (General.bas) , 6

D

Device_OnComm (Comm.bas) , 40
DisplayCommError (Comm.bas) , 30
DisplayCommOk (Comm.bas) , 31
DisplayCurrentTip (frmTip.frm) , 99

Function Index

DisplayErrorMessage (Comm.bas) , 31
DoNextTip (frmTip.frm) , 98
DrawAllCompliedDosesTaken (Calendar.bas) , 62
DrawAllDoseSizeChanges (Calendar.bas) , 61
DrawAllDosesMissed (Calendar.bas) , 63
DrawAllNonCompliedDosesTaken (Calendar.bas) , 63
DrawHorizontalLine (Printing.bas) , 47
DrawSingleCompliedDoseTaken (Calendar.bas) , 73
DrawSingleDoseMissed (Calendar.bas) , 74
DrawSingleDoseSizeChange (Calendar.bas) , 64
DrawSingleNonCompliedDoseTaken (Calendar.bas) , 65

E

EraseDataInMemory (General.bas) , 5
EstablishDeviceComm (Comm.bas) , 39
EventDelete (General.bas) , 2
EventInsert (General.bas) , 3

F

FaxMan1_ConfigurationDone (frmMain.frm) , 76
FaxMan1_FaxStatus (frmMain.frm) , 76
FileExists (General.bas) , 7
FindClosestDateInArray (General.bas) , 23
FindFirstMatchingDateInArray (General.bas) , 23
FindPrescribedDoseSizeForSpecificDay (General.bas) , 3
Form_Activate (frmAllPatients.frm) , 102
Form_Activate (frmDeviceDiagnostics.frm) , 130
Form_Activate (frmDeviceInitialize.frm) , 148
Form_Activate (frmDosingCalendar.frm) , 113
Form_Activate (frmFaxEditGroups.frm) , 142
Form_Activate (frmFaxEditLocations.frm) , 144
Form_Activate (frmFaxLog.frm) , 141
Form_Activate (frmFaxSend.frm) , 139
Form_Activate (frmFaxStatus.frm) , 133
Form_Activate (frmGetDateTime.frm) , 151
Form_Activate (frmOptions.frm) , 89

Function Index

Form_Activate (frmPatientDosingRpt.frm) , 120
Form_Activate (frmReadDeviceData.frm) , 123
Form_Activate (frmRecentDosingGraph.frm) , 109
Form_Activate (frmTip.frm) , 99
Form_Initialize (frmDeviceDiagnostics.frm) , 130
Form_Load (frmAbout.frm) , 92
Form_Load (frmAllPatients.frm) , 103
Form_Load (frmBrowser.frm) , 95
Form_Load (frmDeviceDiagnostics.frm) , 130
Form_Load (frmDeviceInitialize.frm) , 148
Form_Load (frmDosingCalendar.frm) , 113
Form_Load (frmFaxEditGroups.frm) , 142
Form_Load (frmFaxEditLocations.frm) , 144
Form_Load (frmFaxLog.frm) , 141
Form_Load (frmFaxSend.frm) , 139
Form_Load (frmFaxStatus.frm) , 133
Form_Load (frmGetDateTime.frm) , 151
Form_Load (frmLogin.frm) , 86
Form_Load (frmOptions.frm) , 89
Form_Load (frmPatientDosingRpt.frm) , 120
Form_Load (frmPrint.frm) , 124
Form_Load (frmReadDeviceData.frm) , 123
Form_Load (frmRecentDosingGraph.frm) , 109
Form_Load (frmSplash.frm) , 85
Form_Load (frmTip.frm) , 99
Form_MouseMove (frmDosingCalendar.frm) , 115
Form_QueryUnload (frmDeviceInitialize.frm) , 148
Form_QueryUnload (frmPrint.frm) , 125
Form_Resize (frmAllPatients.frm) , 103
Form_Resize (frmBrowser.frm) , 96
Form_Resize (frmDosingCalendar.frm) , 115
Form_Resize (frmFaxLog.frm) , 141
Form_Resize (frmPatientDosingRpt.frm) , 120
Form_Resize (frmPrint.frm) , 125
Form_Resize (frmRecentDosingGraph.frm) , 109

Function Index

Form_Unload (frmDeviceDiagnostics.frm) , 130
Form_Unload (frmDeviceInitialize.frm) , 148
Form_Unload (frmDosingCalendar.frm) , 116
Form_Unload (frmFaxEditGroups.frm) , 142
Form_Unload (frmFaxEditLocations.frm) , 144
Form_Unload (frmFaxSend.frm) , 139
Form_Unload (frmReadDeviceData.frm) , 123
Form_Unload (frmRecentDosingGraph.frm) , 110
frameView_MouseMove (frmDosingCalendar.frm) , 116

G

GetDrugRefNumber (Comm.bas) , 32
GetFaxLocations (Fax.bas) , 57
GetFileNameFromSpec (General.bas) , 17
GetIndexToFaxGroupName (Fax.bas) , 58
GetIndexToFaxLocName (Fax.bas) , 58
GetNISetting (General.bas) , 7
GetKeyValue (frmAbout.frm) , 93
GetOrganRefNumber (Comm.bas) , 32
GetPatientDataFromDisk (General.bas) , 7
GetProgramPreferences (General.bas) , 9
grid_AfterEdit (frmPatientDosingRpt.frm) , 120
grid_Click (frmAllPatients.frm) , 104
grid_DblClick (frmAllPatients.frm) , 104
grid_KeyDown (frmPatientDosingRpt.frm) , 121
grid_RowColChange (frmPatientDosingRpt.frm) , 121

H

HScroll1_Change (frmPrint.frm) , 125
HScroll1_Scroll (frmPrint.frm) , 125

I

InitPageMargins (Printing.bas) , 52
InitPageProperties (Printing.bas) , 53
InterpretDosingData (Comm.bas) , 32
InterpretErrorFlags (Comm.bas) , 33

Function Index

InterpretScoreData (Comm.bas) , 41

IsDoseWithinPrescribedTimeRange (Calendar.bas) , 66

L

lblActivePrinter_Click (frmPrint.frm) , 126

LoadPictureToPrinterControl (Printing.bas) , 52

LoadTips (frmTip.frm) , 98

IstLocations_Click (frmFaxSend.frm) , 140

IstLocations_OblClick (frmFaxSend.frm) , 140

IstLocations_KeyPress (frmFaxEditGroups.frm) , 143

M

Main (General.bas) , 11

MDIForm_Load (frmMain.frm) , 77

MDIForm_Unload (frmMain.frm) , 77

mnuAccessWebSite_Click (frmMain.frm) , 77

mnuFaxConfigure_Click (frmMain.frm) , 77

mnuFaxSend_Click (frmMain.frm) , 77

mnuFaxViewLogs_Click (frmMain.frm) , 77

mnuFileExit_Click (frmMain.frm) , 84

mnuFileMRU_Click (frmMain.frm) , 84

mnuFileOpen_Click (frmMain.frm) , 83

mnuFilePageSetup_Click (frmMain.frm) , 84

mnuFilePrint_Click (frmMain.frm) , 84

mnuFileProperties_Click (frmMain.frm) , 78

mnuFileSave_Click (frmMain.frm) , 78

mnuFileSaveAs_Click (frmMain.frm) , 83

mnuFileSend_Click (frmMain.frm) , 84

mnuGenError_Click (frmMain.frm) , 78

mnuHelpAbout_Click (frmMain.frm) , 79

mnuHelpContents_Click (frmMain.frm) , 82

mnuHelpDeviceDiag_Click (frmMain.frm) , 78

mnuHelpSearch_Click (frmMain.frm) , 82

mnuHelpTips_Click (frmMain.frm) , 78

mnuReadDeviceData_Click (frmMain.frm) , 79

mnuSendDeviceData_Click (frmMain.frm) , 79

Function Index

mnuViewAllPatients_Click (frmMain.frm) , 79
mnuViewCalendar_Click (frmMain.frm) , 79
mnuViewExplorer_Click (frmMain.frm) , 79
mnuViewOptions_Click (frmMain.frm) , 79
mnuViewPatientDosingReport_Click (frmMain.frm) , 80
mnuViewPatientSummary_Click (frmMain.frm) , 80
mnuViewStatusBar_Click (frmMain.frm) , 80
mnuViewToolBar_Click (frmMain.frm) , 80
mnuWindowArrangeIcons_Click (frmMain.frm) , 83
mnuWindowCascade_Click (frmMain.frm) , 83
mnuWindowTileHorizontal_Click (frmMain.frm) , 83
mnuWindowTileVertical_Click (frmMain.frm) , 83
MoveFormObjects (Calendar.bas) , 72

O

OpenPatientData (General.bas) , 13
optViewFaxes_Click (frmFaxLog.frm) , 141
optZoom_Click (frmPrint.frm) , 126

P

ParseDelimString (General.bas) , 18
ParseMemoryContents (Comm.bas) , 41
pnlControls_MouseMove (frmDosingCalendar.frm) , 116
PollDeviceContinually (Comm.bas) , 44
PopulateDeviceCommDialog (General.bas) , 21
PopulateDeviceDiagDialog (General.bas) , 14
PrintAllPatientsSummary (Printing.bas) , 47
PrintCalendar (Calendar.bas) , 66
PrintDosingEventsHeader (Printing.bas) , 51
PrintPageDate (Printing.bas) , 53
PrintPageNumber (Printing.bas) , 53
PrintPatientDosingReport (Printing.bas) , 49

R

RefreshAllOpenForms (General.bas) , 16
RefreshPreview (Printing.bas) , 54

Function Index

ReloadGroupsList (frmFaxSend.frm) , 134
ReloadLocationsList (frmFaxSend.frm) , 134
RemoveAllObjects (Calendar.bas) , 75
RemoveCompliedDosesTaken (Calendar.bas) , 70
RemoveDoseSizeChanges (Calendar.bas) , 70
RemoveDosesMissed (Calendar.bas) , 70
RemoveGroupFromFaxList (Fax.bas) , 58
RemoveNameFromFaxList (Fax.bas) , 59
RemoveNonCompliedDosesTaken (Calendar.bas) , 70
RescaleGrid (frmPatientDosingRpt.frm) , 117

S

SaveDataToNewFile (General.bas) , 18
SaveINISetting (General.bas) , 24
SavePatientData (General.bas) , 19
SaveProgramPreferences (General.bas) , 22
SetCommTimer (Comm.bas) , 45
SetFaxDeviceLabel (Fax.bas) , 60
SetPreviewSize (Printing.bas) , 55
SetPrinterIcon (General.bas) , 16
Slider1_SlideChange (frmAllPatients.frm) , 104
SSLstBar1_ListItemClick (frmMain.frm) , 80
sstab1_Click (frmOptions.frm) , 91
StartSysInfo (frmAbout.frm) , 92

T

tbToolBar_ButtonClick (frmBrowser.frm) , 96
tbToolBar_ButtonClick (frmMain.frm) , 81
timTimer_Timer (frmBrowser.frm) , 96
txtDoseLockoutHours_Change (frmDeviceDiagnostics.frm) , 131
txtDoseLockoutHours_Change (frmDeviceInitialize.frm) , 148
txtDoseSize_Change (frmDeviceDiagnostics.frm) , 131
txtDoseSize_Change (frmDeviceInitialize.frm) , 148
txtDosesPerDay_Change (frmDeviceDiagnostics.frm) , 131
txtDosesPerDay_Change (frmDeviceInitialize.frm) , 148
txtDoseTime_Change (frmDeviceDiagnostics.frm) , 130

Function Index

txtDrug_Click (frmDeviceDiagnostics.frm) , 131
 txtDrug_Click (frmDeviceInitialize.frm) , 149
 txtEndDate_Change (frmRecentDosingGraph.frm) , 110
 txtEndDate_HideDropDown (frmAllPatients.frm) , 105
 txtFax_KeyPress (frmFaxEditLocations.frm) , 144
 txtFaxNumber_GotFocus (frmOptions.frm) , 91
 txtFaxNumber_KeyPress (frmFaxEditLocations.frm) , 145
 txtName_KeyPress (frmFaxEditGroups.frm) , 143
 txtName_KeyPress (frmFaxEditLocations.frm) , 145
 txtOrgan_Click (frmDeviceDiagnostics.frm) , 131
 txtOrgan_Click (frmDeviceInitialize.frm) , 149
 txtPatientFirstName_Change (frmDeviceDiagnostics.frm) , 131
 txtPatientFirstName_Change (frmDeviceInitialize.frm) , 149
 txtPatientID_Change (frmDeviceDiagnostics.frm) , 131
 txtPatientID_Change (frmDeviceInitialize.frm) , 149
 txtPatientLastName_Change (frmDeviceDiagnostics.frm) , 131
 txtPatientLastName_Change (frmDeviceInitialize.frm) , 149
 txtRetries_GotFocus (frmOptions.frm) , 91
 txtRetryInterval_GotFocus (frmOptions.frm) , 91
 txtSerialNumber_Change (frmDeviceDiagnostics.frm) , 131
 txtSerialNumber_Change (frmDeviceInitialize.frm) , 149
 txtStartDate_Change (frmRecentDosingGraph.frm) , 110
 txtStartDate_HideDropDown (frmAllPatients.frm) , 105
 txtTelephone_KeyPress (frmFaxEditLocations.frm) , 145
 txtTxCenter_Change (frmDeviceDiagnostics.frm) , 132
 txtTxCenter_Change (frmDeviceInitialize.frm) , 149
 txtVoiceNumber_GotFocus (frmOptions.frm) , 91
 txtVoiceNumber_KeyPress (frmFaxEditLocations.frm) , 145

U

UpdateCalendar (Calendar.bas) , 74
 UpdatefrmPatientDosingReportHeader (frmPatientDosingRpt.frm) , 117
 UpdatefrmPatientSummaryHeader (frmRecentDosingGraph.frm) , 106
 UpdateListBoxSelections (frmFaxSend.frm) , 139
 UpdatePageButtons (Printing.bas) , 56

Function Index

UpdatePatientDosingGraph (frmRecentDosingGraph.frm) , 106
UpdatePatientGridDisplay (frmPatientDosingRpt.frm) , 117
UpDateRecentFileMenu (General.bas) , 17
UpdateZoomBox (Calendar.bas) , 71
UpDownDoseTime_DownClick (frmDeviceDiagnostics.frm) , 132
UpDownDoseTime_DownClick (frmDeviceInitialize.frm) , 149
UpDownDoseTime_UpClick (frmDeviceDiagnostics.frm) , 132
UpDownDoseTime_UpClick (frmDeviceInitialize.frm) , 150

V

ValidateChecksum (Comm.bas) , 33
ValidateDoseNumbers (General.bas) , 24
ValidatePatientDataSaved (General.bas) , 24
vsPrinter1_EndPage (frmPrint.frm) , 126
vsPrinter1_Error (frmPrint.frm) , 126
vsPrinter1_NewPage (frmPrint.frm) , 127

W

Wait (General.bas) , 25

B

b LogonToWebSite (General.bas) , 8046736

F

File Declarations (Calendar.bas) , 61
File Declarations (Comm.bas) , 26
File Declarations (Fax.bas) , 57
File Declarations (frmAbout.frm) , 92
File Declarations (frmAllPatients.frm) , 100
File Declarations (frmBrowser.frm) , 95
File Declarations (frmDeviceDiagnostics.frm) , 128
File Declarations (frmDeviceInitialize.frm) , 146
File Declarations (frmDosingCalendar.frm) , 111
File Declarations (frmFaxEditGroups.frm) , 142
File Declarations (frmFaxEditLocations.frm) , 144
File Declarations (frmFaxLog.frm) , 141

Function Index

File Declarations (frmFaxSend.frm) , 134
File Declarations (frmFaxStatus.frm) , 133
File Declarations (frmGetDateTime.frm) , 151
File Declarations (frmLogin.frm) , 86
File Declarations (frmMain.frm) , 76
File Declarations (frmOptions.frm) , 88
File Declarations (frmPatientDosingRpt.frm) , 117
File Declarations (frmPrint.frm) , 124
File Declarations (frmReadDeviceData.frm) , 122
File Declarations (frmRecentDosingGraph.frm) , 106
File Declarations (frmSplash.frm) , 85
File Declarations (frmTip.frm) , 98
File Declarations (General.bas) , 1
File Declarations (Printing.bas) , 46

CLAIMS

What is claimed is:

1. A computer-implemented method for monitoring medication dosing by a patient, comprising:
storing patient data, including a medication name and amounts of the medication prescribed for a patient;
5 retrieving patient data, including times and amounts of the medication delivered to the patient;
evaluating data by analyzing drug dispensing data and patient data to determine compliance of the delivered medication to the prescribed medication; and displaying the evaluated data.
- 10 2. The method of claim 1, wherein storing patient data further comprises storing information from a remote device over a communications line.
3. The method of claim 1, wherein storing patient data further comprises storing data from local memory.
- 15 4. The method of claim 1, wherein storing patient data further comprises storing user input.
5. The method of claim 1, wherein retrieving patient data further comprises retrieving data from local memory.
6. The method of claim 1, wherein retrieving patient data further comprises retrieving user input.
- 20 7. The method of claim 1, wherein displaying the evaluated data displays the evaluated data in a patient summary report.
8. The method of claim 1, further comprising printing the evaluated data.
9. The method of claim 1, wherein dosages of multiple patients are monitored, the method comprising:

storing patient data for a plurality of patients, including the medication name and amounts of the medication prescribed for the plurality of patients;
retrieving patient data for the plurality of patients, including times and amounts of medication delivered to the plurality of patients;
5 evaluating data by analyzing the stored patient data for the plurality of patients to determine overall compliance of the delivered medication to the prescribed medication; and displaying the evaluated data.

10. The method of Claim 1, wherein:
storing patient data includes storing amounts of an immunosuppressive medication
10 prescribed for a patient; and
retrieving patient data includes retrieving times and amounts of the immunosuppressive medication delivered to the patient.

11. The method of Claim 1, wherein:
storing patient data includes storing amounts of an analgesic drug prescribed for a
15 patient; and
retrieving patient data includes retrieving times and amounts of the analgesic drug delivered to the patient.

12. The method of Claim 1, wherein:
storing patient data includes storing amounts of an opiate agonist prescribed for a
20 patient; and
retrieving patient data includes retrieving times and amounts of the opiate agonist delivered to the patient.

13. The method of Claim 1, wherein:
storing patient data includes storing amounts of an opiate antagonist prescribed for a
25 patient; and
retrieving patient data includes retrieving times and amounts of the opiate antagonist delivered to the patient.

14. The method of Claim 1, wherein:
storing patient data includes storing amounts of a liquid drug prescribed for a patient;
30 and
retrieving patient data includes retrieving times and amounts of the liquid drug delivered to the patient.

15. The method of Claim 1, wherein the step of retrieving patient data includes retrieving data transmitted via a carrier wave.
16. A computer-implemented method for monitoring patient dosages, comprising:
5 retrieving dosing data, including times and amounts of medication prescribed for a patient;
retrieving patient data, including times and amounts of medication delivered to the patient;
determining evaluation data by analyzing the retrieved dosing and patient data to
10 determine compliance of the delivered medication to the prescribed medication;
and
displaying the evaluation data.
17. A memory device storing computer readable instructions for aiding a computer to monitor patient dosages of a medicine, comprising:
15 instructions for storing patient data, including the medication name and amounts of the medication prescribed for a patient;
instructions for retrieving patient data, including times and amounts of the medication delivered to the patient;
instructions for evaluating data by analyzing drug dispensing data and the patient data
20 to determine compliance of the delivered medication to the prescribed medication; and
instructions for displaying the evaluated data.
18. A computer system for monitoring patient dosages, comprising:
a processor for storing patient data, including a name of a medication and amounts of
25 the medication prescribed for a patient and for retrieving patient data, including times and amounts of the medication delivered to the patient, and evaluating data by
analyzing drug dispensing data and the patient data to determine compliance of the
delivered medication to the prescribed medication; and
a monitor for displaying the evaluated data.
19. The computer system of claim 18, further comprising a communications link linking the
30 processor to a remote device, wherein the retrieved patient data may be received from the remote device over the communications link.

20. The computer system of claim 19, wherein the retrieved patient data is received from the remote device over the communications link via a carrier wave.
21. The computer system of claim 20, further comprising an input device coupled to the processor, wherein the retrieved patient data may be received through the input device.
- 5 22. A method of graphically displaying drug compliance information, the method comprising the computer-implemented steps of:
receiving dosage data representing one or more quantities and one or more
administration times for delivering a drug;
10 receiving administration data representing one or more times when each of a plurality
of doses of the drug was delivered;
generating a graphical display of the drug compliance information on a display device,
wherein the graphical display comprises one or more elements that each
correspond to a time period;
15 displaying, within a first element among the elements, one or more icons that represent
each dose due within said first period; and
rendering each of the icons in one of a plurality of formats based on said dosage data
and said administration data.
23. The method of Claim 22, wherein the step of rendering includes the steps of:
20 determining whether a particular dose due within the first period was correctly
delivered based on said scheduling data and said administration data;
rendering a particular icon in a first format when the particular dose was incorrectly
delivered; and
rendering the particular icon in a second format when the particular dose was correctly
delivered.
- 25 24. The method of Claim 22, wherein the step of receiving administration data includes the
step of receiving data indicating an administration time for said particular dose, and
wherein the method further includes the steps of:
receiving data indicating a time period in which said drug should be delivered;
30 determining whether the particular dose was delivered within the time period; and
rendering the icon in a third format when the particular dose was delivered
within the time period.

25. The method of Claim 22, wherein the step of receiving administration data includes receiving data indicating an administration time for said particular dose, and wherein the method further includes the steps of:
receiving data indicating a time period within the administration time in which said drug
5 should be delivered;
determining whether the particular dose was delivered within the time period; and
rendering the icon as a particular format when the particular dose was delivered within the time period.
26. The method of Claim 22, further including the steps of:
10 displaying a graphical object;
displaying a second set of icons along an axis of the graphical object, in which the second set includes an icon for each dose of the drug delivered within a first period, and the position of each of the second set along the axis identifies when the respective dose was delivered.
- 15 27. The method of Claim 26, further including the step of said user selecting said first grid element associated with said first period.
28. The method of Claim 27, further including the steps of:
selecting a first icon of the second set of icons, wherein the first icon is associated with a first dose, wherein the first dose is associated with a first administration time; and
20 displaying additional information about the first dose, including the administration time.
29. The method of Claim 22, wherein the step of receiving dosage data includes the step of receiving dosage data from a dosage dispensing device.
30. The method of Claim 29, wherein the step of receiving dosage data from a dosage dispensing device includes the step of receiving dosage data from a portable medication
25 administration device.
31. The method of Claim 22, wherein the step of receiving administration data includes the step of receiving administration data from a portable medication administration device.
32. A method of generating data representing patient medication administration compliance, the method including the steps of:

receiving dosage data indicating parameter values for delivering a drug to a patient, wherein the parameter values specify one or more quantities and one or more administration times for delivering doses of the drug to the patient;

receiving administration data that indicates when each of a plurality of doses of the drug was administrated to the patient; and

generating data indicating a portion of the plurality of doses that was delivered according to the parameter values.

33. The method of Claim 32, wherein the step of generating data includes the steps of generating one or more values specifying a portion of said plurality of doses that was delivered within a specified time period of said administration times.

34. The method of Claim 32, further including the step of receiving data specifying a time period for which to generate said compliance data, the time period containing a plurality of days; and

wherein the step of generating data includes the steps of generating, for each day within a period, one or more values specifying a portion of the plurality of doses scheduled for the day that were delivered.

35. A method of managing the administration of drugs to a patient, the method comprising the steps of:

receiving dosage data that represents one or more administration quantities and one or more administration times for delivering doses of a drug to a patient;

transmitting, to a dosage dispensing device, data that specifies said one or more administration quantities and one or more administration times;

receiving administration data that indicates how each of a plurality of doses of the drug was administrated to the patient; and

storing the administration data in a memory device.

36. The method of Claim 35, further including the steps of:
receiving data specifying a lockout period that must elapse after delivering a dose before another dose is delivered to the patient; and
transmitting, to the dosage dispensing device, data that specifies the lockout period.

37. The method of Claim 36, wherein:
the step of receiving dosage data includes the steps of receiving data indicating a volume to deliver; and

the step of transmitting includes the steps of transmitting, to a dosage dispensing device, data that specifies said volume.

38. The method of Claim 35, wherein:
the step of receiving dosage data includes the steps of receiving data specifying a dose
of a liquid drug; and
the step of transmitting includes the steps of transmitting, to a dosage
dispensing device, data that specifies a dose of said liquid drug.

39. The method of Claim 38, further including the steps of:
receiving data indicating that a drug container has been removed from the dosage
dispensing device;
storing the data indicating that said drug container has been removed; and
reporting the data in a report of medication events.

40. The method of Claim 38, further including the steps of:
receiving data indicating that a drug container has been inserted into the dosage
dispensing device;
storing the data indicating that said drug container has been inserted; and
reporting the data in a report of medication events.

41. A method of managing administration of drugs to a patient, the method comprising the
steps of:
receiving data indicating administration times for a drug to be delivered to a
patient and a lockout period that must elapse after delivering a dose before another
dose is delivered to said patient; and
transmitting, to a dosage dispensing device, data that specifies the lockout period.

42. The method of Claim 41, wherein the step of transmitting includes the steps of
transmitting to a dosage dispensing device that dispenses a liquid.

FIG. 1




FIG. 2





FIG. 3



4 / 8




Fig. 4




Fig 5A




Fig. 5B





Fig. 5C



INTERNATIONAL SEARCH REPORT

 International application No.
 PCT/US98/22830

A. CLASSIFICATION OF SUBJECT MATTER

IPC(6) : G06F 15/42

US CL : 364/413.02

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS-SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

U.S. : 364/413.02

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practicable, search terms used)

APS

C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	US 5,347,453 A (MAESTRE) 13 September 1994, col. 5, lines 53-68 to col. 8, lines 1-7, col. 9, lines 47-68 to col. 12, lines 1-63, col. 13, lines 30-68 to col. 16, lines 1-50.	1-41
X	US 5,016,172 A (DESSERTINE) 14 May 1991, col. 2, lines 23-68 to col. 4, lines 1-6	1-17
X	US 4,839,806 A (GOLDFISCHER et al) 13 June 1989, col. 8, lines 17-68 to col. 9, lines 1-6.	17

☐ Further documents are listed in the continuation of Box C.
 ☐ See patent family annex.

* Special categories of cited documents:	*T* later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
A document defining the general state of the art which is not considered to be of particular relevance	*X* document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
B earlier document published on or after the international filing date	*Y* document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art
L document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)	*A* document member of the same patent family
O document referring to an oral disclosure, use, exhibition or other means	
P document published prior to the international filing date but later than the priority date claimed	

Date of the actual completion of the international search

06 JANUARY 1999

Date of mailing of the international search report

07 MAY 1999

 Name and mailing address of the ISA/US
 Commissioner of Patents and Trademarks
 Box PCT
 Washington, D.C. 20231

Facsimile No. (703) 305-3230

Authorized officer

JAMES TRAMMELL

Telephone No. (703) 308-9293

Joni Hill